

Invacare® 9000™ Topaz™ Wheelchair





Yes, you can.™

The Invacare® 9000™ Topaz™ wheelchair is designed specifically to suit the unique and growing needs of the bariatric population of individuals. The Topaz combines the strength and durability of a triple chrome-plated, carbon steel frame with the comfort of heavy-duty front riggings and arms for proper positioning. With a variety of seat widths and depths, frame reinforcements that include gusseted sides, caster journals and dual heavy-duty crossbraces, this chair provides not only ease in maneuverability but stability that will be trusted for years to come. Invacare 9000 Topaz wheelchairs meet RESNA standards for manual wheelchairs in product performance.



**A** The 9000 Topaz is available with a triple chrome-plated, carbon steel frame, along with a variety of frame finishes for customization.

## FEATURES

- Available in vinyl or nylon upholstery
- Durable low-maintenance, triple chrome-plated, carbon steel frame is long-lasting
- Reinforced frame includes gusseted sides, caster journals and dual heavy-duty crossbraces
- Multiple frame color variations available for customizing
- Heavy-duty front riggings provide proper positioning for feet with footrest and elevating legrest
- Heavy-duty arms designed to accommodate weight capacity and provide proper positioning
- Dual-axle positions allow for quick and simple seat-to-floor height variations, 17.5" and 19.5"
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty triple inner liners resist stretching of seat and back
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage

## Specifications

Dimensions	
Seat sizes	
Widths	20", 22", 24", 26", 28", 30"
Depths	18", 20"
Seat-to-floor heights	17.5", 19.5" 15.5" with different wheel and caster packages
Back heights	16", 17", 18"
Front riggings	
Swingaway footrest	
Swingaway padded elevating legrest	
Arm styles	
Full-length fixed height, padded, removable	
Desk-length fixed height, padded, removable	
Full-length adjustable height, padded, removable	
Desk-length adjustable height, padded, removable	
(All arm styles available conventional)	
Weight capacity	700 lb.
Product weight	82 lb.
Warranty	lifetime on frame and crossbraces

Invacare Corporation  
www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44036-2125  
(800) 333-6900

### Canada

5970 Chedworth Way  
Mississauga, Ontario  
L5R 3T9 Canada  
(800) 668-5324

Invacare is a registered trademark and Yes, you can., 9000 and Topaz are trademarks of Invacare Corporation. Specifications are subject to change without notification.