



Invacare® Compass™ XE
Adult Positioning Folding Wheelchair
Suggested Canadian Price List
Effective March 1st, 2005
Revised June 1st, 2007
***UPDATED - DENOTES CHANGES**

NAME _____ PHONE NO _____

ACCOUNT# _____ P.O. _____ DATE _____

ADDRESS _____

FRAME SELECTION

FRAME TYPE

CXE Folding, Depth Adjustable Extended Frame..... \$2,768

CXE CROSSBRACE TYPE

CONV Conventional Crossbrace..... No Charge

LSR Less Seat Rails..... No Charge

1 TRANSPORT TIE DOWN

TRRO Transport Ready Option..... No Charge

TRBKTS Transport Brackets & Pelvic Belt..... No Charge

TROMIT Less Transport Ready Option..... Less net (\$150)

FRAME DEPTH

MEDM Medium Side Frames (14" - 16")..... No Charge

LONG Long Side Frames (14" - 18")..... No Charge

XLONG Extra Long Side Frames (16" - 20")..... No Charge

SEAT SELECTIONS

SEAT WIDTH

W12 12" Wide..... No Charge

W13 13" Wide..... No Charge

W14 14" Wide..... No Charge

W15 15" Wide..... No Charge

W16 16" Wide..... No Charge

W17 17" Wide..... No Charge

W18 18" Wide..... No Charge

W19 19" Wide..... \$149

W20 20" Wide..... \$149

SEAT DEPTH

D12 12" Seat Depth..... No Charge

D13 13" Seat Depth..... No Charge

D14 14" Seat Depth..... No Charge

D15 15" Seat Depth..... No Charge

D16 16" Seat Depth..... No Charge

D17 17" Seat Depth..... No Charge

D18 18" Seat Depth..... No Charge

D19 19" Seat Depth..... No Charge

D20 20" Seat Depth..... No Charge

ARM SELECTIONS

ARM TYPE

2 AT99N Cantilever Non Locking Bolt On Arms \$244

3 FXFULL Conventional Fixed Ht., Full Lgth. Arms \$244

3 FXDESK Conventional Fixed Ht., Desk Lgth. Arms \$244

3 DESKADJCV Conventional Adjust. Ht., Desk Lgth. Arms \$348

3 FULLADJCV Conventional Adjust. Ht., Full Lgth. Arms \$348

2 AT55D Cantilever Locking Desk Lgth. Bolt On Arms \$348

2 AT55F Cantilever Locking Full Lgth. Bolt On Arms \$348

4 TAD T-Arm Adjustable Height, Desk Length \$348

4 TAF T-Arm Adjustable Height, Full Length \$348

Transport Ready Option	
*TRRO Must be ordered with conventional crossbraces, seat cushion and PERSONAL BACK™. TRRO includes installed transport brackets and wheelchair anchored pelvic belt.	
*TRBKTS includes installed transport brackets and wheelchair anchored pelvic belt.	

*** UPDATED - 06/01/2007 - Description Corrected**

Available Seat Widths / Depths for Compass XE							
	D14	D15	D16	D17	D18	D19	D20
W14	X	X	X	X	X	S	S
W15	X	X	X	X	X	S	S
W16	X	X	X	X	X	S	S
W17	X	X	X	X	X	S	S
W18	X	X	X	X	X	X	X
W19	S	S	X	X	X	X	X
W20	S	S	X	X	X	X	X

S = Must be ordered thru specials

*** UPDATED - 06/01/2007 - New Option**

*** UPDATED - 06/01/2007 - New Option**

*** UPDATED - 06/01/2007 - New Option**

*** UPDATED - 06/01/2007 - New Option**

Arm Heights	
AT99N	5" - 10"
Range adjusts in 1" increments	
AT55D	5" - 13"
AT55F	5" - 13"
Range adjusts in 1" increments	
DESKADJCV, FULLADJCV	10" - 14"
Range adjusts in 1" increments	
TAD	Low Range 6" - 10"
Range adjusts in 1/2" increments	
TAF	High Range 8" - 14"
Range adjusts in 1" increments	
All measurements from top of seat rail to top of armrest. On AT55D and AT55F, add 1" for upholstered armrests. For less seat rail chairs, add 1-1/2" to above measurements (top of upper frame rail to top of armrest).	

FOOTNOTES

- Transport Tie Down options not retrofittable to existing models, not field serviceable.
- Only available with FBH18, FBH20, FBH24, fixed back height.
- Only available with CONV, Conventional Crossbraces, and Seat Depths greater than 15". *** UPDATED - 06/01/2007**
- Bolt-On Arm Attachment.

T-ARM HEIGHT

- L6-10 Low Range 6" - 10" (T-Arm Height)..... No Charge
- H8-14 High Range 8" - 14" (T-Arm Height)..... No Charge

BACK SELECTIONS

BACK TYPE

- A1418 10° Bend Back w/Push Handles (14"-18")..... No Charge
- 5 1014 Sportster Back w/o Push Handles (10"-14")..... No Charge
- 5 S1216 Sportster Back w/o Push Handles (12"-16")..... No Charge
- A1216 Straight Back w/ Push Handles (12"-16")..... No Charge
- S1418 Straight Back w/ Push Handles (14"-18")..... No Charge
- FBH18 Fixed Back Height 18"..... No Charge
- FBH20 Fixed Back Height 20"..... No Charge

BACK UPHOLSTERY

- SLIPBACK Slip On Back Upholstery..... No Charge
- ADJTENS Adjustable Tension Back Upholstery..... \$133
- L8601 Omit Back Upholstery..... Less net (\$20)

6 RECLINER BACK

- REC18 Reclining Back 18" High..... \$743
- REC20 Reclining Back 20" High..... \$743
- REC24 Reclining Back 24" High..... \$743
- 7 RECEXT Reclining Headrest Pillow and Extension..... \$122

CASTER & WHEEL SELECTIONS

9 FRONT SEAT-TO-FLOOR HEIGHT

- FSTF1525 15.25" Front Seat-To-Floor Height..... No Charge
- FSTF1575 15.75" Front Seat-To-Floor Height..... No Charge
- FSTF1625 16.25" Front Seat-To-Floor Height..... No Charge
- FSTF165 16.5" Front Seat-To-Floor Height..... No Charge
- FSTF1675 16.75" Front Seat-To-Floor Height..... No Charge
- FSTF17 17" Front Seat-To-Floor Height..... No Charge
- FSTF1725 17.25" Front Seat-To-Floor Height..... No Charge
- FSTF175 17.5" Front Seat-To-Floor Height..... No Charge
- FSTF1775 17.75" Front Seat-To-Floor Height..... No Charge
- FSTF18 18" Front Seat-To-Floor Height..... No Charge
- FSTF1825 18.25" Front Seat-To-Floor Height..... No Charge
- FSTF185 18.5" Front Seat-To-Floor Height..... No Charge
- FSTF1875 18.75" Front Seat-To-Floor Height..... No Charge
- FSTF19 19" Front Seat-To-Floor Height..... No Charge
- FSTF1925 19.25" Front Seat-To-Floor Height..... No Charge
- FSTF1950 19.50" Front Seat-To-Floor Height..... No Charge
- FSTF1975 19.75" Front Seat-To-Floor Height..... No Charge
- FSTF20 20" Front Seat-To-Floor Height..... No Charge
- FSTF2025 20.25" Front Seat-To-Floor Height..... No Charge
- FSTF2050 20.50" Front Seat-To-Floor Height..... No Charge
- FSTF2075 20.75" Front Seat-To-Floor Height..... No Charge
- FSTF2125 21.25" Front Seat-To-Floor Height..... No Charge

CASTER OPTIONS

- AT900 3" Roller Blade Casters..... \$80
- 5ALUM 5" Aluminum Urethane Casters \$143
- AT901 5" Composite Urethane Casters No Charge
- AT902 6" Composite Urethane 3 Spoke Casters No Charge
- 6SKYWAY 6" Composite Urethane 6 Spoke Casters..... No Charge
- 1259 6" X 2" Semi Pneumatic Casters \$122
- AT903 8" X 1 1/4" Composite Urethane Casters No Charge
- 1228 8" X 1 3/4" Semi Pneumatic Casters Lt Grey..... \$140

FORK TYPE

- A1103 3" Fork (Single Hole Fork)..... No Charge
- A1105 5" Fork (3 Hole Fork)..... No Charge
- A1108 8" Fork (2 Hole Fork)..... No Charge

Front and Rear Seat-to-Floor Heights					
*Measured From Floor To Top Of Seat Rail					
Fork	3"	5"	6"	6"	8"
	AT900	5ALUM AT901	AT902 6SKYWAY	1259	1228 AT903
A1103 3" Fork	15.25" 15.75" 16.50" 17.00"	16.25" 16.75" 17.50" 18.00"	N/A	N/A	N/A
A1105 5" Fork	16.25" 16.75" 17.25" 17.50" 18.00" 18.50"	16.75" 17.25" 17.75" 18.00" 18.25" 18.50"	17.25" 17.75" 18.25" 18.50" 18.75" 19.00" 19.50"	N/A	N/A
A1108 8" Fork	N/A	N/A	N/A	18.00" 18.50" 19.00" 19.25" 19.50" 19.75" 20.25" 20.75"	19.00" 19.50" 20.00" 20.25" 20.50" 20.75" 21.25" 21.75"
All Dimensions Are Approximate					

FOOTNOTES

- 5 For push handles see Back Accessories section.
- 6 1" diameter anodized straight canes with push handles comes standard with back upholstery.
- 7 RECEXT must be chosen with reclining back selection.
- 8 For folding frame w/o seat rails, subtract 1 1/2" from the seat-to-floor height to obtain distance from floor to top of side frame.
- 9 Rear seat-to-floor height is set approximately 1" lower than the front seat-to-floor height.

10 WHEEL OPTIONS

- 20COM 20" Composite Wheels..... No Charge
- 22COM 22" Composite Wheels No Charge
- 11 24COM7 24" Composite 7-Spoke Wheels No Charge
- 20SP 20" Spoke Wheels No Charge
- 11 22SP 22" Spoke Wheels..... No Charge
- 24SP 24" Spoke Wheels No Charge

TIRE OPTIONS

- PNEU Pneumatic Tires No Charge
- PNEUFF Flat Free Inserts..... \$119
- 12 URETH Urethane (Treaded) Tires No Charge
- 13 PRIMO Primo Tires \$98
- 14 KIKB KIK® Black Tires..... \$98

MOUNTING HANDRIMS (SPOKE WHEELS ONLY)

- 15 IN Mounting Inside Handrims..... No Charge
- OUT Mounting Outside Handrims..... No Charge

HANDRIM OPTIONS

- ALUMPC Aluminum Anodized Handrims..... No Charge
- A1512 Black Plastic Coated Handrims \$107
- 16 AT1225-8 Plastic Coated w/8 Vertical Projections..... \$286
- AT1224-8 Plastic Coated w/8 Oblique Projections..... \$286
- 14 AT1225-12 Plastic Coated w/12 Vertical Projections..... \$286
- 14 NFITSD Natural Fit Handrims Standard Grip (pair)..... \$355
- 14 NFITSP Natural Fit Handrims Super Grip (pair)..... \$415
- 14 L8702 Omit Handrims..... Less net (\$15)

AXLE & CAMBER SELECTIONS

AXLE PLATE OPTIONS

- STDPT Standard Axle Plates..... No Charge
- MPOSAP Multi-Position Offset Axle Plate..... No Charge
- 17 EXTPT Extended/Amputee Axle Plates \$107

AXLES

- A1577 Quick Release Axles..... No Charge
- PAXLE Permanent Axles Standard
- A1573 Quad Release Axles..... \$152

CAMBER ANGLE

- C00 0° Camber..... No Charge
- C03 3° Camber..... No Charge
- C06 6° Camber..... No Charge

FRAME FINISH SELECTIONS

FRAME FINISH

CONTEMPORARY

- 24P Wet Black..... No Charge
- 60P Silver Metallic..... No Charge
- 71P Silver Vein..... No Charge
- 105P Midnight Blue..... No Charge
- 115P Black Prism..... No Charge

TOO HOT

- 30P Sunny Yellow..... No Charge
- 61P Electric Red..... No Charge
- 119P Bubblegum Pink..... No Charge
- 120P Grape Madness..... No Charge

WAY COOL

- 62P Electric Blue..... No Charge
- 104P Electric Teal..... No Charge
- 121P Lolly Pop Blue..... No Charge
- 122P Cosmic Blue..... No Charge
- 125P Emerald Green..... No Charge
- 127P Grasshopper Green..... No Charge

Knee-to-Heel Range							
For folding frames with seat rails							
70° Tapered		70°		70°MFX		Elevating	
MIN	MAX	MIN	MAX	MIN	MAX	MIN	MAX
16 1/4"	20 1/4"	14"	19"	8 1/4"	14"	15 3/4"	21 3/4"
90°		60°		60°MFX			
MIN	MAX	MIN	MAX	MIN	MAX		
5 3/4"	10 3/4"	15 1/4"	19"	11 1/2"	15 1/2"		

Range adjusts in 1" increments.
*Measurement taken from top of seat rail to back of footrest.

FOOTNOTES

- 10 Spoke wheels are standard with dome button axle.
- 11 Not available on medium frame set at 14" deep.
- 12 Not available with 20SP, 22SP or 24SP.
- 13 Only available on 24SP.
- 14 Only available on 24COM7 and 24SP.
- 15 Decrease overall width by 3/4".
- 16 Only available on 22COM, 24COM7 and 24SP.
- 17 Standard (no charge) with Reclining Back Option.

POSITIONING STRAPS

POSITIONING STRAPS

<input type="checkbox"/>	A1515	Push Button Buckle Seat Positioning Strap.....	\$79
<input type="checkbox"/>	A1321	Hook & Loop Seat Positioning Strap	\$61
<input type="checkbox"/>	9515	Push Button Chest Positioning Strap.....	\$79
<input type="checkbox"/>	7321	Hook & Loop Chest Positioning Strap.....	\$61

BACK ACCESSORIES

18	<input type="checkbox"/>	BP02	Back Pack.....	\$98
	<input type="checkbox"/>	1483	Non-Adjustable Split Handle Bars.....	\$146
	<input type="checkbox"/>	AASH	Adjustable Angle Rigidizing Handle Bar.....	\$424
19	<input type="checkbox"/>	RIGB	Flip Down Rigidizer Bar.....	\$146

SEAT ACCESSORIES

<input type="checkbox"/>	SP-01	Seat Pouch.....	\$78
<input type="checkbox"/>	L8651	Omit Seat Upholstery.....	Less net (\$15)

20 ARM ACCESSORIES

<input type="checkbox"/>	A1571	Fabric Clothing Protector.....	\$110
<input type="checkbox"/>	1671-SM	Rigid Side Guard, Small.....	\$189
<input type="checkbox"/>	1671-LG	Rigid Side Guard, Large.....	\$189
<input type="checkbox"/>	1671-SMH	Quad Rigid Side Guard, Small	\$189
<input type="checkbox"/>	1671-LGH	Quad Rigid Side Guard, Large.....	\$189

WHEEL LOCKS

21	<input type="checkbox"/>	HUMW	Hideaway Undermount Wheel Locks.....	\$139
	<input type="checkbox"/>	9691	Push to Lock Clamp On Wheel Locks.....	No Charge
	<input type="checkbox"/>	A1567	Pull to Lock Clamp On Wheel Locks.....	No Charge
22	<input type="checkbox"/>	A1598	Hill Holder Clamp On Wheel Locks.....	\$139

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/>	3421	3/4" To 1" Clamp Assembly.....	\$46	
<input type="checkbox"/>	A1569	Crutch & Cane Holder.....	\$61	
<input type="checkbox"/>	A1209	Tire Pump.....	\$58	
<input type="checkbox"/>	A1307	Spoke Protectors (for 20", 22", 24", 26").....	\$143	
<input type="checkbox"/>	A1580	Rear Anti-Tippers.....	\$102	
23	<input type="checkbox"/>	A1541	Caster Swivel Locks.....	\$123
	<input type="checkbox"/>	A1576	Wheel Lock Extensions (6")	\$59
	<input type="checkbox"/>		Touch Up Paint.....	\$43

FRONT RIGGING SELECTIONS

24 FRONT RIGGING

<input type="checkbox"/>	SWING1319	60° Swingaway Footrests.....	\$225	
<input type="checkbox"/>	SWING914	60° Swingaway Footrests MFX.....	\$351	
25	<input type="checkbox"/>	70MFX	70° Swingaway Footrests MFX	\$351
	<input type="checkbox"/>	70NTAPER	70° Swingaway Footrests.....	\$225
	<input type="checkbox"/>	70STAPER	70° Tapered Footrests.....	\$262
26	<input type="checkbox"/>	90SWING	90° Swingaway Footrests.....	\$231
	<input type="checkbox"/>	70LIFTNT	70° Lift Off Non-Taper Footrests.....	\$225
	<input type="checkbox"/>	70LIFTMFX	70° Lift Off MFX Footrests.....	\$351
	<input type="checkbox"/>	90LIFTMFX	90° Lift Off MFX Footrests.....	\$231
	<input type="checkbox"/>	AT5044	Elevating Legrests.....	\$529
	<input type="checkbox"/>	ATL4	Center Pivot Articulating Legrests.....	\$366

FOOTNOTES

- 18 Not available with Reclining Back Option.
- 19 Not available on 19" or 20" widths.
- 20 Only available with AT99N, AT55D or AT55F.
- 21 Not available with MEDM and/or D12-D16. *NEW- 06/01/2007
- 22 Not available with URETH or KIKB tires.
- 23 Only available on tubular footplates, AT900, 5ALUM, AT901, AT902, 6SKYWAY.
- 24 All front riggings are center pivot style. Footplates are attached directly to rigging tube when MFX style is selected.
- 25 Black finish only, standard with pull-to-lock wheel locks.
- 26 16" - 20" wide frames only.
- 27 Not available with 70MFX swingaway footrests.
- 28 AT5543 footplates are dorsi-plantar adjustable.
- 29 Not available with 70STAPER or 12" - 15" wide. * UPDATED - 06/01/2007

FOOTPLATE OPTION

27	<input type="checkbox"/>	1651	Flip-Up Composite Footplates.....	No Charge
28	<input type="checkbox"/>	AT5543	Adjustable Angle Flip-Up Footplates.....	\$122
29	<input type="checkbox"/>	A1564	Extra Large Aluminum Footplates.....	\$99
	<input type="checkbox"/>	TUBULAR	Individual Tubular Footplates.....	No Charge

MAX KNEE TO HEEL LENGTH

<input type="checkbox"/>	KTH14	14" Knee to Heel.....	No Charge
<input type="checkbox"/>	KTH1516	15"-16" Knee to Heel.....	No Charge
<input type="checkbox"/>	KTH1618	16"-18" Knee to Heel.....	No Charge

FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/>	A1337	Calf Strap.....	\$59	
30	<input type="checkbox"/>	IMPACTG	Impact Guards.....	No Charge
31	<input type="checkbox"/>	ALPT	Longer Pivot And Slide Tubes.....	\$123
	<input type="checkbox"/>	A1600BK	Heel Loops with Ankle Straps.....	\$79
	<input type="checkbox"/>	HLOOP	Heel Loops without Ankle Straps.....	\$74
	<input type="checkbox"/>	9342	'H' Strap.....	\$87

FOOTNOTES

- 30 Not available with elevating legrests.
- 31 2" longer slide tube used on composite, adjustable angle or tubular footplates only, not available with 70MFX swingaway footrests.

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEAT AND BACK CUSHIONS
SEE PRICE LIST FORM # 05-338C.**

- Basic Chair Includes:**
- Folding, Depth Adjustable Extended Frame
 - Crossbrace Type Selection
 - Transport Tie Down Options
 - Seat Width (W12 - W20)
 - Frame Depth Selection
 - Seat Depth (D12 - D20)
 - Arm Type (AT99N, FXDESK, FXFULL)
 - Back Type Selection (A1418, 1014, S1216, FBH18, FBH20)
 - Slip-On Back Upholstery
 - Frame Finish
 - Caster Options (AT902, AT903 or 6SKYWAY)
 - Seat to Floor Height Selection (FSTF1725 - FSTF2075)
 - Wheel Options
 - Tire Options (URETH)
 - Handrim Options (ALUMPC)
 - Axle Type (A1577 or PAXLE)
 - Axle Brackets (STDPT or MPOSAP)
 - Camber Angle Selection
 - Front Rigging (SWING1319, SWING914, 70NTAPER, 70MFX, 70STAPER OR 90SWING)
 - Flip-Up Composite Footplates (1651)
 - Seat Upholstery
 - Wheel Locks (9691 or A1567)

Specifications	
• Frame Style	Folding
• Seat Widths	12" - 20"
• Seat Depths	12" - 20"
• Overall Chair Width (closed)	Seat width plus 11.5"*
• Overall Chair Width (open)	Seat width plus 8.5"*
• Adjustable Height Back	14" - 18"
• Back Angle Adjustment	80° to 110° in 5° increments
• Camber Angles	0°, 3°, 6°
• Tubing Outside Diameter	1"
• Weight of Chair	19 lbs**/36 lbs***
• Warranty	Lifetime on frames and crossbraces
*0° Angles	
**18"x16", AT902, no rear wheels, no footrests, no anti-tippers, no wheel locks, no arms.	
***18"x16", AT902, 24SP, PNEU, 70NTAPER, seat upholstery, back upholstery, wheel locks, no arms.	

Weight Limitation 250 Pounds

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

Invacare Canada reserves the right to change prices and specifications without notice.
 HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4G 4Z4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)
 Regional Centre - 16769 Boul. Hymus, Kirkland, P.Q. H9H 3L4 Phones: (514) 630-6080 (800) 363-0137(P.Q.) Fax: (514) 630-6120
 DISTRIBUTION CENTRE - #10 -1595 Cliveden Avenue, Delta, BC V3M 6M2 Phone:(800)668-5324(Canada) Fax: (800) 668-5478(Canada)
 Invacare®, and Yes, you can., are trademarks of Invacare Corporation.
 WWW.INVACARE.CA