

BREEZY[®]

S E R I E S



BREEZY



WE LISTEN:

The Breezy philosophy is simple—develop products that don't just meet your needs, but exceed them. Status quo is not even in our vocabulary. Every aspect of our product is designed, built and sold with the people using it in mind. We even integrate our designs around family members, caregivers and suppliers, who have just as much concern over the quality, dependability and operation of the product as you do.

ATTENTION TO DETAIL:

Every concept, design and proto-type is highly scrutinized by our seasoned engineering staff. This ensures you receive a product that is easy to operate, increases your mobility and improves your quality of life.

APPEARANCE:

You take pride in how you look. The clothes, colors and styles you choose all reflect your personality and mood. Why should your chair be any different? We take pride in making sure our Breezy line of chairs meets those needs. Our color selection, ergonomic style and clean appearance don't happen by accident. We conscientiously incorporate these into everything we do.

ROBOTIC WELDING:

Our robotic welders eliminate all chances of "human error" with their precision performance and ability to produce reliable wheelchair frames and components. These precision machines are the most technologically advanced welders available and guarantee every component meets our stringent standards.

PARTS COMMONALITY:

Sunrise produces the highest quality components at an affordable price. The key to succeeding at this has been our parts commonality. Breezy wheelchairs use key components and superior designs from our elite line of Quickie® chairs. Caster journals, armrests, back canes and footrests are just some of the components that are integrated into a family of products designed to exceed your expectations.

BECAUSE LIFE MATTERS:

Visiting with family and friends, finishing that special project, helping your neighbor, playing with your children or grandchildren—this is what life's all about. We are committed to providing you with products that meet any need at any time. Because life matters.



600

Breezy 600

Breezy 600 is the lightest weight feature-packed chair in its class, making it the ideal solution for foot propellers and those needing environmental access. Ten different seat height adjustments allow this chair to fit every client.

BREEZY 600

High-Strength Lightweight



OPTIONS



One-Arm Drive positive lock strut system increases propelling efficiency and is easily removed for folding. Fits Breezy 503, 510 and Breezy 600 models.



Standard, hemi and ultra-hemi configurations (shown above).

Breezy

Breezy 500/510

With its seemingly limitless number of adjustments, configurations and features, the Breezy 510 is the ultimate high-strength, lightweight wheelchair.

BREEZY 510



High-Strength Lightweight

BREEZY 500



Economic value and outstanding quality make the Breezy 500 an excellent choice for virtually any situation.

FEATURES



Our **New Offset Axle** shifts the rear wheel forward to improve handrim access and traction and to reduce caster load. These changes increase efficiency, function and independence.

Our **New optional Amputee Axle** (not shown) offsets the rear wheel almost 3" back from the standard position providing an extremely stable chair.

OPTION



The Breezy 510 Recliner reclines to a full 180° for maximum pressure relief, positioning and comfort. The offset headrest is ergonomically designed for proper support.

Breezy

Breezy 503

The Breezy 503 is a value-packed, lightweight wheelchair offering unsurpassed quality and performance. With its 17" fixed back height, offset axle plate and parts commonality with the Breezy 500, 510, and 600, the Breezy 503 is the perfect choice.

BREEZY 503



Lightweight

Breezy Features & Options

- High Strength/
Lightweight
- Lightweight



Description	Breezy 600	Breezy 510	Breezy 500	Breezy 503
• FRAME				
Width	14" - 20"	14" - 20"	16", 18"	16" - 18"
Depth	16", 18"	16", 18"	16"	16"
Front Seat Ht	13" - 20.5"	13" - 20.5"	17.5" & 18.5"	17.5" & 18.5"
Backrest Ht	16" - 19"	17" - 19"	17" - 19"	17" Fixed
Reclining Backrest (180°)		Opt. (w/Adj. Pivot Pt.)		
Seat Sling	Padded & Non-Padded Nylon	Padded & Non-Padded Nylon	Padded & Non-Padded Nylon	Padded & Non-Padded Nylon
Frame Colors	Forest Green, Silver Texture & Blue	Silver Texture, Blue & Chrome	Silver Texture	Silver Texture
• ARMREST				
Desk Length	Padded	Padded	Padded	Padded
Full Length	Padded	Padded	Padded	Padded
Fixed Non-Removable				
Dual Post Removable		Recliner Only		
Flip Back	Standard	Standard	Standard	Standard
Height Adjustable	Option	Option	Option	Option
Tall Fixed Ht (+ 2.25")				
• HANGERS				
Fixed Non-Removable				
60 Degree Swing-Away	Option	Option	Option	Option
70 Degree Swing-Away	Option	Option	Option	Option
Elevating Legrest	Option	Option	Option	Option
Articulating Legrest	Option	Option	Option	Option
Footplate Style	Composite, Aluminum & Angle Adjustable	Composite & Aluminum	Composite	Composite
Heel Loops	Standard	Standard	Standard	Standard
• REAR AXLE				
Fixed Position Axle				
Forward Offset Axle	Standard	Standard	Standard	Standard
Amputee Axle		Option	Option	Option
One Arm Drive	Option	Option	Option	Option
Quick-Release Axle	Option	Option	Option	Option
• REAR WHEEL				
Size	20", 22" & 24"	20", 22" & 24"	24"	24"
Tire Type	Low Profile, Full Profile, Pneumatic, Pneumatic/ Airless	Low Profile, Full Profile, Pneumatic, Pneumatic/ Airless	Low Profile	Low Profile
Handrim Options	Composite, Aluminum Painted, Aluminum Anodized, Plastic Coated & Projection	Composite, Aluminum Painted, Aluminum Anodized, Plastic Coated & Projection	Composite, Aluminum Painted & Plastic Coated	Composite, Aluminum Painted, Aluminum Anodized, Plastic Coated
• WHEEL LOCKS	Push-to-Lock, Pull-to-Lock	Push-to-Lock, Pull-to-Lock	Push-to-Lock, Pull-to-Lock	Push-to-Lock, Pull-to-Lock
• WEIGHT CAPACITY	250 lbs.	250 lbs.	250 lbs.	250 lbs.
• WEIGHT OF CHAIR (LBS.)	28	32	<34	34.5
• ACCESSORIES				
6" Wheel Lock Ext. Handles	Option	Option	Option	Option
Anti-Fold Device				
Anti-Tippers	Option	Option	Option	Option
Calf Pad Cover	Standard	Standard	Standard	Standard
Chart Pocket				
Crutch Holder	Option	Option	Option	Option
Grade Aids	Option	Option		
Telescoping IV Pole				
Oxygen Tank Holder	Option	Option	Option	Option
Oxygen Tank Holder/IV Pole				
Positioning Belt	Option	Option	Option	Option
Wheelchair Tray	Option	Option	Option	Option

Accessories



Breezy 6" wheel lock extension with bungee cord



Breezy 500/510 Oxygen tank holder

BREEZY WHEELS



1

2

3

4

Wheel Description

- 1 Breezy 510/600 Low Profile Poly Tire with painted handrim
- 2 Breezy 510/600 Pneumatic Tire with plastic coated handrim
- 3 Breezy 500/510/600 Low Profile Poly Tire with composite handrim
- 4 Breezy 510/600 Full Profile Poly Tire with aluminum handrim

* Not all accessories shown. Please consult order form for complete list.