



Quickie

Accessorize Your Life!



www.sunrisemedical.com

Table of Contents

- Miscellaneous AccessoriesFRONT
- Swing-Away Hanger Options1
- Footplates2
- Armrest Options3
- Seat and Back Options4
- Caster Tires and Forks5
- Wheels, Tires and Rims6
- Axels and Axle Plates7
- Wheel Locks and Handrims8
- Miscellaneous AccessoriesBACK



▲ RAGGS MUSIC CD

◀ B. MAX

▼ BACK PACK



▲ TRAVEL WHEELS



▲ SEAT POUCH



▲ H LEG STRAP



▲ FRONT AND BACK IMPACT GUARDS



▲ SINGLE LEG STRAP

SWING-AWAY HANGER OPTIONS



▲ **70°, 80°, 90° FRONT MOUNT**
Popular option for kid's chairs. Designed for use with our flip-up footplates and angle adjustable footplates.

▲ **70°, 80°, 90° HEAVY-DUTY LIFT-OFF HANGER**
This heavy-duty hanger securely locks into place to meet the needs of the most demanding individuals.

▲ **60°, 70° SWING-AWAY HANGER**
Swings out of the way for transfers. Available on both rigid and folding chairs. 70° gives you a shorter overall chair turning radius.



◀ **ELEVATING LEGREST HANGER**
This low-profile legrest minimizes interference and provides comfortable positioning at an economical price. Available on manual and power chairs.



◀ **KIDS ELR LEGREST HANGER**
Great fit for kids and anyone requiring a shorter leg length.



◀ **ARTICULATING LEGREST HANGER**
This legrest assists the movement of your leg and provides support throughout the range of motion. Available on manual and power chairs.




◀ **70° V HANGER**
This sleek "V" swing-away has a streamlined design for maximum performance and style.

FOOT PLATES

PLATFORM FOOTPLATE

The angle-adjustable footplate is specially designed for chairs with a front frame bend of 70° or more.

**▲ 90° FOOTBOARD**

Footboard provides a large flat 90° foot surface.

**▲ FLIP-UP FOOTPLATES WITH FOAM COVERS**

These lightweight footplates with foam covers provide traction and keep feet from slipping.

**▲ ADAPTIVE PLATFORM SPACERS**

Available in 2", 4" and 6" spacers accommodates shorter, lower lengths.

**▲ ALUMINUM FOOTPLATES**

These machined aluminum footplates add 4" of adjustments in .5" increments.

**▲ KIDS' ANGLE-ADJUSTABLE FOOTPLATES**

Designed for use with front mount hangers, these flip-up footplates are angle-adjustable.

**▲ EXTENDED FOOTPLATES**

These machined aluminum footplates add 3" of forward extension to Quickie composite footplates.

**▲ COMPOSITE FLIP-UP FOOTPLATES**

Lightweight with a large, flat surface, these footplates easily flip-up and out of the way. Available on manual and power chairs.

**▲ ONE-PIECE FLUP-UP FOOTPLATE**

This footplate allows folding chairs to accept a rigid footplate and rigid frames to have flip-up footplates. Just unlatch and swing-away to fold.

**▲ A - RIGID FOOTREST
B - Ti FOOTPLATES with COVER****▲ CONTRACTUAL FOOTREST**

Ultimately adjustable to accommodate asymmetrical positioning of the legs and feet. Multiple axis adjustment allows the footplates to be positioned just where you need them.

**▲ ADULT ANGLE-ADJUSTABLE FOOTPLATES**

These footplates are angle- and depth-adjustable to accommodate multiple foot positions. Available on manual and power chairs.


**▲ LOCKING ANGLE-ADJUSTABLE FOOTPLATES**

Angle adjustable hinge provides depth and angle adjustment. Footplates lock together for additional strength.

ARMREST OPTIONS



- ◀ **SINGLE POST HEIGHT-ADJUSTABLE ARMREST**
Designed for easy adjustment, this detachable padded armrest can be moved up and down with ease. Available with either a 10" desk-length pad or 14" full-length pad. Available on manual and power chairs in either pediatric, kids or adult sizes.



- ▶ **FLIP-UP ARMREST**
This unique armrest flips up and down smoothly and easily (Zippie product line).



- ▶ **SWING-AWAY ARMREST**
This padded tubular armrest lifts up and swings out of the way during transfers.



- ▶ **KIDZ ARMREST**
With a smaller receiver, this armrest can be placed in a multitude of positions on the chair.

- ▶ **LENGTH ADJUSTABLE LOCKING FLIP-UP ARMREST**



- ▶ **ADJUSTABLE LOCKING FLIP-UP ARMREST**
This height- and angle-adjustable armrest can be flipped up and back easily, and features a trigger lock to secure the arm in place.



- ▶ **DUAL-POST HEIGHT-ADJUSTABLE ARMREST**
Dual-post for added stability. Height-adjustable for a good fit. Available with desk-length (10") or full-length (14") pads on some manual and most power chairs.




- ▶ **FLIP-BACK HEIGHT-ADJUSTABLE ARMREST**
This armrest is removable and height-adjustable for a perfect fit. Available in desk- and full-length pad.

- ▶ **FLIP-BACK ARMREST**
For easier transfers, this armrest pivots out of the way and locks back in place. Available in desk and full-length pads.



SEAT AND BACK OPTIONS



- ▲ **FOAM CUSHIONS**
Cushions have washable covers and are available in 2" thickness.



- ▲ **LIFTING STRAPS**
To assist when going up stairs or over curbs.



- ▲ **ADJUSTABLE UPHOLSTERY**
Tension-adjustable back support provides personalized comfort.



STANDARD GP UPHOLSTERY

◀ QUICKIE GP PUSH HANDLES



◀ QUICKIE 2 STYLE PUSH HANDLES

STANDARD Q2 UPHOLSTERY

- ◀ **SINGLE STROLLER HANDLE**
A convenient option for Quickie Kidz chair. This single stroller handle easily snaps into place.



- ▲ **SWING-AWAY ADJUSTABLE STROLLER HANDLES**
Height-adjustable within a 6" range, these stroller handles fold flat for easy transportability.



- ▲ **STROLLER HANDLES**
These extended push handles are the ideal compliment to Quickie kids' chairs. Stroller handles come in standard, medium and tall heights.



- ◀ **STROLLER HANDLE EXTENSION**
Stroller handle extensions easily snap in and out of place to provide a comfortable grip at any angle.

CASTER TIRES AND FORKS



- ◀ **3" & 4" MICRO CASTERS**
These lightweight casters provide increased on-court performance and maximum foot clearance.



- ◀ **3" & 4" LIGHTED CASTERS**
Lighted for added safety and increased visibility. And don't forget the fun value!



- ◀ **4" & 5" LOW-PROFILE CASTERS**
These lightweight casters offer a responsive ride and small turning radius. Available with aluminum or composite rims.



- ◀ **POLYURETHANE CASTERS**
Maintenance-free polyurethane casters are available in 6" and 8" sizes and aluminum or composite rims.



- ◀ **6" & 8" PNEUMATIC CASTERS**
These casters have air-filled inner tubes for a softer, smoother ride. Available in 6" and 8" sizes, and aluminum or composite rims.



- ◀ **A - Polyurethane Casters**
(front view)
- ◀ **B - 1.5" wide Poly Casters**
(front view)
Available sizes:
4", 5" and 6"
- ◀ **C - 5" Soft Roll**
(front view)



- ◀ **8" x 2" PNEUMATIC CASTERS**
8"x2" pneumatic casters for added flotation over rough terrain.



- ◀ **SLOTTED FORK**
Available sizes -
3.25", 4.25", 5.25",
6", and 7".



- ◀ **SLOTTED QUICK RELEASE FORK**
Available sizes -
5.25", 6", and 7".



- ◀ **CASTER PIN LOCK**



- ◀ **FROG LEGS SUSPENSION FORK**



- ◀ **TI CHAIR FORK**
Available in
short or tall.



- ◀ **XTR FORK**
One fork with
two stem bolt
options.



- ◀ **MULTI-POSITIONED CASTER FORK**
Vertically aligned
adjustment
holes. This
fork provides greater
range of seat height.



- ◀ **STEM BOLT with 1/2" SPACER**
Also available
with 3/4"
spacer option.

**12" WHEEL**

These wheels are half the size of most wheels and are positioned below the chair's seat out of the child's way. (Zippie TS, Quickie TS, Quickie Kidz.) Also available: sizes 16' and 18"

A. STANDARD WHEELRIMS

Lightweight and rigid for real performance. Available sizes: 20", 22", 24", and 26"

B. BASKETBALL WHEEL

This 24" or 26" high performance wheel is the choice of top Team Quickie athletes. Heat treated or standard rim, 2.0 mm spokes and a hi-lo flange hub.

C. TENNIS WHEEL

This 24" or 26" high performance wheel is the choice of top Team Quickie athletes. Heat treated or standard rim, 1.8 mm spokes and a hi flange hub.

D. SPINERGY WHEEL


The Ultimate Sports and Everyday Wheel! Ultra light weight, one pound lighter per set than other wheels. Spoke colors available: Black, Yellow, and White.

E. 22" MOUNTAIN WHEEL

Using our 22" wheel and wide knobbie tires you can change your wheels without effecting the seat height. (Wheel with knobbie is the same height as our and 24" std. and 24" heat treated.)

F. SIX SPOKE MAG WHEEL

Low-maintenance molded one-piece wheel. Available sizes: 20", 22", and 24"

**A. LATEX TUBULAR TIRES**

These ultra-lightweight high-pressure tubulars are available in 24" and 26" sizes.

B. LOW-PROFILE POLYURETHANE TIRE

This airless low-maintenance tire offers flat-free convenience. Available in 20", 22" and 24" sizes.

C. KEVLAR TIRE

Lightweight, puncture-resistant tire.

D. 24" KNOBBIE TIRE FOR STANDARD RIMS

This all-terrain tire features deep tread for additional traction. For use with standard rims.

E. PNEUMATIC TIRES

Available in 20", 22", 24" and 26" sizes, these air-filled tires offer a cushioned ride.

F. FULL-PROFILE POLYURETHANE TIRES

Available in 22" and 24" sizes, these airless tires provide flat-free convenience in a standard size tire.

G. 24" KNOBBIE TIRE FOR HEAT-TREATED RIMS

This all-terrain tire features deep tread for additional traction. For use with heat-treated rims.

H. HIGH-PRESSURE CLINCHER TIRES

Available in 24" and 26" sizes, this air-filled tire offers high performance for the everyday riders.

I. TURBO TREAD TIRES

Super-lightweight, high-pressure clincher tires, for use with replaceable tubes. Available in 24" and 26" sizes.

J. AIRLESS INSERT

This flat-proof solid insert replaces the innertube. Available in 12", 20", 22" and 24" sizes. Also available for power chairs with 8"x2" or 10" casters and 20" drive wheels.

A B C D E F G H I J




AXLES & AXLE PLATES



- A - AXLE PLATE M6
- B - AXLE PLATE QTS
- C - AXLE PLATE QLX
- D - AXLE PLATE AMPUTEE DOUBLE SLOT INDEX QGP
- E - AXLE PLATE Q2 - Q2HP AND ZIPPIE2
- F - AXLE PLATE DOUBLE SLOT INDEX QGP
- G - CURVED AXLE PLATE Q2 - Q2HP
- H - AXLE PLATE QGP
- I - P AXLE BRACKET QLX

- ▼ D - STEEL, QUICK-RELEASE AXLES
Quick-release mechanism allows wheels to be removed quickly and easily.
- ▼ E - BIG BUTTON AXLES
Protected from impact and a larger surface to press for quick release.
- ▼ F - AXLE SLEEVE

- ▶ A - THREADED AXLE
For riders who don't need removable wheels.
- ▶ B - QUAD-RELEASE AXLES
A special looped quad-release nut is attached to quick-release axles. Ideal for riders with limited hand function.
- ▶ C - TITANIUM QUICK-RELEASE AXLES
These lightweight durable axles are perfect for high performance chairs.



ONE-ARM DRIVE
This dual handrim unit can be mounted on chair's left or right side, is easily folded and adds less than 5 pounds to the weight of the chair. ▼



- ▶ Ti CAMBER TUBES
Camber options: 0°, 2°, 4°, and 8°
- ▶ CAMBER TUBES
0°, 4°, 8°, and 12°
- ▶ R2 CAMBER TUBE
Available with single or dual camber inserts.



- A - AXLE BRACKET (12") ZIPPIE TS
- ▶ B - AXLE BRACKET ZIPPIE TS
- C - EXTENDED AXLE BRACKET ZIPPIE TS





▲ FOOT LOCK

This foot lock easily locks both wheels since it requires only one foot to depress the lever. Accessible from either the rear or the side of the chair.



▲ HIGH MOUNT LOCK

Available in push-to-lock or pull-to-lock.



▲ HIGH-LOCK MOUNT EXTENSION HANDLES

These extension handles are available in 6" or 9" lengths.



▲ HIGH-MOUNT SCISSOR LOCK

This lock is easy to reach and swings out of the way.



▲ SHORT-THROW WHEEL LOCK

Provides maximum locking power with minimum force applied.



▲ GRADE AID

Prevents the chair from rolling back when going up an incline.



▲ UNILATERAL WHEEL LOCK

Allows rider to easily engage both wheel locks with one hand.



▲ SPRING LOADED WHEEL LOCK

This design requires less force to engage.

▼ PLASTIC-COATED HANDRIM

Coated with polyurethane for better grip.

▼ ALUMINUM HANDRIM

This standard handrim features seamless weld construction for strength and durability.

▼ VERTICAL PROJECTIONS

Available with 6, 8, 10 or 12 projections depending on wheel size.

▼ OBLIQUE PROJECTIONS

Available with 6, 8, 10 or 12 projections depending on wheel size.



MISCELLANEOUS ACCESSORIES



◀ SIDE GUARDS



▲ CADDY



▲ BACK POST STABILIZER BAR



▲ TOOL KIT



◀ CATCH ALL



▶ CRUTCH HOLDER



▲ TI TOOL KIT



▲ STRAPS AND POSITIONING BELTS

▼ TRAY TABLE



◀ SPOKE GUARDS



