

QUICKIE 2
FORK/CASTER HEIGHTS
(From floor to top of seat sling)

		CASTER	4" POLY	4" LIGHTED	5" SOFT	5" POLY	5"x 1-1/2" POLY	6" PNEU	6" SEMI PNEU (POLY)	6" x 1-1/2" POLY	8" POLY	8" PNEU	8" X 2" POLY	8" X 2" PNEU
		P/N TIRE W/CASTER	400115	400108	400095	400155	400148	400156	400154	400157	400204	400206	400702	400700
		ACTUAL DIAMETER	3.97	3.95	5.00	4.88	4.98	6.06	5.69	5.73	7.70	7.90	7.50	7.50
FORK	SPACER	FORK HOLE												
5-1/4"	STD	Top	16.69	16.68	17.20	17.14	17.19	17.73	17.55	17.57	N/G	N/G	N/G	N/G
5-1/4"	STD	Bottom	17.61	17.60	18.13	18.07	18.12	18.66	18.47	18.49	19.48	19.58	N/G	N/G
5-1/4"	3/4"	Top	17.44	17.43	17.95	17.89	17.94	18.48	18.30	18.32	N/G	N/G	N/G	N/G
5-1/4"	3/4"	Bottom	18.36	18.35	18.88	18.82	18.87	19.41	19.22	19.24	20.23	20.33	N/G	N/G
5-1/4"	1-1/2"	Top	18.19	18.18	18.70	18.64	18.69	19.23	19.05	19.07	N/G	N/G	N/G	N/G
5-1/4"	1-1/2"	Bottom	19.11	19.10	19.63	19.57	19.62	20.16	19.97	19.99	20.98	21.08	N/G	N/G
5-1/4"	QR	Top	16.81	16.80	17.33	17.27	17.32	17.86	17.67	17.69	N/G	N/G	N/G	N/G
5-1/4"	QR	Bottom	17.74	17.73	18.25	18.19	18.24	18.78	18.60	18.62	19.60	19.70	N/G	N/G
6"	STD	Top	17.44	17.43	17.95	17.89	17.94	18.48	18.30	18.32	19.30	19.40	N/G	N/G
6"	STD	Bottom	18.36	18.35	18.88	18.82	18.87	19.41	19.22	19.24	20.23	20.33	N/G	N/G
6"	3/4"	Top	18.19	18.18	18.70	18.64	18.69	19.23	19.05	19.07	20.05	20.15	N/G	N/G
6"	3/4"	Bottom	19.11	19.10	19.63	19.57	19.62	20.16	19.97	19.99	20.98	21.08	N/G	N/G
6"	1-1/2"	Top	18.94	18.93	19.45	19.39	19.44	19.98	19.80	19.82	20.80	20.90	N/G	N/G
6"	1-1/2"	Bottom	19.86	19.85	20.38	20.32	20.37	20.91	20.72	20.74	21.73	21.83	N/G	N/G
6"	QR	Top	17.56	17.55	18.08	18.02	18.07	18.61	18.42	18.44	19.43	19.53	N/G	N/G
6"	QR	Bottom	18.49	18.48	19.00	18.94	18.99	19.53	19.35	19.37	20.35	20.45	N/G	N/G
7"	STD	Top	18.44	18.43	18.95	18.89	18.94	19.48	19.30	19.32	20.30	20.40	N/G	N/G
7"	STD	Bottom	19.36	19.35	19.88	19.82	19.87	20.41	20.22	20.24	21.23	21.33	N/G	N/G
7"	3/4"	Top	19.19	19.18	19.70	19.64	19.69	20.23	20.05	20.07	21.05	21.15	N/G	N/G
7"	3/4"	Bottom	20.11	20.10	20.63	20.57	20.62	21.16	20.97	20.99	21.98	22.08	N/G	N/G
7"	1-1/2"	Top	19.94	19.93	20.45	20.39	20.44	20.98	20.80	20.82	21.80	21.90	N/G	N/G
7"	1-1/2"	Bottom	20.86	20.85	21.38	21.32	21.37	21.91	21.72	21.74	22.73	22.83	N/G	N/G
7"	QR	Top	18.56	18.55	19.08	19.02	19.07	19.61	19.42	19.44	20.43	20.53	N/G	N/G
7"	QR	Bottom	19.49	19.48	20.00	19.94	19.99	20.53	20.35	20.37	21.35	21.45	N/G	N/G
8X2"	STD	TOP	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G	19.32	19.32
8X2"	STD	Bottom	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G	20.25	20.25
Multi-Short	STD	Short 4	N/G	N/G	16.38	16.32	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G
Multi-Short	STD	Short 3	16.36	16.35	16.88	16.82	N/G	17.41	17.22	N/G	N/G	N/G	N/G	N/G
Multi-Short	STD	Short 2	16.86	16.85	17.38	17.32	17.37	17.91	17.72	N/G	N/G	N/G	N/G	N/G
Multi-Short	STD	Short 1	17.36	17.35	17.88	17.82	17.87	18.41	18.22	18.24	N/G	N/G	N/G	N/G
Multi-Short	3/4"	Short 4	N/G	N/G	17.13	17.07	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G
Multi-Short	3/4"	Short 3	17.11	17.10	17.63	17.57	N/G	18.16	17.97	N/G	N/G	N/G	N/G	N/G
Multi-Short	3/4"	Short 2	17.61	17.60	18.13	18.07	18.12	18.66	18.47	N/G	N/G	N/G	N/G	N/G
Multi-Short	3/4"	Short 1	18.11	18.10	18.63	18.57	18.62	19.16	18.97	18.99	N/G	N/G	N/G	N/G
Multi-Short	1-1/2"	Short 4	N/G	N/G	17.88	17.82	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G
Multi-Short	1-1/2"	Short 3	17.86	17.85	18.38	18.32	N/G	18.91	18.72	N/G	N/G	N/G	N/G	N/G
Multi-Short	1-1/2"	Short 2	18.36	18.35	18.88	18.82	18.87	19.41	19.22	N/G	N/G	N/G	N/G	N/G
Multi-Short	1-1/2"	Short 1	18.86	18.85	19.38	19.32	19.37	19.91	19.72	19.74	N/G	N/G	N/G	N/G
Multi-Tall	STD	Tall 5	N/G	N/G	N/G	N/G	N/G	18.41	18.22	N/G	19.23	19.33	N/G	N/G
Multi-Tall	STD	Tall 4	N/G	N/G	18.38	18.32	N/G	18.91	18.72	N/G	19.73	19.83	N/G	N/G
Multi-Tall	STD	Tall 3	18.36	18.35	18.88	18.82	18.87	19.41	19.22	N/G	20.23	20.33	N/G	N/G
Multi-Tall	STD	Tall 2	18.86	18.85	19.38	19.32	19.37	19.91	19.72	N/G	20.73	20.83	N/G	N/G
Multi-Tall	STD	Tall 1	19.36	19.35	19.88	19.82	19.87	20.41	20.22	20.24	21.23	21.33	N/G	N/G
Multi-Tall	3/4"	Tall 5	N/G	N/G	N/G	N/G	N/G	19.16	18.97	N/G	19.98	20.08	N/G	N/G
Multi-Tall	3/4"	Tall 4	N/G	N/G	19.13	19.07	N/G	19.66	19.47	N/G	20.48	20.58	N/G	N/G
Multi-Tall	3/4"	Tall 3	19.11	19.10	19.63	19.57	19.62	20.16	19.97	N/G	20.98	21.08	N/G	N/G
Multi-Tall	3/4"	Tall 2	19.61	19.60	20.13	20.07	20.12	20.66	20.47	N/G	21.48	21.58	N/G	N/G
Multi-Tall	3/4"	Tall 1	20.11	20.10	20.63	20.57	20.62	21.16	20.97	20.99	21.98	22.08	N/G	N/G
Multi-Tall	1-1/2"	Tall 5	N/G	N/G	N/G	N/G	N/G	19.91	19.72	N/G	20.73	20.83	N/G	N/G
Multi-Tall	1-1/2"	Tall 4	N/G	N/G	19.88	19.82	N/G	20.41	20.22	N/G	21.23	21.33	N/G	N/G
Multi-Tall	1-1/2"	Tall 3	19.86	19.85	20.38	20.32	20.37	20.91	20.72	N/G	21.73	21.83	N/G	N/G
Multi-Tall	1-1/2"	Tall 2	20.36	20.35	20.88	20.82	20.87	21.41	21.22	N/G	22.23	22.33	N/G	N/G
Multi-Tall	1-1/2"	Tall 1	20.86	20.85	21.38	21.32	21.37	21.91	21.72	21.74	22.73	22.83	N/G	N/G
Frog-Leg	2.5	Frog 3	16.79	16.78	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G
Frog-Leg	2.5	Frog 2	17.29	17.28	17.80	17.74	17.79	N/G	N/G	N/G	N/G	N/G	N/G	N/G
Frog-Leg	2.5	Frog 1	17.79	17.78	18.30	18.24	18.29	18.83	18.65	18.67	N/G	N/G	N/G	N/G
Frog-Leg	4	Frog 3	18.29	18.28	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G	N/G
Frog-Leg	4	Frog 2	18.79	18.78	19.30	19.24	19.29	N/G	N/G	N/G	N/G	N/G	N/G	N/G
Frog-Leg	4	Frog 1	19.29	19.28	19.80	19.74	19.79	20.33	20.15	20.17	N/G	N/G	N/G	N/G

For Hemi, subtract 2.125"

Xtender Power Assist Wheel-Caster & Fork option only available with Dimensions that are highlighted in white

N/G = no go