



XTENDER



"Pushing forward in life is easy with the Quickie Xtender power assist. No reasonable obstacle is too tough. The Quickie Xtender will help you make the grade! Game On!"

Marilyn Hamilton

Marilyn Hamilton, Founder of Quickie Wheelchairs



A successful design will accomplish a clearly defined mission which, in the case of the Quickie Xtender, is to offer xtended power when pushing your chair. Whether you are going to the mall, going to the theater or just taking a long stroll in the evening with your family and friends, the Xtender will truly extend life's little pleasures.



www.sunrisemedical.com

XTENDER

With the Quickie Xtender you can xtend the distance you travel, xtend the efficiency of your push and xtend your transportation options.

Quick Disconnect Axle Release

Wheel Synchronization

Transportability



XTENDER SPECIFICATIONS:

Torque multiplier	1:1.5 and 1:3
Overall weight	35 lbs
Largest weight component	17.25 lbs
Available handrims	Plastic coated (1:3), stainless steel (1:1.5)
Battery capacity	7 Ahr
Maximum range	9.5 miles
Available tire size	24"
Available tire type	Pneumatic and pneumatic with inserts
Weight capacity	250 lbs
Warranty	2 years
Time to full charge	2.5 hr to 3.0 hr



Quickie® 2 New Features



Multi-position Fork
Provides greater range of seat height and a tighter turning radius.



Uni-lateral Wheel Lock
Allows rider to easily engage both wheel locks with one hand.



Jay® Precision™ Back
Features adjustable upholstery and height adjustable laterals for ideal positioning support.

