



Invacare® Allegro® Wheelchair



Yes, you can.™



Yes, you can.™

The Invacare® Allegro® chair is the ultimate in versatility. Designed from the ground up as a complete mobility system, the Allegro addresses multiple positioning needs and growth requirements. Its adjustable-depth and adjustable-angle back and modular cross-brace system ensure a perfect fit and adapt to any seating and positioning system. Choose between the folding or rigid model; a version without seat rails for a lower seat-to-floor height; or the reclining back, which is interchangeable with the fold-down back.

## Highlights

- Adjustable-depth frame
- Modular crossbrace system
- Grows from 12" x 12" to 20" x 18"
- Reclining back available
- Rigid model available
- Available without seat rails
- Reinforced caster headtube mounting system
- Several axle plate choices
- Numerous seat-to-floor height positions



Hydraulic cylinders offer smooth recline of up to 170°.



The Allegro without seat rails lowers seat-to-floor height by 1.75" and is ideal for use with seating systems.



The Allegro's unique modular crossbrace system allows a host of configuration possibilities.



The rigid Allegro model provides a stiffer ride and easy integration of seating systems.



The Invacare Transport Ready option consists of four factory-installed tie down brackets that meet today's ANSI/RESNA WC19 vol.1 - section 19.5.3 standard for the frontal impact test.\*



Adjustable-angle stroller handles provide maximum flexibility in adjusting the chair's geometry.

\*As tested with the Invacare® KSS Seating System. Invacare has not performed testing with any other seating system. Invacare recommends that wheelchair users NOT be transported in vehicles of any kind while in wheelchairs. As of this date, the Department of Transportation has not approved any tie-down system for transportation of a user while in a wheelchair in a moving vehicle of any type.

## Specifications

### Dimensions

Overall*	seat width plus 8.5"			
Width	12"-20" in 1" increments			
Depth	12"-18" in 1" increments			
Front seat-to-floor	3"	5"	6"	8"
Fork Type	Caster	Caster	Caster	Caster
3"	15", 16.25"	16", 17.25"	n/a	n/a
3" fork/1" spacer	16", 17.25"	17", 18.25"	n/a	n/a
3" fork/2" spacer	17", 18.25"	18", 19.25"	n/a	n/a
5"	16.5", 17.75"	17.5", 18.75"	18", 19.25"	n/a
5" fork/1" spacer	17.5", 18.75"	18.5", 19.75"	19", 20.25"	n/a
6" upper mounting	16.25", 17"	16.5", 17.25", 18"	17", 17.75", 18.5"	n/a
6" lower mounting	17.5", 18.25"	17.75", 18.5", 19.25"	18.25", 19", 19.75"	n/a
8" 2 holes	n/a	n/a	n/a	18.75", 20", 21.25"

### Rear seat-to-floor\*\*

20"	14.75" - 19.25"
22"	15.75" - 20.25"
24"	16.5" - 21"

### Knee-to-heel\*\*\*

	min.	max.
70° tapered footrest	16.25"	20.25"
70° swingaway footrest	14"	19"
70° MFX footrest	8.25"	14"
Elevating legrest	15.75"	21.75"
90° swingaway footrest	5.75"	10.75"
60° swingaway footrest	15.25"	19"
60° MFX footrest	11.5"	15.5"

### Cane options

Sportster	10"-14" adjustable	n/a fixed
10" with push handles	16"-18" adjustable	n/a fixed
Straight with push handles	12"-16" adjustable	16"-20" fixed

### Arm height†

	Folding frame with seat rails	Rigid frame without seat rails
Adjustable, dual-point	10"-14" conventional	n/a
Fixed, dual-point	10" conventional	n/a
Cantilever	5"-10"	7"-12"
T-Arm	n/a	6"-10" .5 increments 8"-14" 1" increments

Weight 250 lb.

Limited warranty Lifetime on frame and crossbraces

\*Measurement taken from outside of handrim to outside of handrim with 0° camber.

\*\* Rigid frames and folding frames without rails, subtract 1.75" from seat-to-floor height to obtain distance from top of side frame. Range adjusts in 5/8" increments.

\*\*\*Range adjusts in 1" increments. Measurement taken from top of seat rail to back of footrest.

†Measurement taken from top of seat rail to top of armrest. Range adjusts in 1" increments.

Invacare Corporation  
www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

Invacare and Allegro are registered trademarks and Yes, you can. and the blue medallion design are trademarks of Invacare Corporation. Specifications are subject to change without notification.

© 2004

Invacare Corporation

Form No. 01-209  
Rev 05/04



"Invacare says 'Yes, you can'; and I couldn't have said it better myself."