



**Quickie®**

**FOLDING**

**Family**







At Quickie our focus has always been on delivering the best product in terms of quality, performance and reliability. To attain this goal, we are committed to bringing design innovations and the use of quality materials together in a superior product that improves the life of the user.

The Quickie Folding Family offers the superior style, performance and excellence in design that have made the Quickie 2 the leading product in the market. You can always be sure that...

***...Nothing Beats a Quickie***

The Quickie 2 is already the most popular ultra-light folding frame wheelchair available today, and we've found a way to make it even better. With twelve new features and options, the Quickie 2 now offers more fine tune adjustments for improved versatility, customization and ease of use. With the Quickie 2 Family the possibilities are endless.



### Quickie 2 Highlights:

- Fully adjustable axle plate with camber, center of gravity, wheel base width and seat height adjustments
- Even more caster adjustability to accommodate any modifications of the frame that affect alignment
- Angle-adjustable backrest
- Fully modular design for complete customization
- The Quickie True Fit Program allows the chair to grow
- High quality Quickie components and hardware ensure durability and performance
- Widest variety of hanger, footplate, frame, caster, wheel, backrest and armrest options
- New Jay Precision™ Back option
- Travel wheel option
- One-arm drive option (see Quickie Kids picture)
- Uni-lateral (one-hand) wheel lock system option



# Quickie® 2HP

Ultra-Light

**The Quickie 2HP** is the lightest folding chair in the Quickie 2 Family. In addition to all the features of the Quickie 2, the Quickie 2HP also includes a fixed front frame for a more solid ride. The Quickie 2HP is ideal for those who demand high performance, low maintenance and everyday manual mobility.



#### Quickie 2HP Highlights:

- Fixed front frame provides a shorter, lighter, more rigid chair
- Performance wheel and caster options
- Increased maneuverability and propulsion





**The Quickie 2 Kids** is the perfect choice for pediatric riders with changing mobility needs and lifestyles. The Quickie 2 Kids incorporates all of the Quickie 2 features into a smaller package to facilitate environmental accessibility and exploration.

#### Quickie 2 Kids Highlights:

- Kid-sized armrests increase wheel access and positioning options
- Modular frame accommodates multiple pediatric seating solutions
- Shorter footplates and legrests for smaller riders
- The Quickie True Fit Program allows the chair to grow with the rider



# Quickie® 2 HD

## Heavy Duty

The Quickie 2 with the Heavy Duty Package offers all the features of the Quickie 2 with a 350 pound weight capacity. Reinforced features increase the rigidity and performance of the Quickie 2HD making it the perfect choice for those who require a more solid frame and all the flexibility of the Quickie 2.



### Quickie 2HD Highlights:

- Stabilizer bar standard for increased performance and rigidity
- Tri-cross brace improves overall durability
- Reinforced back upholstery





# Quickie® M6

Heavy Duty

Finally, a wider, stronger folding chair. **The Quickie M6** offers seat widths up to 30" and a weight capacity of 650 pounds. The Quickie M6, made of high-strength steel and modeled on the Quickie 2 frame, offers interchangeable frame components that allow the chair to be adjusted for the rider's changing needs.

#### Quickie M6 Highlights:

- 650 pound weight capacity
- Lightweight
- Adjustable back height
- 17" seat-to-floor height
- Stabilizer bar
- Stylish frame available in all Quickie colors



**The Quickie Recliner** has all of the features of a Quickie 2 plus a reclining backrest. The Quickie 2 can be converted to a recliner and back again. The ease of use and effectiveness of this chair make it a big hit with caregivers.

**Quickie Recliner Highlights:**

- Reclines 90°-160°
- Ease of use ensured by the dual trigger release with pin lock on the stabilizer bar
- 10" removable headrest extension



**The Quickie LX** is a high-strength, lightweight chair that offers the benefit of the quality components of the Quickie 2 in a more limited design. The Quickie LX is a price-competitive solution when quality and reliability are still a must!

**Quickie LX Highlights:**

- One-piece side frame with limited adjustability and options
- Angle-adjustable front casters
- Finely-tuned seat positioning and back height adjustability



**The Quickie LXI** features the same one-piece side frame as the Quickie LX, but offers slightly more adjustability through the axle plate. The Quickie LXI— an ultra-lightweight chair with limited adjustability made with quality components to ensure Quickie reliability.



**Quickie LXI Highlights:**

- Axle plate adjustability
- Camber adjustability
- Center of gravity adjustability



# REVOLUTION

Ultra-Light & Maneuverable

**The Quickie Revolution** feels and performs like a rigid frame wheelchair with all the convenience of a folding chair. The unique folding system of the Quickie Revolution allows it to fold up so small and compact that it stows almost anywhere.



**Quickie Revolution Highlights:**

- Innovative V-front folding frame
- Camber tubes available in 0°, 4°, 8° and 12°
- Patented center of gravity adjustment
- Angle-adjustable backrest





# FOLDING FAMILY

Quickie Folding Chairs – A complete folding line to suit your mobility needs.

## Ultra-Light



## Lightweight



## Positioning



## Pediatrics



### Quickie 2 & Quickie 2 Kids & Quickie 2 HD

#### Front Seat Height (For Standard Frames Only)

CASTER FORK MOUNTING HOLE	5" LP POLY 5 1/4" FORK	5" LP POLY 7" FORK	6" PNEU 6" FORK	6" POLY 6" FORK	8" PN 5 1/4" FORK
Top	16 15/16"	18 3/4"	18 5/16"	18 3/16"	18 1/2"
Bottom	17 15/16"	19 3/4"	19 5/16"	19 3/16"	19 1/2"

#### Rear Seat Height (For Standard Frames and Standard Axle Plate Position Only)

REAR AXLE/HOLE POSITION	PNEU	FP	PNEU	LP	FP	HP	PNEU/ KEVLAR	HP CLNCHR
1 AND 2	18 7/8"	19 1/2"	19 3/4"	20"	20 1/4"	20 3/8"	20 1/2"	20 1/4"
2 AND 3	17 3/4"	18 3/8"	18 5/8"	18 7/8"	19 1/8"	19 1/4"	19 3/8"	20 1/8"
3 AND 4	16 5/8"	17 1/4"	17 1/2"	17 3/4"	18"	18 1/8"	18 1/4"	19"
4 AND 5	15 1/2"	16 1/8"	16 3/8"	16 5/8"	16 7/8"	17"	17 1/8"	17 7/8"

#### Overall Width (With Standard Camber and Padded Swing-Away Arm)

FRAME WIDTH	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
SPOKE WHEEL W/HANDRIMS	20 1/2"	21 1/2"	22 1/2"	23 1/2"	24 1/2"	25 1/2"	26 1/2"	27 1/2"	28 1/2"	29 1/2"	30 1/2"

Quickie 2	Quickie 2 HD	Quickie 2 Kids
Weight: 27 lbs		
Weight Capacity: 250 lbs	Weight Capacity: 350 lbs	Weight Capacity: 250 lbs

### Quickie 2HP

#### Front Seat Height (For Standard Frames Only)

CASTER FORK MOUNTING HOLE	3"MICROS 5 1/4" FORK	3"MICROS 6" FORK	3"MICROS 7" FORK	5"LP POLY 5 1/4" FORK	5"LP POLY 6" FORK	5"LP PC 7" FORK
TOP	16 5/16"	17 1/8"	18 1/8"	17 3/16"	18"	19"
BOTTOM	17 5/16"	18 1/8"	19 1/8"	18 3/16"	19"	20"

#### Rear Seat Height (For Standard Frames and Standard Axle Plate Position Only)

REAR AXLE/HOLE POSITION	PNEU	LATEX	HP	PNEU	TEN/B-BALL	LATEX	PNEU
1 AND 2	19 3/4"	20 1/4"	20 3/8"	20 1/2"	20 5/8"	21 1/4"	21 5/8"
2 AND 3	18 5/8"	19 1/8"	19 1/4"	19 3/8"	19 1/2"	20 1/8"	20 1/2"
3 AND 4	17 1/2"	18"	18 1/8"	18 1/4"	18 3/8"	19"	19 3/8"
4 AND 5	16 3/8"	16 7/8"	17"	17 1/8"	17 1/4"	17 1/8"	18 1/4"

#### Overall Width (With Standard Camber and Padded Swing-Away Arm)

FRAME WIDTH	12"	13"	14"	15"	16"	17"	18"	19"	20"
SPOKE/MAG WHEEL W/HANDRIMS	21 1/2"	22 1/2"	23 1/2"	24 1/2"	25 1/2"	26 1/2"	27 1/2"	28 1/2"	29 1/2"

Quickie 2 HP
Weight: 22 lbs
Weight Capacity: 250 lbs

### Quickie M6

FRAME WIDTH:	24"	26"	28"	30"
OVERALL WIDTH:	32.5"	34.5"	36.5"	38.5"
FOLDED WIDTH:	14"			

### Quickie M6

Weight:	53 lbs. w/legrests (26 x 20)
Weight Capacity:	650 lbs

### Quickie Recliner

#### Front Seat Height

CASTER FORK MOUNTING HOLE	5" LP POLY 6" FORK	6" PNEU 6" FORK	6" POLY 6" FORK	8" PNEU 5 1/4" FORK	8" POLY 5 1/4" FORK
TOP	17 3/4"	18 5/16"	18 3/16"	N/A	N/A
BOTTOM	18"	19 5/16"	19 3/16"	19 1/2"	19 5/16"

#### Rear Seat Height

REAR AXLE HOLE POSITION	22" FP Poly	24" FP Poly	24" HP Clincher	24" Pneu/ Kevlar
Top	15"	15 3/4"	15 7/8"	16"
Middle	16"	16 3/4"	16 7/8"	17"
Bottom	17"	17 3/4"	17 7/8"	18"

#### Armrest Height

##### Height-Adjustable

9" - 14"

#### Footrest Height

SWING-AWAY	ELR
14" - 21"	16 1/2" - 19"

#### Overall Width (Without Armrests)

Frame Width	14"	15"	16"	17"	18"	19"	20"	Folded
SPOKE/MAG WHEEL W/HANDRIMS & QUICK-RELEASE AXLES	22 1/4"	23 1/4"	24 1/4"	25 1/4"	26 1/4"	27 1/4"	28 1/4"	12 1/2"
SPOKE/MAG WHEEL W/HANDRIMS & THREADED AXLES	22"	23"	24"	25"	26"	27"	28"	12 1/2"

### Quickie Recliner

Weight:	42 lbs. (w/ swing-away footrests, height-adjustable armrests, headrest extension and mag. wheels)
Weight Capacity:	250 lbs

### Quickie LX & LXI

#### Front Seat Height

CASTER FORK MOUNTING HOLE	5"LP Poly 4 3/4" Fork	5"LP Poly 5 1/4" Fork	6"Poly 6" Fork	6"Pneu 6" Fork	8"Poly 5 1/4" Fork	8"Pneu 5 1/4" Fork	8"x2" Caster 8"x2" Fork
TOP	N/A	N/A	17"	17"	N/A	N/A	N/A
BOTTOM	16"	16 1/2"	18"	18"	18"	18"	N/A

#### Rear Seat Height

REAR AXLE HOLE POSITION (LX ONLY)	20" LP POLY	20" PNEU	22" LP POLY	22" PNEU	24" LP POLY	24" PNEU
TOP	15"	15 1/2"	16"	16 1/2"	17"	17 1/2"
BOTTOM	16"	16 1/2"	17"	17 1/2"	18"	18 1/2"

REAR AXLE HOLE POSITION (ADJ.AXLE)	20" LP POLY	20" PNEU	22" LP POLY	22" PNEU	24" LP POLY	24" PNEU
TOP	N/A	15"	15 1/2"	16"	16"	17"
MIDDLE	15 1/2"	16"	16 1/2"	17 1/2"	17 1/2"	18"
BOTTOM	17"	17 1/2"	18"	18 1/2"	18 1/2"	19 1/2"

#### Armrest Height

#### Footrest Height

FLIP-BACK HEIGHT ADJ.	FLIP-BACK/ PADDED S/A HEIGHT ADJ.	HT. ADJ.	SWING-AWAY	ELR/ALR
9 1/4"	8 3/4"-12 5/8"	8 3/4"	7 1/4"-12 1/4"	12 3/4"-15 1/4"
			8 3/4"-13 3/4"	14"-16 1/2"
				15 1/4"-17 3/4"
				16 1/2"-19"
				17 3/4"-20 1/4"

#### Overall Width (Without Armrests)

FRAME WIDTH	12"	14"	16"	18"	20"	FOLDED
MAG WHEEL W/HANDRIMS	19"	21"	23"	25"	27"	12 1/2"
SPOKEWHEEL W/HANDRIMS	20"	22"	24"	26"	28"	12 1/2"

### Quickie Recliner

Weight:	28 lbs
Weight Capacity:	250 lbs

# SPECIFICATIONS

<b>Quickie Revolution</b>							
<b>Front Seat Height (For Standard Frames Only)</b>							
CASTER FORK	5" MICRO	5" POLY		6" PNEU		6" POLY	
MOUNTING HOLE	5 1/4"	5 1/4"		5 1/4"		5 1/4"	
TOP	N/A	17 1/4"		17 1/8"		17 3/4"	
BOTTOM	17 1/4"	18 1/4"		18 7/8"		18 3/4"	
<b>REAR SEAT HEIGHT (FOR STANDARD FRAMES AND AXLE POSITION ONLY)</b>							
CAMBER/REAR	22"	24"	24"	24" HT	26"	26"	26" HT
AXLE HOLE	PNEU	PNEU	LATEX	TEN/BBALL	HP	PNEU	TEN/BBALL
POSITION				TURBO TRD	CLINCHER		TURBO TRD
0°, 4°/NO SPACER	15 1/4"	16 1/8"	15 3/4"	16 1/4"	16 7/8"	17 1/4"	17 5/8"
8°, 12°/NO SPACER	14 1/2"	15 3/8"	15"	15 1/2"	16 1/8"	16 1/2"	16 7/8"
8°, 4° + 3/4°	16"	16 7/8"	16 1/2"	17"	17 5/8"	18"	18 3/8"
8°, 12° + 3/4°	15 1/4"	16 1/8"	15 3/4"	16 1/4"	16 7/8"	17 1/4"	17 5/8"
0°, 4° + 1 1/2°	16 3/4"	17 5/8"	17 1/4"	17 3/4"	18 3/8"	18 3/4"	19 1/8"
8°, 12° + 1 1/2°	15"	16 7/8"	16 1/2"	17"	17 5/8"	18"	18 3/8"
<b>Overall Width (With 4° of camber)</b>							
Frame Width	14"	15"	16"	17"	18"	19"	20"
SpokeWheel*							
with Handrims	24"	25"	26"	27"	28"	29"	30"
<b>Quickie M6</b>							
Weight:	23 lbs.						
Weight Capacity:	250 lbs						

## Quickie 2 True Fit Program

Physical changes can cause a chair not to fit properly. To ensure proper fit, Sunrise offers one-time replacement for growth of the seat sling, back upholstery, footplates, hangers, back tubes, frame and crossbrace at no extra charge. (For up to three years from delivery on all Quickie 2 family wheelchairs shipped after November 15, 1995.)

ACTUAL MEASUREMENTS MAY VARY ACCORDING TO CHAIR CONFIGURATION AND ACCESSORIES. PLEASE CONSULT ORDER FORM FOR COMBINATIONS AND RESTRICTIONS. SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE.

LIFETIME FRAME AND AXLE WARRANTY

We highly recommend distribution through a certified Rehab Professional when fitting and ordering any Quickie Products.

