



Yes, you can.®

# Absolute™ Series

REHAB  
POWER  
CUSTOM MADE



*Comfort comes standard™*

# Invacare® Absolute™ Series Seating Products



Yes, you can.®

**Invacare® Absolute™ Series: Simple seating solutions for comfort, fit and function.**

For clients with mild to moderate seating needs, the Invacare® Absolute™ Series Cushions offer a wide range of simple but highly functional solutions. Absolute™ and Single Density cushions are designed to provide stability and comfort at an affordable price. The Invacare® Comfort-Mate Extra cushion goes one step further with a comfortable, molded design and a durable moisture-proof coating.

## Accessories and Options



### Invacare® Absolute™ Cushion


Model no. EC

- Molded polyurethane foam provides comfort and durability
- Waterproof cover protects foam from moisture
- Non-slip bottom with hook and loop fasteners keeps cushion securely in place

Cushion Thickness

A	Height at adductor	3.25"
B	Height at abductor	2.25"
C	Height at rear	2"

Cushion weight approximately 2 lb.



### Invacare® Single Density Cushion


Model no. SD

- Contoured deformation cut foam provides comfort and stability
- Available in over 100 custom sizes
- Waterproof cover with non-slip bottom

Cushion Thickness

A	Height at rear (adductor)	2"
B	Height at rear (center)	1.5"
C	Height at front	2.5"

Cushion weight approximately 2 lb.



### Invacare® Comfort-Mate Extra Cushion


Model no. CMEX

- Comfortable, molded dual-firmness polyurethane foam promotes good postural position
- Moisture-proof coating seals and protects foam from moisture; easily wipes clean
- Waterproof cover removes for laundering
- Beveled sides allow cushion to be used easily with Infinity™ Drop Base

Cushion Thickness

A	Height at adductor	2.625"
B	Height at abductor	3.0"
C	Height at rear	1.75"
D	Height at leg trough	2.25"

Cushion weight approximately 2.3 lb.



**Invacare Corporation**  
www.invacare.com

#### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

#### Canada

570 Matheson Blvd. E.  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

Invacare, the Blue Medallion Design and Yes, you can. are registered trademarks and Absolute, Infinity and Comfort comes standard are trademarks of Invacare Corporation. Specifications are subject to change without notification.

© 2005

Invacare Corporation

Form No. 05-142



**"Invacare says 'Yes, you can.®,' and I couldn't have said it better myself."**