

QUANTUM[®]
REHAB

INNOVATIVE REHAB SOLUTIONS

A Division of Pride Mobility Products[®] Corp.



SYNERGY[®] SHAPE

The new Synergy[®] SHAPE back offers superior comfort and pressure management utilizing dual layer foam construction with Sensus foam layer. It's unique one-hand, one-motion release feature makes removal extremely convenient. The contoured design assists with maintaining proper positioning.

FEATURES...

- 4 inches of depth adjustment
- Recline angle adjustment 57 to -30 degrees
- Lightweight
- High-density, high-resiliency base for support and to prevent bottom out
- Lateral supports available
- Tilt rated to 300 lbs.