



## SYNERGY SHAPE

### Back Cushion

The new Synergy Shape back cushion offers clients superior back support with pressure management. The SHAPE's dual layer foam construction features Sensus foam and its contoured design assists with proper positioning. Its unique one-hand, one-motion quick release feature makes removal extremely convenient.



**MSRP: \$799**

**14"×19" & 22"×19" - SEBX10005**


**ADP Price: \$640**

**16"×19" , 18"×19" & 20"×19" - SEBX10010**

**ADP Price: \$600**

## SOLUTION™ Seating Cushion

The Solution seating cushion has a contoured, wrap-and-lift base designed with high density, high resiliency foam that redistributes weight away from areas vulnerable to skin breakdown. A twin-cell gel insert and dual foam layers combine with a fluid resistant inner liner and 4-way stretch cover to provide increased pressure management and positioning along with durability and comfort.



**MSRP: \$399-\$449**


**ADP Price: \$319-\$359**

**SESX10015 - Standard Size**

**SESX10020 - Custom Sizes\***

## SPECTRUM Seating Cushion

The Spectrum seating cushion takes our wrap-and-lift design to the next level by providing a choice of deep twin-cell gel or Sensus form-retaining foam insert. The choice of insert allows the Spectrum to offer more personalized pressure management and positioning to match the individual client's needs while providing enhanced comfort.



**MSRP: Visco/Gel \$549-\$599**

**ADP Price: \$439-\$479**

**Choice of Inserts:**

**SESX10025 - Standard Size Visco**


**SESX10030 - Custom Sizes Visco\***

**SESX10035 - Standard Size Gel**

**SESX10040 - Custom Sizes Gel\***

## Structure Seating Cushion

The Structure seating cushion expands even further on our wrap-and-lift design by combining a choice of twin-cell gel or Sensus form-retaining foam insert with modular positioning components. This allows for client-specific customization for maximum pressure management, unparalleled positioning and incredible comfort.



**MSRP: Visco/Gel \$749-\$799**

**ADP Price: \$599-\$639**

**Choice of Inserts:**

**SESX10060**

**Standard Size Visco/Gel**

**SESX10065**

**Custom Sizes Visco\*/Gel\***