

Concept 45





1. Adjustable angle back from 90 to 110 degrees
2. Seat to floor heights from 14.5" to 18.75"



New back posts to ensure client comfort. Allows for better headrest adjustability.



Extended arms to allow for better tray support and to assist in transferring.



1. Standard padded headrest.
2. Optional comfort headrest with adjustable Multi Axis mount.

The **Concept 45** is one of Invacare Canada's most affordable compact tilt-n-space chairs available. With 45° of tilt on a compact frame, the Concept 45 provides the user with excellent support and greater flexibility for overall functioning.

Standard Features:

- Full Adjustable Height Arms
- Adjustable Angle Back
- Permanent Rear Axle
- 22" Urethane Rear Wheels
- Push to Lock Wheel Locks
- ABS Seat Pan
- Chrome Handrims
- 7" Urethane Casters
- Black Metallic Frame
- Heel Loops
- Anti Tippers
- Standard Padded Headrest
- Lower Frame Impact Guard
- 250 lb. Weight Limit

No Charge Options:

- Choice of Composite or Aluminum Footplates
- Choice of Chrome or Aluminum Handrims
- Choice of 8 Seat to Floor Heights
- Pull to Lock Brakes

Upcharge Options:

- 19"/20" Seat Widths
- 19" Seat Depth
- 22"/26" Fixed Back Heights
- 20"/24" Recliner Back Heights
- Comfort Headrest with Multi Axis Mount
- Extra Large Footplates
- Adjustable Angle Footplate
- Choice of 6 Front Rigging Options
- Personal Back

Specifications:

- Overall Width: 23 1/4"
- Overall Depth: 44 3/4"
- Overall Height: 40"
- Weight: 55 lbs. (16"x16" chair) w/front rigging. w/o seating
- Seat Width: 15" - 20"
- Seat Depth: 15" - 19"
- Seat-to-Floor: 14.5" - 18.75"
- Seat-to-Back Height: 18", 20"
- Tilt: 45 degrees
- Recline 30 degree

Note: Above specs are measured on a 16x16 chair.

Frame Colour



Metallic Black



Invacare Canada
www.invacare.ca

Canada
570 Matheson Blvd. E.
Unit #8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Invacare is a registered trademark and Yes, you can. is a trademark of Invacare Corporation. Specifications are subject to change without notification.

Made in Canada.



"Invacare says 'Yes, you can,' and I couldn't have said it better myself."

© 2002
Invacare Corporation
Form No. 02-252C
Rev. September 2004