



The New Generation

Transport Chair – TC-100

The AMS TC-100 folding transport chair features loop-style hand brakes that lock down on large 12" rear tires. The TC-100 includes removable padded armrests, swing-away footrests and seat belt.



Built with the highest standards in quality and reliability

Since 1989 AMS HOMECARE INC. has built an enviable reputation by delivering high quality, reliable and comfortable self-care equipment addressing the needs of thousands of customers seeking to maintain their freedom and independence. We are dedicated to providing innovative, reasonably priced solutions specially developed with discriminating consumers in mind.

The Shoprider® brand represented exclusively in Canada by AMS HOMECARE INC. is recognized worldwide for superior design, technical sophistication as well as functionality in mobility products. Through warehouse-service centers conveniently located across the country, AMS HOMECARE INC. is the leader in the distribution of Shoprider® scooters and accessories.

In addition, we will be offering an expanding product line of AMS-branded products and accessories such as our current scooter canopy & cane holders and bed rails.

AMS HOMECARE INC. is a growing company committed to bringing to market the best products, devices and solutions to satisfy mobility and homecare needs.



Specifications Transport Chair (model TC-100):

Seat width	inch	19
Seat depth	inch	16
Back height	inch	18
Weight capacity	kg (lb)	136 (300)
Unit weight	kg (lb)	14 (30)
Wheel sizes	Front (in)	8
	Rear (in)	12
Removable padded armrests	Yes/No	Yes
Swing-away footrests	Yes/No	Yes
Handle with loop hand brake	Yes/No	Yes
Rear tire locks	Yes/No	Yes
Steel chrome plated finish	Yes/No	Yes
Fold-down back	Yes/No	Yes

All specification subject to change without prior notice.

Contact your AMS Homecare Dealer today:

Features:

- Loop style hand brakes
- 12" rear wheels
- Removable padded armrests

