



Shown: Manual Tote

MANUAL WHEELCHAIR CARRIERS

CADDY

The easy-to-install Silver Star™ Caddy wheelchair carrier from Pride® makes transporting your manual wheelchair simple and inexpensive. Plus, like all Silver Star products, the Caddy features an all-steel frame and quality construction.

- Removable sections for easy storage
- Mounts to Class I hitch
- 100 lbs. carrying capacity

MANUAL TOTE

The Silver Star™ Manual Tote wheelchair carrier from Pride® is an easy, inexpensive way to transport your manual wheelchair. Plus, like all Silver Star products, the Manual Tote features an all-steel frame and quality construction.

- Spring loaded platform for easy lifting
- Mounts to Class I hitch
- 100 lbs. carrying capacity

HITCHHIKER

Lightweight, simple to install, and easy to use, the Silver Star™ Hitchhiker manual wheelchair carrier from Pride® is the easy way to transport your manual wheelchair. Just push your wheelchair into position and the Hitchhiker locks it in place. Plus, like all Silver Star products, the Hitchhiker features an all-steel frame and quality construction.

SilverStar™ Lifts ★

A Division of **Pride**® Mobility Products Corp.

MANUAL WHEELCHAIR CARRIERS



Caddy



Manual Tote



Hitchhiker

The Caddy, Manual Tote and Hitchhiker can easily lift the Pride® Stylus™ and folding wheelchairs.

CADDY

Base weight lift	13 lbs. ¹
Maximum manual chair weight	100 lbs.
Maximum product width	N/A
Maximum travel distance from the ground	N/A
Hitch class	Class I ²
Vehicle platform	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover
Warranty	3-year transferable

MANUAL TOTE

Base weight lift	40 lbs. ¹
Maximum manual chair weight	100 lbs.
Maximum product width	N/A
Maximum travel distance from the ground	19"
Hitch class	Class I ²
Vehicle platform	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover
Warranty	3-year transferable

HITCHHIKER

Base weight lift	20 lbs. ¹
Maximum manual chair weight	100 lbs.
Maximum product width	14"
Maximum travel distance from the ground	19"
Hitch class	Class I ²
Vehicle platform	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover
Warranty	3-year transferable

¹ Weight does not include any accessories

² If vehicle is equipped with a Class II or III hitch, it requires a hitch reducer; additional charges apply

SilverStar™ Lifts 

A Division of **Pride** Mobility Products Corp.

182 Susquehanna Ave. • Exeter, PA 18643 • (800) 800-8586 • www.silverstarmobility.com • www.pridemobility.com

SPECIFICATIONS