

The new Invacare® Lynx™ L-3 and L-4 micro-portable scooters allow you to accomplish daily activities with ease! You'll be amazed at how simple they are to assemble and disassemble for transport. Their compact size and seamless functionality allow them to easily transition from the tight turns of the home to the wide-open spaces of the outdoors. With many of the same attributes on both models, you can choose between two stylish scooters packed with great features at a competitive price.

Lynx™ L-3

The three-wheeled Lynx L-3 scooter takes portability to a new level. The L-3 model conveniently disassembles into four compact pieces for easy transport. Pick the destination, stow it in the trunk and go... the L-3 model will do the rest. With an estimated seven miles of range, a plug-in battery charger and flat-free tires, consumers can plan their schedule around what they want to do, not around the limitations of their scooter.

Features

- Easy 4-piece disassembly
- Estimated 7 mile travel range
- 0 - 5 mph variable speed control
- 250 lb. weight capacity
- Connector-free, drop-in battery pack
- On or off-board battery charging
- Foam-filled, flat-free tires
- Comfortable 360° rotating chair
- Ample leg room and foot clearance
- 3-position adjustable tiller

Shroud Colors



Blue



Red



360° Rotating Swivel Chair



On/Off Chair Battery Charging



Lynx™ L-4



The Lynx L-4 offers the features of the L-3 model with the performance and bold look of a four-wheel scooter. The stable four-wheel base of the L-4 model handles varying surfaces with ease. Customize the L-4 model to fit the consumer by simply adjusting the seat height, armrest width and tiller angle. Get comfortable, get set and go! With a top speed of five mph, you can feel confident the Lynx L-4 scooter will get you to your destination on time.

Features

- Easy 4-piece disassembly
- Estimated 7 mile travel range
- 0 - 5 mph variable speed control
- 250 lb. weight capacity
- Connector-free, drop-in battery pack
- On or off-board battery charging
- Foam-filled, flat-free tires
- Comfortable 360° rotating chair
- Ample leg room and foot clearance
- 3-position adjustable tiller

Shroud Colors



Blue

Red

Intuitive User Controls



Flat-free Tires



Specifications

Overall dimensions	L-3	L-4
Base length	39.8"	39.8"
Base width	26"	26"
Step height (floor pan to ground)	4"	4"
Seat height (floor pan to seat cushion)	14"-16" padded seat	14"-16" padded seat
Seat dimensions		
Width	16.5"	16.5"
Depth	14.5"	14.5"
Back height	16"	16"
Arm dimensions		
Length	11.6"	11.6"
Height between seat and arm	9"	9"
Width between arms	15"	15"
Speed range		
	0-5 mph	0-5 mph
Estimated travel range	7 miles	7 miles
Weight limit		
	250 lb.	250 lb.
Maximum incline capability		
	8° grade	8° grade
Tires		
Front	8" x 2"	8" x 2"
Rear	8" x 2"	8" x 2"
Turning radius		
	39.4"	51.2"
Weight		
Front base section	30 lb.	39 lb.
Rear base section (with shrouds)	25 lb.	25 lb.
Entire seat	24 lb.	24 lb.
Batteries	23 lb.	23 lb.
Basket	1 lb.	1 lb.
Total weight		
	103 lb.	112 lb.
Limited warranty - frame, electronics, transaxle		
	1 year	1 year

Four Easy Steps for Portability



1. *Remove Seat*



2. *Remove Battery*



3. *Fold Tiller Down*



4. *Remove Back Wheels*



Disassembles for Easy Transport

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Invacare, the Blue Medallion Design and Yes, you can. are registered trademarks and Lynx is a trademark of Invacare Corporation. Specifications are subject to change without notification.

© 2006
Invacare Corporation
Form No. 06-047



"Invacare says 'Yes, you can,' and I couldn't have said it better myself."



Yes, you can.™

Lynx™ Scooters



Anything else is a compromise.™