



Yes, you can.®

# Meteor Scooter





Yes, you can.®



Model no. S8465-ARIN - Red  
Model no. S8465-SBIN - Blue

The Invacare® Meteor is a comfortable, safe ride providing complete freedom to enjoy total independence. This robust scooter incorporates a powerful motor designed to carry a maximum user weight of 440 lbs. A convenient driving position can be obtained by adjusting seat height and horizontal seat position. The armrests are width and height adjustable, with a reclining backrest that can be angled to your requirements. Finally, the swivel seat and swing-up armrests allow for easy transfer when getting on and off of the Meteor. The Meteor is available in a 4-wheel version only and emphasizes safety, user friendliness and unique driving capabilities.

## Features

- Handle Bar - an ergonomic handle bar enables light and easy turning
- Reclining Backrest - adjustment for individual comfort
- Swivel Seat - enables easy access getting in and out of the seat
- Adjustability - optimal driving position can be obtained by adjusting seat height, horizontal seat position and the width and height of the armrests
- Light System - stylish lighting system
- Hand Brake - use hand brake as both park brake and traditional hand brake
- Throttle - can be operated by a light "pull action"
- Safety Belt - provides support while driving
- Fully programmable for individual client needs

**Invacare Canada**  
www.invacare.ca

**Canada**  
570 Matheson Blvd. E.  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324



### TILLER

Independent of the Tiller angle the scooter headlight always shines directly ahead.



### CONTROL PANEL

The spacious control panel includes:

- "Easy touch" buttons
- Illuminating indicators when button is activated
- Accurate easy to read battery indicator
- Handy "Speed dial"



### SEAT POSITIONING

7" sliding rails allows adjustment of the seat position.



### ADJUSTABLE ARMRESTS


The "flip-up" armrests are width, angle and height adjustable.

## Specifications

Width	25 1/2"
Height	48"
Length (w/o bumpers and anti-tippers)	57 1/2"
Seat width	17-21 1/2"
Seat depth	17 1/2"
Seat height (measured from the chassis, manually adjustable)	18-20"
Backrest height (w/o headrest)	17"
Backrest height (with headrest)	25"
Backrest angle (manually adjustable)	0°, 10°, 15°, 30°
Armrest height	9-11"
Weight, including batteries (lbs.)	292
Maximum weight capacity (lbs.)	440
Speed	12 km/h
Batteries	2 x 70 Ah
Tire pressure	40 psi
Minimum turning radius	63"
Safe climbing capability 265 lb load	N/A
Safe climbing capability up to 330 lb load	12°
Maximum obstacle height	4 1/2"

© 2006  
Invacare Canada

Form No. 05-451C



*"Invacare says 'Yes, you can'.  
and I couldn't have said it  
better myself."*