

GO-GO[®] TRAVEL LIFT

The new Go-Go[®] Travel Lift is the ideal enhancement for people on the go who want to take along their travel vehicle or small scooter. Featuring a unique rigid steel skeleton and color enhanced body, and available at a low cost, the Go-Go Travel Lift offers function, style and value.

- Perfectly complements the Go-Go[®] Ultra Travel Scooter and Go-Chair
- Adjustable base legs, boom length and height for easy vehicle installation
- Standard corded remote for easy operation
- 130 lbs. lifting capacity
- End cover protects automobile while hoisting
- Upper removable boom arm weighs less than 14 lbs.
- Base cover enhances appearance while protecting cargo
- Boom end cover helps protect inside of vehicle while hoisting



Go-Go Travel Lift shown with optional Trailer Hitch Mount available March 2005.

800-800-8586 (U.S.) • 888-570-1113 (Canada) • www.pridemobility.com

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

GOGO TRAVEL LIFT 11-04

Pride[®]
Mobility Products Corp.