



Pride Mobility Products Corp.	BACKPACKER INTERIOR LIFT⁴	DOUBLE ENTRY EXTERIOR SCOOTER PLATFORM	FULL EXTERIOR PLATFORM	POWER CHAIR EXTERIOR PLATFORM WITH DOCKING STATION¹	QUAD XL SCOOTER EXTERIOR PLATFORM WITH EXTENDED RAILS	POWER CHAIR EXTERIOR PLATFORM WITH STRAPS	SINGLE ENTRY EXTERIOR SCOOTER PLATFORM	POWER TOTE MINI EXTERIOR LIFT	CADDY WHEELCHAIR EXTERIOR CARRIER	HITCHHIKER WHEELCHAIR EXTERIOR CARRIER	MANUAL TOTE WHEELCHAIR EXTERIOR CARRIER	POWER TOTE WHEELCHAIR EXTERIOR LIFT
BASE WEIGHT	106 lbs.	106 lbs.	107 lbs.	108 lbs.	116 lbs.	105 lbs.	112 lbs.	47 lbs.	13 lbs.	20 lbs.	40 lbs.	47 lbs.
MAXIMUM SCOOTER/ POWER CHAIR OVERALL WEIGHT	325 lbs.	325 lbs.	325 lbs.	325 lbs.	325 lbs.	325 lbs.	325 lbs.	150 lbs.	100 lbs.	100 lbs.	100 lbs.	100 lbs.
MAXIMUM PRODUCT WIDTH	26"	27"	27"	27"	27"	27"	27"	22"	N/A	14"	N/A	14" Folded
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	26"	20"	20"	20"	20"	20"	20"	20"	N/A	19"	19"	20"
MOTOR	Sealed 12 Volt DC	Sealed 12 Volt DC	Sealed 12 Volt DC	Sealed 12 Volt DC	Sealed 12 Volt DC	Sealed 12 Volt DC	Sealed 12 Volt DC	Sealed 12 Volt DC	Must manually lift 50% of the weight of the wheelchair to store	Uses handles to tilt and store platform	Uses a spring loaded platform for lifting	Sealed 12 Volt DC
WARRANTY	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable
HITCH CLASS	N/A	Class III Class II (190 lbs. max unit weight)	Class III Class II (190 lbs. max unit weight)	Class III Class II (190 lbs. max unit weight)	Class III Class II (190 lbs. max unit weight)	Class III Class II (190 lbs. max unit weight)	Class III Class II (190 lbs. max unit weight)	Class II or Class III	Class I ²	Class I ²	Class I ²	Class I ²
OPTIONS	Battery pack w/charger	License plate kit, Battery pack w/charger, Swing-away mount	License plate kit, Battery pack w/charger, Swing-away mount	License plate kit, Battery pack w/charger, Swing-away mount	License plate kit, Battery pack w/charger, Swing-away mount	License plate kit, Battery pack w/charger, Swing-away mount	License plate kit, Battery pack w/charger, Swing-away mount	License plate kit, Battery pack w/charger, Swing-away mount ³	N/A	N/A	N/A	License plate kit, Swing-away mount ³
VEHICLE PLATFORM	Mini Van	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover
LIFTS PRIDE® PRODUCTS	All Pride® Scooters, except: Hurricane®, Legend® XL, Wrangler® All Jazzy® PCs with contoured, Specialty & Synergy® seating except: 1170, 1420, 1650 All Jet™ PCs with contoured, Specialty & Synergy® seating All Quantum® PCs with contoured, Specialty & Synergy® seating, except: 1170, 1420, 1650 TRU-Balance® seating on the 1113ATS, 1103 Ultra, Q1107, Q600, Q610	Rally® 3 & 4, Sundancer® 3, Legend® 3 & 4, Celebrity® X 3 & 4, Victory® 3 & 4, Maxima 3 & 4	All Jazzy® PCs with contoured, Specialty & Synergy® seating except: 1170, 1420, 1650 All Jet™ PCs with contoured, Specialty & Synergy® seating All Quantum® PCs with contoured, Specialty & Synergy® seating, except: 1170, 1420, 1650 TRU-Balance® seating on the 1113ATS, 1103 Ultra, Q1107, Q600, Q610	All Jazzy® PCs with contoured, Specialty & Synergy® seating, except: 1170, 1420, 1650, 600, 610 All Quantum® PCs with contoured, Specialty & Synergy® seating, except: 1170, 1420, 1650, 600, 610 TRU-Balance® seating on the 1113ATS, 1103 Ultra, Q1107	Legend® XL & Hurricane®	All Jazzy® PCs with contoured Specialty & Synergy® seating, except: 1170, 1420, 1650, 600, 610 All Quantum® PCs with contoured Specialty & Synergy® seating, except: 1170, 1420, 1650, 600, 610 TRU-Balance® seating on the 1113ATS, 1103 Ultra, Q1107	All Pride® Scooters (3-wheel only)	Sonic®, 3-wheel Go-Go® Ultra Travel Scooter & Go-Go Elite Traveller™	Pride® Stylus™ & Folding Wheelchairs	Pride® Stylus™ & Folding Wheelchairs	Pride® Stylus™ & Folding Wheelchairs	Pride® Stylus™ & Folding Wheelchairs

¹ Requires power chair modification; call for assistance

² If vehicle is equipped with a Class II or Class III hitch, it requires a hitch reducer; additional charges apply

³ If using a Class III Hitch

⁴ Side entry also available, max product width 25.5"; max product length 32"