




### Revo<sup>®</sup> 3-wheel

### Revo<sup>®</sup> 4-wheel

<b>Model number:</b>	SC63	SC64
<b>Weight capacity:</b>	300 lbs.	300 lbs.
<b>Maximum speed:</b>	Up to 4.5 mph	Up to 4.5 mph
<b>Ground clearance:</b>	2.25"	2.25"
<b>Turning radius:</b>	38"	50"
<b>Overall length:</b>	43"	44.75"
<b>Overall width:</b>	21.5"	21.5"
<b>Wheels</b>		
<b>Front:</b>	2" x 8" solid	2" x 8" solid
<b>Rear:</b>	3" x 9" solid	3" x 9" solid
<b>Anti-tips:</b>	Rear	Rear
<b>Suspension:</b>	N/A	N/A
<b>Transport weight</b>		
<b>Total w/o batteries:</b>	106.5 lbs.	125 lbs.
<b>Heaviest piece when disassembled:</b>	48 lbs. (rear section)	51.5 lbs. (front section)
<b>Drive train:</b>	Sealed transaxle, 24 volt DC motor	
<b>Dual braking system:</b>	Electronic, regenerative and electromechanical (electronic, regenerative disc brakes)	
<b>Power elevating seat:</b>	N/A	
<b>Std. Seating:</b>	Compact foldable molded plastic	
<b>Battery charger:</b>	Off-board, 5 amp (standard)	
<b>Per-charge range:</b>	Up to 20 miles	
<b>Battery requirements:</b>	(2) 12V deep cycle (AGM or gel cell type recommended)	
<b>Battery size:</b>	U-1	
<b>Warranty:</b>	3-year limited warranty	

#### Colors:



Red

Blue

#### Benefits

- Innovative interlocking drop-in battery packs for easy removal/installation
- Easily disassembles into 5 lightweight, easy-to-manage pieces
- 300 lbs. weight capacity
- Controls at your fingertips
- Easy-to-use external freewheel lever
- State-of-the-art, microprocessor-based motor controller
- Seat offers height adjustment, adjustable 360° seat rotation, and flip-up armrests for optimal user comfort
- Adjustable tiller for maximum comfort
- Headlight and taillight for added safety
- On-board accessory bracket
- Easily accessible charger port and circuit breaker

#### Options

- Forearm crutch holder
- Lap belt
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Heavy-duty dust cover
- Saddlebag
- Quad cane holder
- Cane, clip attachment



- Backpacker
- Double Entry
- Full Platform



Double entry



Backpacker

These lifts are compatible with the Revo<sup>®</sup>.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

800-800-8586 (US) • 888-570-1113 (Canada) • [www.pridemobility.com](http://www.pridemobility.com)

**Pride<sup>®</sup>**  
Mobility Products Corp.

# REVO®

THE REVO PROVIDES SLEEK STYLING  
AND SUPER PORTABILITY.

## REVO®

The Revo® delivers “full-sized” scooter style, performance, and features in a compact, lightweight design that easily disassembles for transport while remaining highly affordable.



**Pride**<sup>®</sup>  
Mobility Products Corp.

800-800-8586 (US) · 888-570-1113 (Canada) · [www.pridemobility.com](http://www.pridemobility.com)

REVOSS06-06