

Footnotes

1. Solid tires available on most models.
2. Varies with user weight, terrain type, battery charge, battery condition and tire pressure.
3. Battery weight not included.
4. This product is coded by the SADMERC as A9270 (non-covered item). Specifications are preliminary.
5. Battery pack total weight is 27 lbs. each.
6. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
7. Battery pack is 20 lbs.
8. Durable solid plastic
9. Deep cycle (AGM or gel cell type recommended)
10. 12 AH battery pack weighs 21 lbs.; 17 AH battery pack weighs 33 lbs.
11. Standard with power elevating seat
12. Up to 13 miles with optional 17 AH batteries
13. All shroud panels included

Pride® Scooters

"END-TO-END" SOLUTIONS



	Go-Go® ULTRA X ⁴ TRAVEL SCOOTER	Go-Go ELITE TRAVELLER™	Revo®	Sonic®	Rally®	Sundancer® ¹¹	Legend®	Celebrity® X	Victory®	Maxima	Legend® XL	Wrangler® PMV	Hurricane® PMV
Front Tires ¹	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2" x 8" SOLID	2" x 7.5" SOLID	2" x 8"	2.5" x 9" SOLID	3.5" x 10"	3.5" x 10"	3.5" x 10"	3" X 10" SOLID	4" x 13"	6" x 13"	4" x 13"
Rear Tires ¹	2.5" x 8" SOLID	2.5" x 8" SOLID	3" x 9" SOLID	2.5" x 7.5" SOLID	2.5" x 9"	2.5" x 9" SOLID	4" x 10"	4" x 10"	3.5" x 10"	4" X 10" SOLID	4" x 13"	6" x 13"	4" x 13"
Maximum Speed ²	UP TO 4.25 MPH	UP TO 4.25 MPH	UP TO 4.5 MPH	UP TO 4.25 MPH	UP TO 4.25 MPH	UP TO 4.5 MPH	UP TO 5.75 MPH	UP TO 6 MPH	UP TO 5 MPH	UP TO 5.25 MPH	UP TO 5.8 MPH	UP TO 5/10 MPH	UP TO 4.5/9 MPH
Ground Clearance ⁶	2.625" (3W) 2.5" (4W)	2.625" (3W) 2.5" (4W)	2.25"	1.75"	3"	3.5"	3.5"	3.5"	4"	4.5"	4"	7"	4"
Turning Radius (3-wheel) ⁶	32.5"	32.5"	38"	34"	39"	38"	41"	41.5"	41"	52"	N/A	N/A	N/A
Turning Radius (4-wheel) ⁶	43.25"	43.25"	50"	N/A	54"	N/A	57"	59"	58"	58"	65"	77"	65"
(3-wheel) Length/Width ⁶	36.5"/19.25"	36.5"/19.25"	43"/21.5"	37"/21"	41"/21"	41.625"/23.25"	45.375"/24"	46"/24"	47"/24"	46.25"/25.5"	N/A	N/A	N/A
(4-wheel) Length/Width ⁶	39.75"/19.25"	39.75"/19.25"	44.75"/21.5"	N/A	44"/21"	N/A	47"/24"	46.75"/24"	49"/24"	46.875"/25.5"	52"/24"	63.5"/31"	52"/24"
Seat to Ground Height Range	20.5"-22"	20.5"-22"	23"-25"	20"-21"	20"-22"	18"-23"	22"-25"	23"-24.25"	21.5"-23.5"	22"-24"	23"-26"	27.75"-31.75"	23"-26"
Seat Type (standard)	COMPACT, FOLDABLE (17"x17")	COMPACT, FOLDABLE (17"x17")	FOLDABLE MOLDED PLASTIC (18"x16")	FOLDABLE MOLDED PLASTIC (18"x16")	FOLDABLE MOLDED PLASTIC (18"x16")	DELUXE, CONTOURED LOW-BACK (17"x16")	CRS DELUXE RECLINING HIGH-BACK W/SLIDERS & HEADREST (18"x18")	DELUXE MOLDED PLASTIC W/ SLIDERS (19"x16")	FOLDABLE MOLDED PLASTIC (18"x16")	DELUXE CONTOURED MEDIUM-BACK (22"x18.5") (18"x18")	CRS DELUXE RECLINING HIGH-BACK W/SLIDERS & HEADREST (18"x18")	DX CONTOURED HIGH-BACK WITH ADJ. HEADREST & SLIDERS (18"x16")	CRS DELUXE RECLINING HIGH-BACK W/SLIDERS & HEADREST (18"x18")
(3/4-wheel) Weight ³	72.5 LBS./76.5 LBS.	75 LBS./79 LBS.	106.5 LBS./125 LBS.	97.5 LBS./N/A	93 LBS./116 LBS.	148 LBS./N/A	155 LBS./166 LBS.	141.5 LBS./151 LBS.	110.5 LBS./116.5 LBS.	175 LBS./184 LBS.	N/A/207 LBS.	N/A/216 LBS.	N/A/209 LBS.
Seat to Deck Height Range	15"-16.5"	15"-16.5"	16.25"-18.25"	14"-15"	16.5"-18.5"	15.5"-20.5"	17"-20"	16.25"-17.5"	17"-19"	16"-18"	16.5"-19.5"	19"-23"	16.5"-19.5"
Battery Weight	9 LBS. EA.	9 LBS. EA./14 LBS. EA. ¹⁰	24.5 LBS. EACH ⁵	9 LBS. EACH	14 LBS. EA./24.5 LBS. EA.	24.5 LBS. EACH	24.5 LBS. EACH	24.5 LBS. EA./32.5 LBS. EA.	24.5 LBS. EACH	24.5 LBS. EA./38 LBS. EA.	38 LBS. EACH	63.5 LBS. EACH	38 LBS. EACH
Battery Requirements ⁹	(2) 12AH	(2) 12AH OR (2) 17AH	(2) U-1	(2) 12AH	(2) 17AH OR (2) U-1	(2) U-1	(2) U-1	(2) U-1 OR (2) 40AH	(2) U-1	(2) U-1 OR (2) NF-22	(2) U-1 OR (2) NF-22	(2) GROUP 27	(2) U-1 OR (2) NF-22
Per-Charge Range (max.) ²	UP TO 10 MILES	UP TO 10 MILES ¹²	UP TO 20 MILES	UP TO 10 MILES	UP TO 25 MILES	UP TO 20 MILES	UP TO 25 MILES	UP TO 25 MILES	UP TO 25 MILES	UP TO 30 MILES	UP TO 30 MILES	UP TO 25 MILES	UP TO 30 MILES
Battery Charger	OFF-BOARD 2A	OFF-BOARD 2A	OFF-BOARD 5A	OFF-BOARD 2A	ON-BOARD 3A	ON-BOARD 3A	ON-BOARD 3A	OFF-BOARD 5A	ON-BOARD 3A	OFF-BOARD 5A	OFF-BOARD 5A	ON-BOARD 5A	OFF-BOARD 5A
Suspension	N/A	N/A	N/A	N/A	N/A	N/A	REAR	(3-WHEEL) REAR (4-WHEEL) FRONT AND REAR	N/A	REAR	FRONT AND REAR	REAR	FRONT & REAR
Maximum Weight Capacity	260 LBS. MAX	260 LBS. MAX	300 LBS. MAX	300 LBS. MAX	300 LBS. MAX	300 LBS. MAX	350 LBS. MAX	350 LBS. MAX	350 LBS. MAX	500 LBS. MAX	400 LBS. MAX	400 LBS. MAX	400 LBS. MAX
Standard Body Colors	RED, BLUE ⁸	RED, BLUE, SILVER ¹³	RED, BLUE	RED, BLUE ⁸	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CHAMPAGNE, CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CHAMPAGNE, CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE
Available Power Elevating Seat	N/A	N/A	N/A	N/A	N/A	STANDARD	YES	YES (SPECIAL ORDER)	N/A	YES (SPECIAL ORDER)	YES (SPECIAL ORDER)	YES (SPECIAL ORDER)	YES (SPECIAL ORDER)
Weight of Heaviest Piece	28 LBS. (3-WHEEL) (FRONT SECTION) 33.5 LBS. (4-WHEEL) (FRONT SECTION)	28 LBS. (3-WHEEL) (FRONT SECTION) 33.5 LBS. (4-WHEEL) (FRONT SECTION)	48 LBS. (3-WHEEL) (REAR SECTION) 51.5 LBS. (4-WHEEL) (FRONT SECTION)	35 LBS. (FRONT OR REAR SECTION)	45 LBS. (3-WHEEL) (REAR SECTION) 45 LBS. (4-WHEEL) (FRONT/REAR SECTION)	71 LBS. (REAR SECTION)	61 LBS. (REAR SECTION)	56 LBS. (3-WHEEL) (FRONT SECTION) 65.5 LBS. (4-WHEEL) (FRONT SECTION)	50 LBS. (REAR SECTION)	64 LBS. (3-WHEEL) (REAR SECTION) 72.5 LBS. (4-WHEEL) (FRONT SECTION)	82 LBS. (REAR SECTION)	N/A	82 LBS. (FRONT SECTION AND REAR SECTION)
	TRAVEL		PORTABLE		MID-SIZE LUXURY		FULL-SIZE LUXURY			HEAVY-DUTY			PMVs

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.