



ENHANCED! *Sundancer*[®]

UNMATCHED MOBILITY WITH AN ADDED DIMENSION
3-WHEEL SCOOTER



Contemporary styling and superb performance complement one another perfectly in the **Sundancer** scooter. The Sundancer offers power, range and comfort plus a wealth of standard convenience features—like a standard power elevating seat—that help you accomplish your daily activities with ease. Highly maneuverable indoors, with the power, durability and balance for outdoor use.

The **Sundancer** delivers mobility plus! Just think of all you'll accomplish with the help of the Sundancer's *standard power elevating seat*:

- Instantly be at eye level with your family, friends and co-workers to engage in conversation and enjoy more intimate social interaction
- Gain easier access to your environment to more readily accomplish all of your daily activities
- Increase your independence—reach the things you need without having to ask for assistance
- Assist with transfers
- ***Take control of your life!***



“Power seat elevators for many of my clients are critical and absolutely medically necessary. They allow my clients to reach and carry out mobility related activities of daily living at different heights, without causing as much strain on the back and upper extremities.”

Mark Schmeler, MS, OTR/L, ATP,
University of Pittsburgh

800-800-8586 (U.S.)
888-570-1113 (Canada)
www.pridemobility.com

**MOBILITY & INDEPENDENCE
WITHIN ONE PRODUCT**

Pride[®]
Mobility Products Corp.

Specifications

Model number(s)	3-wheel: SC2000PS 4-wheel: N/A
Weight capacity	300 lbs.
Maximum speed	Up to 4.5 mph
Ground clearance	3.5"
Turning radius	3-wheel: 38" 4-wheel: N/A
Overall length	3-wheel: 41.625" 4-wheel: N/A
Overall width	23.25"

Tires

Front	2.5" x 9" solid
Rear	2.5" x 9" solid
Anti-Tips	Rear
Standard seating	Deluxe contoured low-back (17"W x 16" D) w/power elevating seat
Suspension	N/A
Drive train	Sealed transaxle, 24-volt DC motor
Dual braking system	Electronic, regenerative and electromechanical
Battery charger	On-board, 3 amp
Per-charge range	Up to 20 miles
Battery requirements	(2) 12V deep-cycle
Compatible battery sizes	U-1 (AGM or gel cell type recommended)

Weight

Batteries (2 required)	U-1: 24.5 lbs.
Total w/o batteries	3-wheel: 148 lbs.
Heaviest piece	71 lbs. (rear)

Colors



Candy
Apple
Red



Viper
Blue

Benefits

- Standard power elevating seat
- Adjustable tiller for driver comfort
- Full directional lighting package
- Contoured low-back seat with adjustable back angle and flip-up armrests for comfortable positioning
- Maximum weight capacity of 300 lbs.
- On-board battery charger — charge it wherever you go; up to 20-mile range per battery charge
- Speed range adjustable up to 4.5 mph
- Easy to maneuver through hallways or tight spaces in the home
- Spacious front basket standard
- Mag wheels for great style; 9-inch front and rear tires for enhanced outdoor performance
- Easy to disassemble for transport
- Extensive range of options and accessories available
- 3-year limited warranty

Options

- Deluxe contoured medium-back seat
- Forearm crutch holder
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Rear basket
- Safety flag
- Walker holder
- Heavy-duty dust cover
- Saddlebag
- Lap belt
- Cup holder
- Weather breaker canopy
- Quad cane holder

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Measurements are typically +/- 1/8"; weights +/- 1 lbs.