



# Victory<sup>TM</sup>

## 3-and 4-wheel scooters



Determined focus led to the ideal balance of hot style, solid performance and exceptional value found in the Victory<sup>®</sup> Scooter—a focus that provides you with a truly fun-loving scooter equipped standard with loads of luxury features at an affordable price. Live large even on a tight budget with a Pride<sup>®</sup> Victory.

*Stylish  
design  
&  
premium  
performance  
at an  
excellent  
value.*



VICTORY<sup>TM</sup> 4  
SC-1700  
(Candy Apple Red)

888-570-1113  
[www.pridemobility.com](http://www.pridemobility.com)

**Pride**  
Mobility Products Co.

## Specifications

Model number(s)	SC1600: (3-wheel) SC1700: (4-wheel)
Weight capacity	350 lbs.
Maximum speed	Up to 5 mph
Ground clearance	4"
Turning radius	3-wheel: 41" 4-wheel: 58"
Overall length	3-wheel: 47" 4-wheel: 49"
Overall width	3-wheel: 24" 4-wheel: 24"

## Tires

Front	3.5" x 10"
Rear	3.5" x 10"
Anti-tips	None
Standard seating	High back with headrest or Medium back
Power elevating seat	Option: No
Suspension	NA
Drive train	Sealed transaxle, 24 volt DC motor
Dual braking system	Electronic, regenerative and electromechanical
Battery charger	On-board, 3 amp
Battery requirements	(2) 12V deep-cycle (AGM or gel cell)
Battery size	U-1, 32AH

## Weight

Battery	32AH: 24.5 lbs.
Total w/o batteries	(3-wheel): 102 lbs. (4-wheel): 108 lbs.
Heaviest piece when disassembled	45 lbs. (rear section)

## Colors



Candy  
Apple  
Red



Viper  
Blue

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

## Features

- Single point lock-up for easy scooter disassembly
- Ultra-sleek, flush-mount oval headlight with mirror finish
- Restyled sporty alloy mag wheels
- Heavy-duty front (4-wheel) and rear polyurethane foam rubber bumpers
- State-of-the-art, microprocessor-based controller for superior performance & efficiency
- Controller features:
  - "Push-too-fast" protection
  - High Pedal Disable (HPD)
  - Internal fault detection
  - Anti-roll back/forward
  - Load compensation
- Water and shock-resistant Guest battery charger
- Convenient, externally mounted charger ammeter
- Standard front basket
- On-board accessory bracket
- Easily accessible charger port and circuit breaker

## Options

- Forearm crutch holder
- Lap belt
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Heavy-duty dust cover
- Saddle bag
- Solid tires
- Quad cane holder
- Clip, cane attachment



VICTORY™  
SC-1600  
(Viper Blue)

**Pride**  
Mobility Products Co.

380 Vansickle Road  
St. Catharines, Ontario  
**1-888-570-1113**  
[www.pridemobility.com](http://www.pridemobility.com)