

# Chameleon

## Chameleon - GK8/GK83

Shoprider's new Chameleon the fully transportable scooter that changes with you. Taking the Chameleon apart is as simple as 1-2-3, making it easy for you to take anywhere or store anyplace. Feel like a different colour for your dinner party tonight, no problem. The Chameleon comes with 3 different body panel colours, which can be quickly & easily changed. Chameleon is the first scooter that lets you change the color scheme of your whole scooter in just seconds.



## Built with the highest standards in quality and reliability

Since 1989 AMS HOMECARE INC. has built an enviable reputation by delivering high quality, reliable and comfortable self-care equipment addressing the needs of thousands of customers seeking to maintain their freedom and independence. We are dedicated to providing innovative, reasonably priced solutions specially developed with discriminating consumers in mind.

The Shoprider® brand represented exclusively in Canada by AMS HOMECARE INC. is recognized worldwide for superior design, technical sophistication as well as functionality in mobility products. Through warehouse-service centers conveniently located across the country, AMS HOMECARE INC. is the leader in the distribution of Shoprider® scooters and accessories.

In addition, we will be offering an expanding product line of AMS-branded products and accessories such as our current scooter canopy & cane holders and bed rails.

AMS HOMECARE INC. is a growing company committed to bringing to market the best products, devices and solutions to satisfy mobility and homecare needs.



Body style may differ from picture.

## Features:

- Portable 3 or 4-wheel scooter
- Ultra lightweight
- Easy as 1-2-3 to dismantle
- Interchangeable body panels

## Colours:



Red



Blue



Champagne

## Specifications Chameleon GK8 and GK83:

		Model GK8	Model GK83
Dimensions: L x W x H	<i>mm</i>	1120 x 550 x 840	1080 x 560 x 889
	<i>inches</i>	44 x 22 x 33	42.5 x 22 x 35
Number, size of tires	<i>Front</i>	2 x 7"	1 x 8"
	<i>Rear</i>	2 x 8"	2 x 8"
Weight capacity	<i>kg (lb)</i>	114 (250)	114 (250)
Total weight (with 12 Ah batteries, seat)	<i>kg (lb)</i>	(43) 94	(41) 89
Battery capacity	<i>_V_Ah x Pcs</i>	12v12ah x 2	12v12ah x 2
Charger (on/off board)	<i>V/A</i>	2 Amp, off board	2 Amp, off board
Output power of motor	<i>hp</i>	0.5	0.5
Maximum speed	<i>km/hr (mph)</i>	6 (3.8)	6 (3.8)
Max. driving distance with each full charge (depending on weight and terrain)	<i>km (miles)</i>	12 (7.5)	12 (7.5)
Turning radius	<i>mm (inches)</i>	940 (37)	838 (33)
Colours included (base is black)		Red, Blue, Champagne	Red, Blue, Champagne

All specification subject to change without prior notice.

Contact your AMS Homecare Dealer today:



AMS Homecare Inc. 1360 Cliveden Avenue, Delta, BC V3M 6K2 Canada  
Telephone: 604.273.5173 • Fax: 604.273.9312 • www.amshomecare.com

**SHOPRIDER**