

# SHOPRIDER®

## The New Generation



### **Cobra - 778EL**

Shoprider's new mid sized Cobra 778EL is truly cutting edge. Sleek new aerodynamic body design along with new features, including electric seat. The Cobra offers 4-wheel stability with 3-wheel maneuverability. Perfect for indoor and outdoor use.

## **Built with the highest standards in quality and reliability**

Since 1989 AMS HOMECARE INC. has built an enviable reputation by delivering high quality, reliable and comfortable self-care equipment addressing the needs of thousands of customers seeking to maintain their freedom and independence. We are dedicated to providing innovative, reasonably priced solutions specially developed with discriminating consumers in mind.

The Shoprider® brand represented exclusively in Canada by AMS HOMECARE INC. is recognized worldwide for superior design, technical sophistication as well as functionality in mobility products. Through warehouse-service centers conveniently located across the country, AMS HOMECARE INC. is the leader in the distribution of Shoprider® scooters and accessories.

In addition, we will be offering an expanding product line of AMS-branded products and accessories such as our current scooter canopy & cane holders and bed rails.

AMS HOMECARE INC. is a growing company committed to bringing to market the best products, devices and solutions to satisfy mobility and homecare needs.



### Specifications Cobra (model 778EL):

Overall dimension: L* x W x H	mm (in)	1117 (44) x 571 (22.5) x 1041 (41)
Number, size of tires (3-wheel chassis)	Front	Front: 2 x 8"
	Rear	Rear: 2 x 10"
Weight capacity	kg (lb)	136.1 (300)
Total weight (with batteries, basket, seat)	kg (lb)	73 (160)
Battery weight	kg (lb)	2 x 11.4 (2 x 25)
Battery capacity	_V_Ah x Pcs	12V34Ah x 2
Charger (on/off board)	V/A	24V/3A, off board
Output power of motor	hp	0.8
Maximum speed	km/hr (mph)	8.0 (5.0)
Safe climbing angle	degree	8° 14%
Max. driving distance with each full charge (depending on weight and terrain)	km (mile)	38 (24)
Turning radius	inch	38 in
Suspension	Yes/No	No
Detachable chassis	Yes/No	Yes (& stretchable)
Electric elevating seat	Yes/No	Yes
Handle bar type		C type delta

*All specification subject to change without prior notice.*

Contact your AMS Homecare Dealer today:


### Features:

- Electric seat lift
- Sliding swivel seat with manual seat recline
- Delta tiller
- Articulating front wheels
- Full light package (headlights, tail lights, signal & hazard)

### Colours:



Silver



Burgundy



AMS Homecare Inc. 1360 Cliveden Avenue, Delta, BC V3M 6K2 Canada  
Telephone: 604.273.5173 • Fax: 604.273.9312 • www.amshomecare.com

**SHOPRIDER**