



SHOPRIDER®

The New Generation

Landcruiser - 889XLS

The Landcruiser by Shoprider is truly a new breed of heavy-duty scooters. The 889XLS offers unprecedented style, with a weight capacity capable of up to 500lbs. With a 2hp motor, deluxe 20" wide high back captain's seat and generous legroom, you now have a scooter that offers incredible power and solid dependability.



Built with the highest standards in quality and reliability

Since 1989 AMS HOMECARE INC. has built an enviable reputation by delivering high quality, reliable and comfortable self-care equipment addressing the needs of thousands of customers seeking to maintain their freedom and independence. We are dedicated to providing innovative, reasonably priced solutions specially developed with discriminating consumers in mind.

The Shoprider® brand represented exclusively in Canada by AMS HOMECARE INC. is recognized worldwide for superior design, technical sophistication as well as functionality in mobility products. Through warehouse-service centers conveniently located across the country, AMS HOMECARE INC. is the leader in the distribution of Shoprider® scooters and accessories.

In addition, we will be offering an expanding product line of AMS-branded products and accessories such as our current scooter canopy & cane holders and bed rails.

AMS HOMECARE INC. is a growing company committed to bringing to market the best products, devices and solutions to satisfy mobility and homecare needs.



Body style may differ from picture.

Specifications Landcruiser (model 889XLS):

Overall dimensions: L x W x H	mm (in)	1430 (56.3) x 700 (27.6) x 1120 (44)
Number, size of tires: (black non-marking tires)	Front Rear	Front: 2, 13" Rear: 2, 13"
Weight capacity	kg (lb)	226.8 (500)
Total weight (w/batteries)	kg (lb)	140 (309)
Battery capacity	A.H.	33 ~ 75A.H.
Charger (on/off board)	A	8A, off board
Output power of motor	hp	2
Maximum speed	km/hr (mph)	13 (8)
Safe climbing angle	degree (%)	8° (14)
Max. driving distance with each full charge (depending on weight and terrain)	km (mile)	55 (34)
Turning radius	inch	59 in
Suspension	Yes/No	Yes (hydraulic & coil)
Detachable chassis	Yes/No	Yes
Emergency hand brake	Yes/No	Yes
20" wide, high back seat	Yes/No	Yes
Handle bar type		Delta type

All specification subject to change without prior notice.

Contact your AMS Homecare Dealer today:

Features:

- Heavy-duty transmission
- Full suspension (heavy-duty)
- 2 hp motor
- 20" wide sliding swivel seat with manual seat recline
- Full light package (headlights, tail lights, signal & hazard)
- Delta tiller

Colours:



Silver



Burgundy



AMS Homecare Inc. 1360 Cliveden Avenue, Delta, BC V3M 6K2 Canada
Telephone: 604.273.5173 • Fax: 604.273.9312 • www.amshomecare.com

SHOPRIDER