

# SHOPRIDER®

## The New Generation

### Sovereign - 888-3

Our value mid-sized three wheeled scooter offers a greater travel range while maintaining its ability to manoeuvre. Perfect for outdoor and moderate indoor use.



### Sovereign - 888-4

Our value mid-sized four wheeled scooter is one of North America's favourites. Perfect for outdoor and moderate indoor use.



## Specifications Sovereign 888-3 and 888-4:

		Model 888-3	Model 888-4
Overall dimensions: L x W x H	mm	1143 x 609 x 965	1244 x 533 x 939
	inches	45 x 24 x 38	49 x 21 x 37
Number, size of tires	Front	1, 10"	2, 10"
	Rear	2, 10"	2, 10"
Weight capacity	kg (lb)	136 (300)	136 (300)
Battery capacity	_V_Ah x Pcs	12v34ah x 2	12v34ah x 2
Charger (on/off board)	V/A	24v/3a, off board	24v/3a, off board
Output power of motor	HP	1.0	0.8
Maximum speed	km/hr (mph)	8 (5)	8 (5)
Safe climbing angle	degree (%)	8° (14%)	8° (14%)
Max. driving distance with each full charge (depending on weight and terrain)	km (miles)	40 (25)	40 (25)
Turning radius	inches	37	50
Suspension	Yes/No	No	No
Detachable chassis	Yes/No	Yes	Yes
Electric seat	Yes/No	No	No
Handle bar type		T bar or Delta tiller	T bar or Delta Tiller

*All specification subject to change without prior notice.*

### Features (Model 888-3):

- Five position locking tiller
- Deluxe Captain's seat
- Swivel seat with sliders
- Height adjustable arms
- Front head light
- Tail lights

### Colours:

Burgundy



### Features (Model 888-4):

- Five position locking tiller
- Deluxe Captain's seat
- Swivel seat with sliders
- Height adjustable arms
- Front head light
- Tail lights

### Colours:

Burgundy



## Built with the highest standards in quality and reliability.

Since 1989 AMS HOMECARE INC. has built an enviable reputation by delivering high quality, reliable and comfortable self-care equipment addressing the needs of thousands of customers seeking to maintain their freedom and independence. We are dedicated to providing innovative, reasonably priced solutions specially developed with discriminating consumers in mind.

The Shoprider® brand represented exclusively in Canada by AMS HOMECARE INC. is recognized worldwide for superior design, technical sophistication as well as functionality in

mobility products. Through warehouse-service centers conveniently located across the country, AMS HOMECARE INC. is the leader in the distribution of Shoprider® scooters and accessories. In addition, we will be offering an expanding product line of AMS-branded products and accessories such as our current scooter canopy & cane holders and bed rails.

AMS HOMECARE INC. is a growing company committed to bringing to market the best products, devices and solutions to satisfy mobility and homecare needs.

Contact your AMS Homecare Dealer today:



AMS Homecare Inc. 1360 Cliveden Avenue, Delta, BC V3M 6K2 Canada  
Telephone: 604.273.5173 • Fax: 604.273.9312 • www.amshomecare.com

# SHOPRIDER