

SHOPRIDER®

The New Generation



Trailblazer - 889SEL/889SL

The Trailblazer one of Shoprider's largest scooter with a 1.2hp motor and 13" wheels, the Trailblazer provides a comfortable ride, even on rough terrain. The Trailblazer models have a weight capacity of up to 350lbs, along with a heavy-duty suspension.

Built with the highest standards in quality and reliability

Since 1989 AMS HOMECARE INC. has built an enviable reputation by delivering high quality, reliable and comfortable self-care equipment addressing the needs of thousands of customers seeking to maintain their freedom and independence. We are dedicated to providing innovative, reasonably priced solutions specially developed with discriminating consumers in mind.

The Shoprider® brand represented exclusively in Canada by AMS HOMECARE INC. is recognized worldwide for superior design, technical sophistication as well as functionality in mobility products. Through warehouse-service centers conveniently located across the country, AMS HOMECARE INC. is the leader in the distribution of Shoprider® scooters and accessories.

In addition, we will be offering an expanding product line of AMS-branded products and accessories such as our current scooter canopy & cane holders and bed rails.

AMS HOMECARE INC. is a growing company committed to bringing to market the best products, devices and solutions to satisfy mobility and homecare needs.



Body style may differ from picture.

Features:

- Full suspension (heavy-duty)
- Electric seat lift (889SEL)
- Sliding swivel seat with manual seat recline
- Delta tiller
- Full light package (headlights, tail lights, signal & hazard)
- New aerodynamic body designs

Colours:



Royal Blue



Burgundy

Specifications Trailblazer (models 889SEL/889SL):

Overall dimensions: L* x W x H	mm (in)	1244 (49) x 640 (25) x 1120 (44)
Number, size of tires (black, non-marking tires)	Front	Front: 2, 13"
	Rear	Rear: 2, 13"
Weight capacity	kg (lb)	158.8 (350)
Total weight (with batteries, basket, seat)	kg (lb)	100 (221)
Battery weight	kg (lb)	2 x 11.4 (2 x 25)
Battery capacity	_V_Ah x Pcs	12V - 34Ah (889SEL) 12V - 50Ah (889SL)
Charger (on/off board)	V/A	24V/3A, off board
Output power of motor	hp	1.2
Maximum speed	km/hr (mph)	13 (8)
Safe climbing angle	degree	8° 14%
Max. driving distance with each full charge (depending on weight and terrain)	km (mile)	31 (20) (889SEL)
		48 (30) (889SL)
Turning radius	inch	55 in
Suspension	Yes/No	Yes
Detachable chassis	Yes/No	Yes
Electric elevating seat	Yes/No	889SEL
Handle bar type		Type delta

All specification subject to change without prior notice.

Contact your AMS Homecare Dealer today:



AMS Homecare Inc. 1360 Cliveden Avenue, Delta, BC V3M 6K2 Canada
Telephone: 604.273.5173 • Fax: 604.273.9312 • www.amshomecare.com

SHOPRIDER