



# SHOPRIDER®

## The New Generation



### **Voyager - Model 888SEL**

Shoprider's new mid-sized Voyager 888SEL are full of industry leading features. Sleek aerodynamic body design, full suspension, and electric seat lift. The Voyager offers 4-wheel stability. Perfect for outdoor and moderate indoor use.

## **Built with the highest standards in quality and reliability**

Since 1989 AMS HOMECARE INC. has built an enviable reputation by delivering high quality, reliable and comfortable self-care equipment addressing the needs of thousands of customers seeking to maintain their freedom and independence. We are dedicated to providing innovative, reasonably priced solutions specially developed with discriminating consumers in mind.

The Shoprider® brand represented exclusively in Canada by AMS HOMECARE INC. is recognized worldwide for superior design, technical sophistication as well as functionality in mobility products. Through warehouse-service centers conveniently located across the country, AMS HOMECARE INC. is the leader in the distribution of Shoprider® scooters and accessories.

In addition, we will be offering an expanding product line of AMS-branded products and accessories such as our current scooter canopy & cane holders and bed rails.

AMS HOMECARE INC. is a growing company committed to bringing to market the best products, devices and solutions to satisfy mobility and homecare needs.



Body style may differ from picture.

Two 34A-12V batteries

### Specifications Voyager (model 888SEL):

|   |             |                   |
|---|-------------|-------------------|
| Overall dimensions: L* x W x H  | inches      | 48 x 22.5 x 43    |
| Number, size of tires   | Front       | Front: 2, 10"     |
|   | Rear        | Rear: 2, 10"      |
| Weight capacity   | kg (lb)     | 136.1 (300)       |
| Total weight (with batteries, basket, seat)                                   | kg (lb)     | 95 (210)          |
| Battery weight  | kg (lb)     | 2 x 11.4 (2 x 25) |
| Battery capacity  | _V_Ah x Pcs | 12V34Ah x 2       |
| Charger (on/off board)  | V/A         | 24V/3A, off board |
| Output power of motor   | hp          | 0.8               |
| Maximum speed   | km/hr (mph) | 8 (5.0)           |
| Safe climbing angle   | degree (%)  | 8° (14%)          |
| Max. driving distance with each full charge (depending on weight and terrain) | km (mile)   | 38 (24)           |
| Turning radius  | inch        | 50 in             |
| Suspension  | Yes/No      | Yes               |
| Detachable chassis  | Yes/No      | Yes               |
| Electric elevating seat   | Yes/No      | Yes               |
| Handle bar type   |             | Type delta        |

All specification subject to change without prior notice.

Contact your AMS Homecare Dealer today:


### Features:

- Full light package
- New Aerodynamic body design
- Delta tiller
- Electric seat lift
- Sliding swivel seat with seat recline
- Full suspension (heavy duty)

### Colours:



Royal Blue



Burgundy



AMS Homecare Inc. 1360 Cliveden Avenue, Delta, BC V3M 6K2 Canada  
Telephone: 604.273.5173 • Fax: 604.273.9312 • www.amshomecare.com

**SHOPRIDER**