



Shock-absorbing "cradle effect" protects power chairs from jarring vibrations.



All-new design dramatically reduces metal to metal contact for enhanced performance and longer product life.



ENHANCED!



Easy to use, three-point power chair tie down system allows for one handed operation and is self-retracting and self-tensioning.



Innovative, lightweight self-leveling platform allows for easy loading and unloading, even on inclines or uneven terrain.

FULL PLATFORM EXTERIOR LIFT

Your on-the-go lifestyle demands that you bring your power chair with you. Now, Pride Mobility Products® Corp., the number one power chair and scooter manufacturer, offers the lift you need to take your mobility product on the road. Backed by the service you have come to expect and trust from Pride®, with the full support of our expert technical service staff, a Silver Star™ Lift is the natural complement to your Pride Power Chair.

The perfect fusion of Pride's engineering achievement and Silver Star Lifts' field-proven expertise brings you the most advanced, durable and feature-filled exterior lift system available.

A Silver Star exterior lift makes it easy to take your Pride Power Chair with you wherever you go, complementing your product perfectly and enhancing your mobility. An all-steel frame and protective coating mean Silver Star Lifts are built to be durable, while quality construction means you can count on your lift to be dependable. From loading to lifting, Silver Star Lifts are designed with ease of use in mind.

- 325 lbs. lifting capacity *(350 lbs. coming soon)*
- 3-year transferable warranty
- Most adaptable exterior lift available

SilverStar™ Lifts

A Division of **Pride**® Mobility Products Corp.

182 Susquehanna Ave. • Exeter, PA 18643 • (800) 800-8586 • www.silverstarmobility.com • www.pridemobility.com

FULL PLATFORM EXTERIOR LIFT



Unoccupied lift stowed.



Lift down, ready for product.



Loaded product secured and ready for lifting.



Product lifted and ready for transport.

SPECIFICATIONS

Base weight lift	107 lbs. ¹
Power chair overall weight	325 lbs.
Maximum product width	27"
Maximum product wheelbase	36"
Maximum travel distance from the ground	20"
Motor	Sealed 12 volt DC
Hitch class	Class II or III ²
Vehicle platform	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover
Warranty	3-year transferable



Can easily lift the following Pride® Power Chairs³

- Jazzy® 600
- Jazzy® 610
- Jazzy® 1103 Ultra
- Jazzy® 1107
- Jazzy® 1113 ATS
- Jazzy® 1121
- Jazzy® 1122
- Jazzy® Select
- Jazzy® Select 14
- Jazzy® Select 14 XL
- Jet™ 3 Ultra
- Jet™ 7
- Jet™ 10 Ultra
- Go-Chair®

Can easily lift the following Quantum® Power Chairs³

- Quantum® 600
- Quantum® 610
- Quantum® 1107
- Quantum® 1121
- Quantum® 1402
- Quantum® 6000
- Quantum® 6000 XL

Options

- License plate kit
- Optional battery pack with charger
- Swing-away joint
- Extended swing-away joint
- Power chair cover

¹ Weight does not include any accessories
² 190 lbs. maximum product weight with a Class II hitch
³ Vehicle dependent

Note: The Full Platform Exterior Lift can lift all of the listed power chairs with contoured, Synergy® or Specialty seating. Please call Customer Service when requesting information about lifting with a power positioning system.

SilverStar™ Lifts

A Division of **Pride** Mobility Products Corp.

182 Susquehanna Ave. • Exeter, PA 18643 • (800) 800-8586 • www.silverstarmobility.com • www.pridemobility.com