



Shown: Power Chair Lift with Straps and Jazzy[®] 1103 Ultra

POWER CHAIR EXTERIOR LIFTS

Silver Star[™] Power Chair Lifts from Pride[®] are the safe, reliable way you can take your power chair on the road. Quality construction means you can count on your lift to be dependable, while an all-steel frame and protective coating mean Silver Star Lifts are built to be durable. From loading to lifting, Silver Star Lifts are designed with ease of use in mind.

- 325 lbs. lift capacity
- 3-year transferable warranty
- Optional with automatic docking station for securing power chair to lift

SilverStar[™] Lifts 

A Division of **Pride** Mobility Products Corp.

POWER CHAIR EXTERIOR LIFTS



Power Chair Lift in lowest position.



Power Chair Lift with Jazzy® 1103 Ultra raised.



Power Chair Lift folded.

SPECIFICATIONS

POWER CHAIR LIFT WITH STRAPS/POWER CHAIR LIFT WITH DOCKING STATION

Base weight lift	105/108 lbs. ¹
Power chair overall weight	325 lbs.
Maximum product width	27"
Maximum product length	36" ²
Maximum travel distance from the ground	20"
Motor	Sealed 12 volt DC
Hitch class	Class II or III ³
Vehicle platform	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover
Warranty	3-year transferable



Power Chair Lift with docking station (passive locking system that secures the power chair to the platform automatically)

Can easily lift the following Jazzy® Power Chairs⁴

- Jazzy® 1103 Ultra
- Jazzy® 1107
- Jazzy® 1113 ATS
- Jazzy® 1121
- Jazzy® 1122

Can easily lift the following Quantum® Power Chairs⁴

- Quantum® 1107
- Quantum® 1121
- Quantum® 1402

Can easily lift the following Jet™ Power Chairs⁴

- Jet™ 2
- Jet™ 2HD
- Jet™ 3 Ultra
- Jet™ 7
- Jet™ 10
- Jet™ 10 Ultra

Pride® Power Dock Adapters

- Jazzy® 1103 Ultra
- Jazzy® 1121
- Jazzy® 1122
- Jet™ 2
- Jet™ 3 Ultra

Note: If a specific Power Dock is not listed, please contact Customer Service at (800) 800-8586 for assistance.

Options

- License plate kit
- Optional battery pack with charger
- Swing-away joint
- Extended swing-away joint
- Power chair cover

¹ Weight does not include any accessories

² Without front riggings

³ 190 lbs. maximum unit weight on a Class II hitch

⁴ Vehicle dependent

Note: These lifts can lift all of the listed power chairs with contoured, Synergy® or Specialty seating. Please call Customer Service when requesting information about lifting with a power positioning system.

Silver Star™ Lifts

A Division of **Pride** Mobility Products Corp.

182 Susquehanna Ave. • Exeter, PA 18643 • (800) 800-8586 • www.silverstarmobility.com • www.pridemobility.com

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Due to manufacturing tolerances and continual product improvement, length and width can be subject to variance of + or - 3%.