



Shown: Power Tote with Stylus™

POWER TOTES

Silver Star™ Power Tote Carriers by Pride® are easy to install and operate. The PT-100 can lift your manual wheelchair with ease and the PT-150 is able to lift your small scooter or travel vehicle. These handy, lightweight lifts feature the same steel frame and quality construction as all Silver Star Lifts by Pride.

- Weigh just 47 lbs.
- Power Tote (PT-100) easily attaches to Class I hitch & lifts 100 lbs.
- Power Tote Mini (PT-150) easily attaches to Class II hitch & lifts 150 lbs.

SilverStar™ Lifts ★

A Division of **Pride**® Mobility Products Corp.

POWER TOTES



Power Tote



PT-100



Power Tote Mini



PT-150

SPECIFICATIONS

POWER TOTE (PT-100)

Base weight lift	47 lbs. ¹
Power chair overall weight	100 lbs.
Maximum product width	14" folded
Maximum travel distance from the ground	20"
Motor	Sealed 12 volt DC
Warranty	3-year transferable
Hitch class	Class I ³
Vehicle platform	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover
Can easily lift	Pride® Stylus™ and folding wheelchairs

POWER TOTE MINI (PT-150)

Base weight lift	47 lbs. ¹
Power chair overall weight	150 lbs.
Maximum product width	22"
Maximum travel distance from the ground	20"
Motor	Sealed 12 volt DC
Warranty	3-year transferable
Hitch class	Class II or Class III
Vehicle platform	Full-size Van, Mini Van, Truck, Sedan, SUV and Crossover
Can easily lift	3-wheel Go-Go® Ultra, Go-Go Elite Traveller™, or Sonic®

Options

- License plate kit
- Optional battery pack with charger ²
- Swing-away joint ⁴
- Extended swing-away joint ⁴
- Scooter cover

¹ Weight does not include any accessories

² Currently only available on the Power Tote Mini

³ If vehicle is equipped with Class II or Class III hitch, it requires a hitch reducer; additional charges apply

⁴ If using a Class III hitch

SilverStar™ Lifts 

A Division of **Pride** Mobility Products Corp.

182 Susquehanna Ave. • Exeter, PA 18643 • (800) 800-8586 • www.silverstarmobility.com • www.pridemobility.com