



Shown: Double Entry Scooter Lift with Revo®

SCOOTER EXTERIOR LIFTS

Add a whole new dimension to your mobility with a Silver Star™ Scooter Lift by Pride®, the reliable solution to taking your scooter on the road. An all-steel frame and protective coating mean all Silver Star Lifts are built to be durable, while quality construction means you can depend on your lift for the long haul.

- Double Entry and Quad XL lifts available
- 3-year transferable warranty
- 325 lbs. lifting capacity

Silver Star™ Lifts ★

A Division of **Pride** Mobility Products Corp.

182 Susquehanna Ave. • Exeter, PA 18643 • www.silverstarmobility.com • www.pridemobility.com
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

SCOOTER EXTERIOR LIFTS

SPECIFICATIONS



Double Entry Scooter Lift



Quad XL Scooter Lift

DOUBLE ENTRY SCOOTER LIFT

The Double Entry Scooter Lift allows you to enter and exit the lift platform from either side with your 3 or 4-wheel Pride® Scooter.

Base weight lift	106 lbs. ¹
Maximum scooter weight	325 lbs.
Maximum product width	27"
Maximum travel distance from the ground	20"
Motor	Sealed 12 volt DC
Hitch class	Class II or III ³
Vehicle platform	Full-size Van, Mini-van, Truck, Sedan, SUV and Crossover
Warranty	3-year transferable

QUAD XL SCOOTER LIFT

The Silver Star™ Quad XL Scooter Lift by Pride® is a sturdy, longer scooter lift for use with Pride's larger products, the Hurricane® and Legend® XL.

Base weight lift	116 lbs. ¹
Maximum scooter weight	325 lbs.
Maximum product width	27"
Maximum travel distance from the ground	20"
Motor	Sealed 12 volt DC
Hitch class	Class II or III ³
Vehicle platform	Full-size Van, Mini-van, Truck, Sedan, SUV and Crossover
Warranty	3-year transferable

¹ Weight does not include any accessories

² Vehicle dependent

³ 190 lbs. maximum scooter weight with Class II hitch

The Double Entry can easily lift the following Pride® Scooters ²

- Celebrity® X 3 & 4-wheel
- Go-Go Elite Traveller™ 3 & 4-wheel
- Go-Go® Ultra X 3 & 4-wheel
- Legend® 3 & 4-wheel
- Maxima 3 & 4-wheel
- Rally® 3 & 4-wheel
- Revo® 3 & 4-wheel
- Sundancer® 3-wheel
- Victory® 3 & 4-wheel

The Quad XL Scooter Lift can easily lift the following Pride® Scooters ²

- Hurricane®
- Legend® XL

Options

- License plate kit
- Optional battery pack with charger
- Swing-away joint
- Extended swing-away joint
- Scooter cover

Note: Silver Star Lifts, while designed specifically for Pride Mobility products, are capable of lifting other scooter brands.

Silver Star™ Lifts

A Division of **Pride** Mobility Products Corp.

182 Susquehanna Ave. • Exeter, PA 18643 • www.silverstarmobility.com • www.pridemobility.com
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)