

scalamobil® iQ

*The mobile
stairclimbing aid*



viamobil®

*Push and brake
assistance for
wheelchairs*



alber

INVACARE
Yes, you can.™

Alber® Power Add-on Products

scalamobil iQ

With the scalamobil iQ, obstructed paths are accessible again. No matter whether the staircases are steep, narrow or spiral, the scalamobil overcomes all types of stairs without any problems and can be operated with very little effort, both simply and safely. The scalamobil's unique design and patented sensor, make stair travel safe and comfortable for both the consumer and the caregiver. The scalamobil iQ, helping you reach new heights.



FEATURES

- The scalamobil iQ is eminently suitable for use outside one's home.
- With its practical quick-release locks it can easily be dismantled into three handy parts.
- The scalamobil iQ fits into any car trunk. It can easily be loaded and transported.
- The scalamobil iQ can be attached to most manual wheelchairs. Usually all that is required is a special bracket on the existing wheelchair to attach the scalamobil (not permanently).
- To climb stairs, the rear wheels are simply removed. The advantage is that the wheelchair occupant need not change seats.

Specifications

Battery pack	2 x 12 V; 3, 4 Ah
Battery charger (input)	90 - 240 V; AC
Battery charger (output)	24 V; DC
Operating voltage	24 V
Climbing speed	6 - 16 steps/min.
Motor rating	176 W
Max. weight of occupant	max. 120 kg; 264.6 lbs
Capacity per battery charger	approx. 300 steps
Drive unit (without battery)	16.6 kg; 36.6 lbs
Battery	2.7 kg; 5.9 lbs
Handle	5.2 kg; 11.5 lbs
Total weight	24.5 kg; 54.0 lbs
Width (with mounted handle)	48 cm

viamobil

Pushing a wheelchair can be a difficult and tiresome task. Let the viamobil do the work. Uphill, downhill, braking, the viamobil's powerful drive system helps you to travel through a variety of terrain. The viamobil's take apart design makes it easy to transport the viamobil unit. The viamobil, the caregivers choice for people on the move.



FEATURES

- The viamobil can be fitted to practically every wheelchair.
- An inconspicuous mounting system guarantees safe attachment of the viamobil to the wheelchair.
- The wheelchair will remain fully operational and the mounting brackets will not inhibit folding once viamobil has been removed.
- The powerful electric unit tucks away under the seat and is held by two quick-lock devices.
- The handy rechargeable battery pack is securely attached to the drive unit and is automatically locked into position.
- Automatic charger is supplied to charge battery.
- The viamobil is operated by an ergonomic control unit that mounts to the wheelchair handles.
- The viamobil relieves strain on the escorting person and increases freedom of movement.
- The drive unit is lifted and lowered at the press of a button.

Specifications

Battery pack	2 x 12 V; 12 Ah
Battery charger (input)	100 - 240 V; AC
Battery charger (output)	29.4 V; DC
Operating voltage	24 V
Climbing speed	1 - 5.5 km/h
Motor rating	150 W
Climbing and braking capacity	max 18%
Range per single charge	ca. 25 km
Maximum load	130 kg; 286 lbs
Battery pack	9 kg; 19.8 lbs
Drive unit	12.9 kg; 28.4 lbs
Control unit	0.9 kg; 2 lbs
Battery charger	0.9 kg; 2 lbs



Invacare Canada
www.invacare.ca

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

© 2005
Invacare Canada
Form No. 05-438C

alber



"Invacare says 'Yes, you can',
and I couldn't have said it
better myself."