

Alber®
e-motion
power assist

Alber® e-motion® power assist



The e•motion power wheels provide assistance to the consumer to maintain or raise their level of active mobility. The revolutionary design of e•motion enables the consumer to tackle upward and downward slopes without the fatigue and joint strain usually associated with propelling a manual wheelchair. The e•motion's self-contained system (noise-free, gearless motors) is discretely mounted in the wheel hubs. Whether for paediatric or adult use, there is no need to change wheelchairs – e•motion fits most manual chairs with a simple bracket installation. The e•motion wheels have a quick release feature, making them easy to interchange with regular wheels and also convenient for travel.

FEATURES

- Consists of 2 programmable power wheels (22" or 24"), 1 pair of batteries, automatic quick charger and mounting brackets
- 2 levels of assisted support (indoor/outdoor)
- 10 driving programs
- Automatic power save feature
- Maintenance-free, leak-proof 3.0 Ah NiMH batteries
- Battery charger included, with automatic charge and shut off feature
- Range: 8 miles approx (depending on assistance level and occupant weight)
- Speed: 0-6 km/h
- Total weight: 53 lbs. (excluding wheelchair)
- Weight capacity: 264 lbs.

ACCESSORIES

- Anti-tip supports
- Spoke covers
- Remote control
- Rubber-coated push rim covers
- Power converter
- Extra battery packs



alber



"Invacare says 'Yes, you can',
and I couldn't have said it
better myself."

Arnold Palmer

Invacare Canada
www.invacare.ca

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

© 2005
Invacare Canada
Form No. 05-436C