

Invacare® Pronto® M71™ Power Wheelchair

POWER • CUSTOM  
REHAB  
SEATING • STAIRS

# Pronto® M71™



Yes, you can.®



*Anything else is a compromise.™*

# Invacare® Pronto® M71™ Power Wheelchair

The Invacare® Pronto® M71™ power wheelchair offers flexibility to meet mild to moderate positioning needs in a compact base that provides excellent indoor maneuverability. Built on the foundation of SureStep® suspension and true center-wheel drive technology, the Pronto M71 offers proven performance and a take apart frame for maximum transportability.

## Center-Wheel Drive



Center-wheel drive offers the tightest turning radius possible, as well as improved tracking and ease of driving, over all other Invacare configurations.

## SureStep®



This patented automatic system ensures that everyday obstacles and thresholds are handled elegantly for a smooth ride.

## MK<sub>5</sub>™ & MK6i™ Electronics



MK<sub>5</sub>™ and MK6i™ electronics, the leading power wheelchair electronics developed and designed by Invacare, include a MK<sub>5</sub> NX™ (non-expandable controller) on the van seat and MK6i 60 amp controller on the rehab seat models with SPJ™+ (Single Purpose Joystick Plus). This combination allows for full turning capability at slow speeds and provides push button on/off, speed selection and horn, as well as a bright LED display.

## Transport



The Pronto M71 can be disassembled without tools for easier transport.

## Other Important Features

- 4-point seat mount
- 2-pole motors
- 3 amp on-board battery charger
- Accepts U1 batteries



