

Invacare® Pronto® M94™ Power Wheelchair



Yes, you can.™

# Pronto® M94™



*Anything else is a compromise.™*



## Pronto M94 with Van Seat

- Width range: 18"-24"
- Depth range: 16"-22"
- Semi-reclining back (standard) for additional comfort
- Width adjustable flip-back arms (standard) with optional height adjustment
- Ergonomic seating contours provide support for longer durations of use
- Adjustable height headrest for head and neck support
- Built-in seat pouch on back
- Standard easy to clean grey vinyl upholstery

## Van Seat Upholstery



Grey vinyl

## Pronto M94 with Rehab Seat

- Width range: 18"-24"
- Depth range: 16"-22"
- Angle adjustable back offers range 80° to 100°
- Sling back upholstery provides comfort and support for mild positioning needs
- Round canes accept other back options for additional positioning

## Shroud Colors



Candy Red



Deep Blue








Magic Black

# Invacare® Pronto® M94™ Power Wheelchair

The Invacare® Pronto® M94™ with SureStep® combines innovative features and exceptional quality to provide a power wheelchair that meets the mobility needs of bariatric consumers. The M94 features the patent-pending SureStep technology that makes negotiating thresholds and transitions of up to three inches high effortless. Also featured is the patent-pending Stability Lock technology that helps ensure that the chair will remain stable with the consumer in control. These two innovative features paired with its sleek styling make the Pronto M94 a truly exceptional power chair for consumers weighing up to 500 pounds.

Invacare® Pronto® M94™ Power Wheelchair

<p>Center-Wheel Drive</p>		<p>Center-wheel drive offers the tightest turning radius, as well as improved tracking and ease of driving, compared to our other configurations.</p>
<p>SureStep®</p>		<p>This patented automatic system ensures that everyday obstacles and thresholds are handled elegantly for a smooth ride.</p>
<p>Stability Lock</p>		<p>An important issue with more maneuverable chairs is stability, but this patented technology helps ensure that the chair remains stable, with the consumer in control.</p>
<p>TCD™ Traction Control Design</p>		<p>TCD improves traction in forward and reverse - especially important for transitioning onto ramps and backing up in grass.</p>
<p>4-Pole Motors</p>		<p>Proven Invacare® 4-pole motors provide high torque and true gear reduction for many rehab needs.</p>



Yes, you can.®

## Specifications

Overall dimensions		
Length (with footboard up)	39"	
Length (with footboard down)	40" - 42"	
Width (excluding joystick)	25.9"	
Height	37" - 48"	
Seat dimensions		
	Van seat	
	Semi-recline	Adjustable Back (ASBA)
Width	18" - 24"	18" - 24"
Depth	16" - 22"	16" - 22"
Overall Height	37" - 40"	35" - 43"
Back height	18"	16" - 20"
with headrest	25"	n/a
Back angle	80° - 115°	80° - 100°
Seat-to-floor height		
Van	21" - 23"	
ASBA	18.5" - 21"	
Weight		
With batteries	273 lb.	
Without batteries	199 lb.	
Wheel and caster sizes		
Front casters	6" x 2"	
Drive wheels	14" x 3"	
Rear casters	6" x 2"	
Maximum speed	4.5 mph	
Incline capability	9°	
Weight capacity	500 lb.	
Turning radius		
Front (with footboard down)	19.5"	
Rear	21.5"	
Ground clearance	3.25"	
Limited warranties		
Base frame	60 months	
Seat frame	60 months	
Electronics	12 months	
Motor / gearbox	18 months	

Invacare® Pronto® M94™ Power Wheelchair

**Invacare Corporation**  
www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

Invacare, the Blue Medallion Design, Yes, you can., Pronto and SureStep are registered trademarks and Anything else is a compromise, M94, MK5, NX, SPJ and TCD are trademarks of Invacare Corporation. Specifications are subject to change without notification.

Patent pending.

© 2006

Invacare Corporation  
Form No. 04-001  
Rev. 01/06



*"Invacare says 'Yes, you can,' and I couldn't have said it better myself."*