



NAME _____ PHONE NO _____

ACCOUNT # _____ P.O. _____ DATE _____

ADDRESS _____

FRAME SELECTION

CHAIR TYPE

- 1 3GTQSP Torque SP Base for Seat w/MK6i Electronics..... \$6,978
- 1 3GTQSPBASE Torque SP Base Only w/MK6i Electronics..... \$6,798
- 2 3GSEAT Seat Only..... See Seat Selection

3 TRANSPORT TIE DOWN

- TRRO Transport Ready Option \$388
- TRBKTS Wheelchair Transport Brackets \$325

USER WEIGHT LIMITS

- U300 User Weight <300 Lbs..... No Charge

FRAME TYPE

- 4 SFS Short Base..... No Charge
- 5 LFS Long Base..... \$403

6 ADJUSTABLE SEAT WIDTH RANGE

- 1216W 12" to 16" Seat Width Range..... No Charge
- 1620W 16" to 20" Seat Width Range..... No Charge
- 2022W 20" to 22" Seat Width Range..... \$429

ASBA FRAME WIDTH

- 7 W12 12" Seat Width..... No Charge
- 7 W13 13" Seat Width..... No Charge
- 7 W14 14" Seat Width..... No Charge
- 7 W15 15" Seat Width..... No Charge
- W16 16" Seat Width..... No Charge
- W17 17" Seat Width..... No Charge
- W18 18" Seat Width..... No Charge
- W19 19" Seat Width..... No Charge
- W20 20" Seat Width..... No Charge
- W21 21" Seat Width..... No Charge
- W22 22" Seat Width..... No Charge

VAN SEAT WIDTH

- VW18 18" Wide (17" Depth)..... No Charge
- VW20 20" Wide (19" Depth)..... \$209
- VW22 22" Wide (19" Depth)..... \$209

PERFORMANCE UPGRADE

- 8 TRUTK TrueTrack GB™ Driving Package..... \$1,465

SEAT SELECTION

- ADJASBA Seat w/Adj. Seat Tilt & Back Angle..... \$774
- VAN Van Style Captain Seat..... \$1,014
- ADJRECL Recliner..... \$1,499
- 9 CTS Contoura Seating \$553
- OMCU Omit Cushion..... Less Net (\$25)

* TRANSPORT TIE DOWN	
<p>TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.</p> <p>TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.</p> <p>As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.</p> <p>* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.</p>	
<ul style="list-style-type: none"> * Torque SE Base & Seat * 300 lb. user weight limit * Short or Long Base Frame * Adjustable Rehab Seat with Adjustable Seat Tilt & Back Angle * Any Seat Width * Any Seat Depth * Low Seat-to-Floor Height * Sectional or Straight Rear Uprights * 18" - 24" ASBA Back Height * 90° or 95° Back Angle * Any Arm Type * 8" x 2.25" Semi Pneumatic Casters * Suspension Shocks * 14" x 3" Tires with Foam Filled Inerts * Rear Wheel Standard Location * Any Color Frame Finish * Black Nylon or Black Vinyl Upholstery * Center Pivot Style Bracket * 60° Swingaway, 70° Swingaway or 70° Tapered Footrests * Group 24 Battery Box 	<ul style="list-style-type: none"> * True Track Upgrade * Fender Option * Flip-Up Composite Footplates, Aluminum Footplates-Deluxe, Adjustable Angle Flip Up Footplates * Calf Strap * Impact Guards * All PSR / PSF Joysticks * Left Hand Mounted Joystick * "T" Handle Flexible Joystick Ext. * Straight Handle Flexible Joystick Ext. * MK6i Remote Programmer * SPJ™+ or MPJ™+ Joysticks * Reinforced Back Upholstery * Group 24 Batteries
<p>Please note that all TRRO chairs are set at a 10° seat angle</p>	

FOOTNOTES

- 1 Torque SP comes standard with 300 lb. user weight capacity, 6.5 MPH top speed, performance 4-pole motor package, MK6i SPJ+ joystick, MK690 Controller. Base only comes standard with MK690ACC Controller. Base only does not come standard with Seat Positioning Strap.
- 2 When ordering a replacement seat and changing depth to 20" or deeper, order new arms.
- 3 Must select Group 24 battery boxes for TRRO or TRBKTS.
- 4 SFS is standard on seat depths 12" - 19".
- 5 LFS is standard on seat depths 20" - 22"
- 6 Jr. seat widths require a straight rear upright and have a user weight limit of 150lbs. Contoura Seating available in even widths from 16"-22".
- 7 Manual reclining back not available.
- 8 TrueTrack driving technology, which enables driver to track straighter regardless of the terrain, includes Monroe Spring Suspension, <300 lb. user weight limit, and 7 MPH maximum speed. Must pick MK6TT TrueTrack controller.
- 9 Selected seat depth on Contoura assumes 1" of back cushion compression. For example: When choosing an 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick. Please see page 7 for seat size offerings.

ADJUSTABLE SEAT DEPTH RANGE

- 10 1218D 12" to 18" Seat Depth Range..... No Charge
- 1619D 16" to 19" Seat Depth Range..... No Charge
- 1922D 19" to 22" Seat Depth Range..... \$355

ASBA SEAT DEPTH

- 11 D12 Seat Depth 12" (Short Base)..... No Charge
- 11 D13 Seat Depth 13" (Short Base)..... No Charge
- 11 D14 Seat Depth 14" (Short Base)..... No Charge
- 11 D15 Seat Depth 15" (Short Base)..... No Charge
- D16 Seat Depth 16" (Short Base)..... No Charge
- D17 Seat Depth 17" (Short Base)..... No Charge
- D18 Seat Depth 18" (Short Base)..... No Charge
- D19 Seat Depth 19" (Short Base)..... No Charge
- D20 Seat Depth 20" (Long Base)..... No Charge
- D21 Seat Depth 21" (Long Base)..... No Charge
- D22 Seat Depth 22" (Long Base)..... No Charge

12 ADJASBA SEAT TO FLOOR HEIGHT

- LOW Seat To Floor 17 1/2"..... No Charge
- MED Seat To Floor 19 3/4"..... No Charge
- TALL Seat To Floor 21 1/2"..... \$252

VAN SEAT STYLE

- SEMIR Semi Recline..... No Charge
- FULLR Full Recline..... \$497

BACK SELECTION

ASBA BACK TYPE

- SECT Sectional Rear Uprights..... No Charge
- SRUP Straight Rear Uprights..... \$283
- 9720P No Push Handles (Straight Canes)..... \$111
- L8661 Omit Rear Uprights..... Less Net (\$25)

13 RECLINER BACK HEIGHT

- BHR18 Back Height 18"..... No Charge
- BHR20 Back Height 20"..... No Charge
- BHR22 Back Height 22"..... No Charge
- BHR24 Back Height 24"..... No Charge
- BHR26 Back Height 26"..... No Charge

BACK HEIGHT

ADJASBA

- BH12 Back Height 12"..... No Charge
- BH13 Back Height 13"..... No Charge
- BH14 Back Height 14"..... No Charge
- BH15 Back Height 15"..... No Charge
- BH16 Back Height 16"..... No Charge
- BH17 Back Height 17"..... No Charge
- BH18 Back Height 18"..... No Charge
- BH19 Back Height 19"..... No Charge
- BH20 Back Height 20"..... No Charge
- BH21 Back Height 21"..... No Charge
- BH22 Back Height 22"..... No Charge
- BH23 Back Height 23"..... No Charge
- BH24 Back Height 24"..... No Charge

CONTOURA

- TSBH25 Back Height 25"..... No Charge
- TSBH295 Back Height 29.5"..... No Charge

FOOTNOTES

- 10 Must select 12" - 18" Depth Range (1218D) when ordering the 12" - 16" Seat Width Range (1216W).
- 11 Manual reclining back not available.
- 12 Based on 5° tilt and 16" seat depth on ADJASBA add 1/4" for additional 2" of seat depth.
- 13 Used only with manual reclining seat selection.

ASBA BACK ANGLE DEGREE

<input type="checkbox"/>	85BA	85 Degree.....	\$66
<input type="checkbox"/>	90BA	90 Degree.....	No Charge
<input type="checkbox"/>	95BA	95 Degree.....	\$66
<input type="checkbox"/>	100BA	100 Degree.....	\$66
<input type="checkbox"/>	105BA	105 Degree.....	\$66

BACK ACCESSORIES

<input type="checkbox"/>	L8601	Omit Back Upholstery.....	Less Net (\$20)
<input type="checkbox"/>	RBP	Reinforced Back Upholstery.....	\$84

ARM SELECTION

14 ARM TYPE

15	<input type="checkbox"/>	18F	Full Length Fixed Hgt. Conv. Flip Up.....	\$244
15	<input type="checkbox"/>	28F	Desk Length Fixed Hgt. Conv. Flip Up.....	\$244
	<input type="checkbox"/>	19F	Full Length Adj. Hgt. Conv. Flip Up.....	\$348
	<input type="checkbox"/>	29F	Desk Length Adj. Hgt. Conv. Flip Up.....	\$348

MK6i SELECTION

CONTROLLER

<input type="checkbox"/>	MK690	90 Amp Controller.....	Standard	
8	<input type="checkbox"/>	MK6TT	True Track Controller.....	No Charge

STANDARD DRIVER CONTROL

<input type="checkbox"/>	SPJPLUS	SPJ+ Single Drive Joystick.....	Standard
--------------------------	---------	---------------------------------	-----------------

16 BASE ONLY SPJ+ OPTIONS

<input type="checkbox"/>	SPJAP	SPJ+ With Powered Seating Switch.....	Standard	
17	<input type="checkbox"/>	SPJAPM	Single Actuator Operation through SPJ+.....	\$1,535

MULTIPLE DRIVE CONTROL OPTIONS

<input type="checkbox"/>	MPJM6U	MPJ+ Multiple Joystick Drive Joystick.....	\$956	
18	<input type="checkbox"/>	PSFM6U	Personal Joystick, Inductive on Front.....	\$956
18	<input type="checkbox"/>	PSRM6U	Personal Joystick, Inductive on Rear.....	\$956
20	<input type="checkbox"/>	PBOD	Push Button On/Off w/PB Drive Select.....	\$84
19,20	<input type="checkbox"/>	PBSS	Push Button On/Off w/Speed Select.....	\$84
20	<input type="checkbox"/>	OR	On/Off Right.....	No Charge
20	<input type="checkbox"/>	OL	On/Off Left.....	No Charge
21	<input type="checkbox"/>	ASLRDYM6	MK6 Display Only (No Driver Control).....	\$226

FOR ALL YOUR ASL DRIVE SYSTEM NEEDS PLEASE REFER TO THE FOLLOWING ASL PRICE LIST FORM 05-220C

<input type="checkbox"/>	1500M6	Proportional RIM Head Control.....	\$2,263
<input type="checkbox"/>	SMHD	Small Headrest.....	\$225
<input type="checkbox"/>	LGHD	Large Headrest.....	\$225
<input type="checkbox"/>	1558MM6	Compact Joystick w/ Display.....	\$1,024

ALTERNATIVE DIGITAL CONTROL SYSTEMS

22	<input type="checkbox"/>	SNPM6	Sip-N-Puff.....	\$1,993
22	<input type="checkbox"/>	5018M6	Wafer Board w/SNP/DI.....	\$2,375
22	<input type="checkbox"/>	5020M6	Mini Tash Joystick & Digital Interface Box.....	\$2,506

FOOTNOTES

- 14 All arms are flip up and detachable.
- 15 Not available with w/seat widths 12"-15". Must use adjustable height.
- 16 Base only Systems for Multiple Actuator Power Seating Systems will require MPJ+, PSR+, PSD+, or MK6i Display & Alternative Driver Control.
- 17 Allows Actuator Operation through Driver Control. Includes Egg Switch for Attendant Control.
- 18 TOD-SS: Toggle On/Off-Drive Select, Speed Potentiometer Standard.
- 19 PSRM6 Only.
- 20 Drive select through mode switch, or mode switch and Joystick command.
- 21 Display only, allows choice from order form.
- 22 Must also use either a Multiple Drive Joystick or MK6 Display.
- 23 Requires SNPM6 Sip-N-Puff Interface Box.

SIP N PUFF KIT

<input type="checkbox"/>	PKG32666	Therapin Sip-N-Puff Breath Tube Kit.....	\$411
--------------------------	----------	--	-------

ADDITIONAL CONTROL CHOICES

<input type="checkbox"/>	1558M6	Compact Joystick - Less Std Mount.....	\$875	
23	<input type="checkbox"/>	5018	Wafer Board.....	\$382
23	<input type="checkbox"/>	5020	Mini Joystick w/Push.....	\$513
<input type="checkbox"/>	PACM6	Proportional Attendant Control.....	\$965	
23	<input type="checkbox"/>	1552M	Attendant Control.....	\$529

JOYSTICK MOUNTING HARDWARE

24	<input type="checkbox"/>	QLAM6	Quad Link Retract Joystick Mount MPJ+.....	\$334
24	<input type="checkbox"/>	QLA	Quad Link Retract Joystick Mount SPJ+.....	\$334
25	<input type="checkbox"/>	QLM	Quad Link Retract. Joystick Mount MPJ+ Van.....	\$334
25	<input type="checkbox"/>	QL	Quad Link Retract. Joystick Mount SPJ+ Van.....	\$334
	<input type="checkbox"/>	LEFT	Left Hand Mounted Joystick.....	No Charge
	<input type="checkbox"/>	RIMHW	Rim Hardware For 1500.....	\$229
	<input type="checkbox"/>	PKG32669	Complete Bib Assembly (Tash/1558).....	\$379
	<input type="checkbox"/>	RCM	Rear Cane Mount for Compact Joystick.....	\$149

<input type="checkbox"/>	CGAT400-427I	Gatlin Midline Mount Horizontal - Invacare MK6i for MPJM6	\$914
<input type="checkbox"/>	CGAT400-426D	Gatlin Midline Mount Horizontal - Invacare for SPJPLUS	\$914
<input type="checkbox"/>	CGAT400V-427I	Gat. Midline Mount Vertical (chin mount) - Invacare MMK6i for MPJM6	\$914
<input type="checkbox"/>	CGAT400V-426D	Gat.Mid. Mnt Vertical (chin mount) - Invacare for SPJPLUS	\$914
<input type="checkbox"/>	CGAT400-480R	Gatlin Midline Mount Horizontal - Invacare 1812	\$914
<input type="checkbox"/>	CGAT400V-480R	Gatlin Midline Mount Vertical (chin mount) - Invacare 1812	\$914

ELECTRONIC ACCESSORIES

	<input type="checkbox"/>	1813M6	MK6i Programmer w/Pro Memory Card.....	\$575	* UPDATED - 06/18/2007 - Increase in Price
	<input type="checkbox"/>	1813SD	Pro Memory Card.....	\$105	
	<input type="checkbox"/>	1560	"T" Handle Flexible Joystick Extension.....	\$86	
	<input type="checkbox"/>	1561	Straight Handle Flexible Joystick Extension.....	\$64	
	<input type="checkbox"/>	1826	Chin Cup.....	\$36	
	<input type="checkbox"/>	A24VPS	24V MK6i Auxiliary Power Source.....	\$235	
	<input type="checkbox"/>	AUX12M6	Auxillary Module For 1 & 2.....	\$659	
	<input type="checkbox"/>	AUX34M6	Auxillary Module For 3 & 4.....	\$659	
26	<input type="checkbox"/>	MPY	Mono Port "Y" Cable for MPJ+ or Display.....	\$20	

FRAME AND UPHOLSTERY COLOUR

FRAME FINISH

CONTEMPORARY

<input type="checkbox"/>	24P	Wet Black.....	No Charge
<input type="checkbox"/>	60P	Silver Metallic.....	No Charge
<input type="checkbox"/>	71P	Silver Vein.....	No Charge
<input type="checkbox"/>	105P	Midnight Blue.....	No Charge
<input type="checkbox"/>	115P	Black Prism.....	No Charge

TOO HOT

<input type="checkbox"/>	30P	Sunny Yellow.....	No Charge
<input type="checkbox"/>	61P	Electric Red.....	No Charge
<input type="checkbox"/>	119P	Bubblegum Pink.....	No Charge
<input type="checkbox"/>	120P	Grape Madness.....	No Charge

WAY COOL

<input type="checkbox"/>	62P	Electric Blue.....	No Charge
<input type="checkbox"/>	104P	Electric Teal.....	No Charge
<input type="checkbox"/>	121P	Lolly Pop Blue.....	No Charge
<input type="checkbox"/>	122P	Cosmic Blue.....	No Charge
<input type="checkbox"/>	125P	Emerald Green.....	No Charge
<input type="checkbox"/>	127P	Grasshopper Green.....	No Charge

UPHOLSTERY COLOURS

<input type="checkbox"/>	U240	Black Nylon Upholstery.....	Standard
<input type="checkbox"/>	U67	Black Vinyl.....	No Charge

VAN UPHOLSTERY

<input type="checkbox"/>	GC	Grey Cloth.....	No Charge
<input type="checkbox"/>	GV	Grey Vinyl.....	No Charge
<input type="checkbox"/>	TV	Tan Vinyl.....	No Charge

CASTER AND WHEEL OPTIONS

CASTER OPTIONS

	<input type="checkbox"/>	1259	6" X 2" Semi-Pneumatic Casters.....	\$122	
	<input type="checkbox"/>	1228-3	8" X 2 1/4" Semi Pneumatic. Casters Lt Grey.....	Standard	FOOTNOTES
	<input type="checkbox"/>	1251-3	8" X 2" Pneumatic Casters.....	No Charge	24 QLA and QLAM6 include ARM250 standard.
27	<input type="checkbox"/>	1220-3	9" X 2 3/4" Pneumatic Casters.....	\$151	25 Only available with Van Seating.
27	<input type="checkbox"/>	1465-3	9" X 2 3/4" Foam Filled Casters.....	\$276	26 For use with MPJ+ or MK6i. Display to provide additional switch site.

28 CASTER SUSPENSION

<input type="checkbox"/>	9213-6	Shock Fork For 6" Casters.....	\$195	
29	<input type="checkbox"/>	9213-8	Shock Fork For 8" Casters.....	\$195

WHEEL OPTION

<input type="checkbox"/>	1430-3	14" X 3" Pneumatic Tires.....	No Charge	
<input type="checkbox"/>	1431-3	14" X 3" Tire w/Foam Filled Inserts.....	\$184	
30,31	<input type="checkbox"/>	1440	14" X 4" Pneumatic Tires.....	\$231
30,31	<input type="checkbox"/>	1441	14" X 4" Tires w/Foam Filled Inserts.....	\$402

31 FREE WHEEL HUB

<input type="checkbox"/>	FHUB	Free Wheel Hub.....	\$172
--------------------------	------	---------------------	-------

REAR WHEEL LOCATION

<input type="checkbox"/>	RWF2	Rear Wheel Forward 2".....	\$65	
33	<input type="checkbox"/>	RWF1	Rear Wheel Forward 1".....	\$65
<input type="checkbox"/>	RWS	Rear Wheel Standard.....	Standard	

34 FENDER OPTION

<input type="checkbox"/>	FDR	Fender.....	\$61
--------------------------	-----	-------------	------

POSITIONING STRAPS

SEAT POSITIONING STRAPS

35	<input type="checkbox"/>	1311BK	Airline Buckle Seat Positioning Strap.....	\$79
<input type="checkbox"/>	1515	Auto Style Buckle Seat Positioning Strap.....	Standard	

CHEST POSITIONING STRAPS

<input type="checkbox"/>	7311BK	Airline Buckle Chest Positioning Strap.....	\$79
<input type="checkbox"/>	9515	Auto Style Chest Positioning Strap.....	\$79
<input type="checkbox"/>	7321BK	Hook & Loop Chest Positioning Strap.....	\$61

MISC ACCESSORIES

MISCELLANEOUS ACCESSORIES

36	<input type="checkbox"/>	1303	Wheel Locks.....	\$81
<input type="checkbox"/>	VANO2	O2 Holder - Van.....	\$236	
<input type="checkbox"/>	SEATO2	Oxygen Holder (ASBA Seat Only).....	\$236	
37	<input type="checkbox"/>	9205-3	Suspension Shocks.....	\$420

BATTERY SELECTION

38 INSTALLATION OF BATTERIES

<input type="checkbox"/>	22NFBATTERY	22 NF Battery (On Chair).....	\$534
<input type="checkbox"/>	24BATTERY	24 Gel Battery (On chair).....	\$654

FOOTNOTES

- 27 Not available with short base.
- 28 Not available with 9" casters.
- 29 Not available with TRRO or TRBKTS.
- 30 14" x 4" wheels increase overall width of the base by 2 1/2". Not available with TrueTrack package (TRUTK).
- 31 Does not include Storm Series® mag wheel.
- 32 Free wheel hub increases overall width 3/4". Not needed or available with TrueTrack package (TRUTK).
- 33 This location changes the center of gravity and performance. Forward wheel position will make turning and climbing easier, but more prone to "wheely".
- 34 Not available with TrueTrack (TRUTK).
- 35 1311BK not available with 12" - 15" width.
- 36 Not available with pneumatic tires.
- 37 Not available with 9" casters.
- 38 Batteries will be installed by Invacare Corporation at no additional charge. Chair will be shipped complete w/batteries.

BATTERY BOXES

<input type="checkbox"/>	22NF	22NF Battery.....	Standard
<input type="checkbox"/>	GP24	Group 24 Battery Boxes	\$351

BATTERY CHARGER OPTIONS

<input type="checkbox"/>	110CHARGER	110 Volt Battery Charger.....	No Charge
<input type="checkbox"/>	LESTR	Lester 8 Amp Battery Charger.....	\$332
<input type="checkbox"/>	L8900	Omit Charger.....	No Charge

HANGER BRACKET SELECTION

HANGER BRACKET TYPE

<input type="checkbox"/>	ASB	Centre Pivot Style Bracket.....	No Charge
<input type="checkbox"/>	PSB	Pin Style Bracket.....	No Charge
<input type="checkbox"/>	NHB	No Hanger Bracket.....	No Charge

FRONT RIGGING OPTIONS

FRONT RIGGING

CENTRE PIVOT STYLE

<input type="checkbox"/>	60HD	60° Swingaway Footrest.....	\$225
<input type="checkbox"/>	70NDH	70° Swingaway Footrest.....	\$225
<input type="checkbox"/>	90SWING	90° Swingaway Footrest.....	\$231
<input type="checkbox"/>	70STAPER	70° Tapered Footrest.....	\$262
<input type="checkbox"/>	AT5544	Elevating Legrests.....	\$529

PIN STYLE

<input type="checkbox"/>	PAS4A	Smart Artic Leg 17.25-21.5 (15" wide & less).....	\$366
<input type="checkbox"/>	P93M	Min Swingaway Footrest 7"-11".....	\$214
<input type="checkbox"/>	P904A	Swingaway Elevating Leg 16-19 3/4".....	\$333
³⁹ <input type="checkbox"/>	PAL4A	Smart Articulating Legrests 17.25-21.5" (16" wide plus).....	\$366

POWER LEGS

⁴⁰ <input type="checkbox"/>	ELRPW	Power ELR with Articulation (pair).....	\$2,856
⁴¹ <input type="checkbox"/>	4WSB	4 Way Switch Box for Power Legs.....	No Charge
⁴² <input type="checkbox"/>	S4WSB	Multiple Actuator Interface Box.....	\$1,973
<input type="checkbox"/>	FWT	4 Way Toggle.....	No Charge
<input type="checkbox"/>	QPB	4 Push Buttons.....	No Charge

FOOTPLATE OPTIONS

<input type="checkbox"/>	1651	Flip-Up Composite Footplates.....	No Charge
<input type="checkbox"/>	1350	Extra Large Aluminum Footplates.....	\$99
<input type="checkbox"/>	1450	Aluminum Footplates-Deluxe.....	No Charge
<input type="checkbox"/>	AT5543	Adjustable Angle Flip Up Footplates.....	\$122
<input type="checkbox"/>	AT5043	Flip Up Footplates.....	No Charge
<input type="checkbox"/>	S1263	One Piece Wooden Footboard.....	\$131

FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/>	LPT2	Longer Pivot Slide Tube For 1350.....	\$123
<input type="checkbox"/>	A1337	Calf Strap.....	\$59
<input type="checkbox"/>	IMPACTG	Impact Guards.....	No Charge
<input type="checkbox"/>	ALPT	Longer Pivot And Slide Tube.....	\$123
<input type="checkbox"/>	A9342	"H" Heel Strap.....	\$87
<input type="checkbox"/>	ALPT3	3" Extension For 90° Footrest	\$123
<input type="checkbox"/>	ALPT4	Longer Pivot And Slide Tube 4"	\$123

FOOTNOTES

- ³⁹ For seat widths greater than or equal to 16".
- ⁴⁰ Independent powered legrests include elevating and articulating legrests with angle adjustable footplates. Extension range of 13"-20". Requires center pivot hanger.
- ⁴¹ Required with all multiple actuator systems. Provides a D-9 port for the 4-way toggle or 4 push buttons.
- ⁴² Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control.

Torque SP Specifications		
* Overall Width	24.75" 27.25"	w/14"x3" wheel w/14"x4" wheel
**Overall Length	29.5" 32.5"	Short Base Long Base
***Seat-to-Floor	17 1/2" Standard 19 3/4" Optional	Low Medium
Battery	22NF deep cycle; gel/sealed (Not supplied)	
Weight w/o batteries 149 lb. (40 - 45) lb. each battery}		
Shipping Weight 266 lb. (18" X 18" seat w/Group 24 batteries)		
*Measurement based on wheel base. Seat width is width picked plus 7.5" for the armrest and joystick.		
**Measured from headtube to anti-tipper.		
***Based on 5° tilt and 16" seat depth; add 1/4" for additional 2" of seat depth		
Note: All specifications and dimensions are approximate.		

Manual Recline Specifications			
Seat Widths	14"-22" (1" increments)		
Seat Depths	16"-22" (1" increments) (Adjustable 16"-19", Fixed 20"-22")		
Back Heights	24" Standard, 18 1/2", 20", 22", 24", 26"		
Weight Limits	300 lb. with Standard Cylinders		
Headrest	10" Height		
Recline Range	90° to 170° (Approximate)		
Overall Back Height (w/o extension)	Low Seat	Med. Seat	High Seat
Standard Base	41 1/2"	43 1/2"	45 1/4"
(Measured with 5° Tilt and 10" for Extension)			

Weight Limit 300 Pounds

Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

NOTE: All weight limits are inclusive of add-on equipment and seat systems

- Chair Features and Options**
- Independent rear wheel suspension with new rubber element design
 - 6.5 mph standard speed range
 - Programmable performance to meet individuals needs
 - 17 1/2" Seat-to-Floor Standard, 19 3/4" and 21 1/2" Optional
 - Seat location system allows for balanced weight distribution
 - Continuous tilt seat 0° - 10°
 - Dynamic anti tippers
 - Rigid seat pan
 - Convenient rear battery access
 - 22NF batteries (standard) In a single 22NF battery box
 - Single battery connector
 - Arms are flip back/removable
 - Adjustable depth rigging attachment 2" (1" increments)
 - Automatic dual mode battery charger (Lester Brand)
 - Compact size: 24 3/4" overall base width w/14" x 3" wheels
 - Weight limit <300 lb.

Contoura Seat Sizes

	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X

*Not available with Contoura Seat Cushion

Back Heights 25" or 29.5"

Footrest Extension Range

	AT5043	AT5543	1651
60HD	N/A	12 1/4"-16 1/2"	13 1/2"-17 1/2"
SWING914	N/A	7 3/4"-12 1/2"	9 1/2"-14 1/2"
70MFX	6 3/4"-12 1/2"	4 1/2"-10"	N/A
70NHD	N/A	10 3/4"-15 3/4"	11 3/4"-17"
70STAPER	N/A	12 1/2"-16 1/2"	13 1/4"-17 1/4"
90SWING	4"-9"	5"-8"	N/A

Measured from Seat Pan to Back of Footplate

Memory Card - BASIC

- * Standard with ALL Expandable Driver Controls. (MPJ+, PSR+, PSF+, MK6i Display)
- * Used to Back Up/Restore Programmed Settings for Only One Chair.
- * May be given to user or kept with provider's files for safe keeping.
- * Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- * Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- * Standard with ALL MK6i Programmers.
- * Available also with USB Card Reader on Order form.
- * Contains Advanced Diagnostics, Help Screens, Software Updates, File Storage/Retrieval.
- * Used to back up multiple chairs and programming settings.
- * Intended for use by qualified providers only.
- * Not compatible with SPJ+ Joysticks

Extended Warranty Information:

Electrical Components: **2 Years**

Motors (2 and 4-pole): **2 Years**

GB Motors: **5 Years**

* Warranty valid to original purchaser who purchases this product within Canada

* For more product warranty information, please refer to the product owner's manual

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEAT AND BACK CUSHIONS
SEE PRICE LIST FORM # 05-338C.**

*** UPDATED - 06/18/2007 - New Warranty Box**

Invacare Canada reserves the right to change prices and specifications without notice.
 HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4G 4Z4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)
 Regional Centre - 16769 Boul. Hymus, Kirkland, P.Q. H9H 3L4 Phones: (514) 630-6080 (800) 363-0137(P.Q.) Fax: (514) 630-6120
 Invacare®, and Yes, you can., are trademarks of Invacare Corporation.
 WWW.INVACARE.CA