



NAME _____ PHONE NO _____

ACCOUNT# _____ P.O. _____ DATE _____

ADDRESS _____

BASE

CHAIR TYPE/POWER SEATING BASE

- 3GAR-CG Arrow for Formula CG with MK6i..... \$8,879
- 3GTQ-CG Torque SP for Formula CG with MK6i..... \$6,978

TRANSPORT TIE DOWN

- ¹ TRBKTS Wheelchair Transport Brackets..... \$325

USER WEIGHT LIMITS

ARROW

- U300 Under 300 lbs..... No Charge

² TORQUE SP

- U300 Under 300 lbs..... No Charge

DRIVE TYPE OTION

- ³ TRUTK True Track Package..... \$1,465

VENT TRAYS

- ⁴ COMP Compact Vent Tray..... \$1,141

FRAME TYPE

- ⁵ SFS Short Base..... No Charge
- LFS Long Base..... \$403

CASTER OPTIONS

- 1259 6" X 2" Semi-Pneumatic Casters..... \$122
- 1228-3 8" X 2 1/4" Semi-Pneumatic Casters Light Grey..... **Standard**
- 1251-3 8" X 2" Pneumatic Casters..... No Charge
- ⁶ 1220-3 9" X 2 3/4" Pneumatic Casters..... \$151
- ⁶ 1465-3 9" X 2 3/4" Foam Filled Casters..... \$276

⁷ CASTER SUSPENSION

- ⁸ 9213-6 Shock Fork For 6" Casters..... \$195
- 9213-8 Shock Fork For 8" Casters..... \$195

WHEEL OPTIONS

- 1430-3 14" X 3" Pneumatic Tires..... **Standard**
- 1431-3 14" X 3" Tire w/Foam Filled Inserts..... \$184
- ⁹ 1440 14" X 4" Pneumatic Tires..... \$231
- ⁹ 1441 14" X 4" Tires w/Foam Filled Inserts..... \$402

FREE WHEEL HUB

- ¹⁰ FHUB Free Wheel Hub..... \$172

FRAME FINISH

CONTEMPORARY

- 24P Wet Black..... No Charge
- 60P Silver Metallic..... No Charge
- 71P Silver Vein..... No Charge
- 105P Midnight Blue..... No Charge
- 115P Black Prism..... No Charge

***Transport Tie Down**

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

FOOTNOTES

- 1 Only Available with MED or TALL Seat to Floor Height.
- 2 Performance 4-pole motors, <300 lb. weight limit and 6.5 mph maximum speed. Standard on Torque SP.
- 3 TrueTrack includes Gearless Brushless GB™ motors which reaches 7 mph and has a 300 lb. weight capacity and Monroe® spring shock suspension. Must upgrade to MK6TT Controller. Available on Torque SP.
- 4 For ventilators with self-contained power supply and back-up. Vent tray on chair reduces chair weight capacity by approximately 50 lbs.
- 5 Not available with seat depths 19"-22". Must pick long base.
- 6 Not available with short base.
- 7 Not available with 9" casters.
- 8 Not available with SFS.
- 9 14 x 4" wheels increase the overall width of the base by 2.5". Does not include 3G design wheel. Not available with True Track.
- 10 Will add .75" overall width. Not available with Arrow GB or TrueTrack system. Does not include 3G designed wheel or fender.

TOO COOL

<input type="checkbox"/> 30P	Sunny Yellow.....	No Charge
<input type="checkbox"/> 61P	Electric Red.....	No Charge
<input type="checkbox"/> 119P	Bubblegum PinkV.....	No Charge
<input type="checkbox"/> 120P	Grape Madness.....	No Charge

WAY COOL

<input type="checkbox"/> 62P	Electric Blue.....	No Charge
<input type="checkbox"/> 104P	Electric Teal.....	No Charge
<input type="checkbox"/> 121P	Lolly Pop Blue.....	No Charge
<input type="checkbox"/> 122P	Cosmic Blue.....	No Charge
<input type="checkbox"/> 126P	Shamrock Green.....	No Charge
<input type="checkbox"/> 127P	Grasshopper Green.....	No Charge

SEAT

SYSTEM TYPE

<input type="checkbox"/> CGE	7" Elevating Seat - up to 300 lbs.....	\$2,147
------------------------------	--	---------

SEAT SELECTION

<input type="checkbox"/> ADJASBA	Rehab Seat w/Adj. Seat Tilt Back Angle.....	\$774
<input type="checkbox"/> BPO	Back Posts Only - Round Canes Only.....	No Charge
<input type="checkbox"/> CTS	Contoura Seating.....	\$553
<input type="checkbox"/> CURVABS	Curved ABS.....	Standard
<input type="checkbox"/> OMCU	Omit Cushion.....	Less Net (\$25)

SEAT WIDTH RANGE

<input type="checkbox"/> 1620W	16" to 20" Width (Adjustable).....	No Charge
<input type="checkbox"/> 2022W	20" to 22" Width (Adjustable).....	No Charge

10 SEAT WIDTH SETTING

<input type="checkbox"/> W16	16" Seat Width.....	No Charge
11 <input type="checkbox"/> W17	17" Seat Width (Not Available w/Contoura).....	No Charge
<input type="checkbox"/> W18	18" Seat Width.....	No Charge
11 <input type="checkbox"/> W19	19" Seat Width (Not Available w/Contoura).....	\$278
<input type="checkbox"/> W20	20" Seat Width.....	\$278
11 <input type="checkbox"/> W21	21" Seat Width ((Not Available w/Contoura).....	\$454
<input type="checkbox"/> W22	22" Seat Width.....	\$454

ADJUSTABLE SEAT DEPTH RANGE

<input type="checkbox"/> 1619D	16" to 19" Deep (Adjustable).....	No Charge
<input type="checkbox"/> 1922D	19" to 22" Deep (Adjustable).....	No Charge

SEAT DEPTH SETTING

13 <input type="checkbox"/> D16	16" Seat Depth.....	No Charge
13 <input type="checkbox"/> D17	17" Seat Depth.....	No Charge
<input type="checkbox"/> D18	18" Seat Depth.....	No Charge
12 <input type="checkbox"/> D19	19" Seat Depth.....	\$278
12 <input type="checkbox"/> D20	20" Seat Depth.....	\$278
12,14 <input type="checkbox"/> D21	21" Seat Depth.....	\$454
12,14 <input type="checkbox"/> D22	22" Seat Depth.....	\$454

15 SEAT TILT

<input type="checkbox"/> 0SA	0°.....	Standard
<input type="checkbox"/> 5SA	5°.....	No Charge

ASBA SEAT BACK HEIGHT

<input type="checkbox"/> SECT	Sectional Rear Uprights.....	No Charge
<input type="checkbox"/> SRUP	Straight Rear Uprights.....	\$283
<input type="checkbox"/> 9720P	No Push Handles (Straight Canes).....	\$111

BACK HEIGHT

CURVED ABS OR ADJUSTABLE ASBA

<input type="checkbox"/> BH16	Back Height 16".....	No Charge
<input type="checkbox"/> BH17	Back Height 17".....	No Charge
<input type="checkbox"/> BH18	Back Height 18".....	No Charge
<input type="checkbox"/> BH19	Back Height 19".....	No Charge
<input type="checkbox"/> BH20	Back Height 20".....	No Charge
<input type="checkbox"/> BH21	Back Height 21".....	No Charge
<input type="checkbox"/> BH22	Back Height 22".....	No Charge
<input type="checkbox"/> BH23	Back Height 23".....	No Charge
<input type="checkbox"/> BH24	Back Height 24".....	No Charge

FOOTNOTES

- 10 Contoura only available in even widths - non-adjustable in width. Please see chart on page 6.
- 11 Not available with Contoura.
- 12 Must pick long base.
- 13 Not available with Contoura seat cushion on 22" wide.
- 14 Not available with Contoura seat cushion on 16" or 18" wide.
- 15 Indicates static tilt of seat. Impacts seat-to-floor height, every 5 degrees raises the seat-to-floor height.

CONTOURA

- TSBH25 Back Height 25"..... No Charge
- TSBH295 Back Height 29.5" No Charge

ADJUSTABLE ARM TYPE (PAIR)

- 19F Full Length Adj. Hgt. Conv. Flip Up..... \$348
- 29F Desk Length Adj. Hgt. Conv. Flip Up..... \$348

MIXED ADJUSTABLE ARM TYPE (PAIR)

- 19FR Full Length Adj. Hgt. Conv. Flip Up Right..... \$174
- 29FR Desk Length Adj. Hgt. Conv. Flip Up Right..... \$174
- 19FL Full Length Adj. Hgt. Conv. Flip Up Left..... \$174
- 29FL Desk Length Adj. Hgt. Conv. Flip Up Left..... \$174

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap \$79
- 1515 Push Button Style Seat Positioning Strap..... Standard

CHEST POSITIONING STRAPS

- 7311BK Airline Buckle Chest Positioning Strap \$79
- 9515 Push Button Style Chest Positioning Strap \$79
- 7321BK Hook & Loop Chest Positioning Strap \$61

ELECTRONICS

CONTROLLER

TORQUE SP

- MK690ACC Controller w/Actuator Control Standard
- MK6TT True Track Controller..... No Charge

ARROW

- MK6TT True Track Controller..... Standard

STANDARD DRIVER CONTROL

- MPJM6 MPJ+ Multiple Drive Joystick..... Standard

ALTERNATIVE PROPORTIONAL CONTROL OPTIONS

- 16 PSFM6 Personal Joystick, Inductive on Front. No Charge
- 16 PSRM6 Personal Joystick, Inductive on Rear No Charge
- 18 PBOD Push Button. On/Off w/Pb Drive Select..... \$84
- 17,18 PBSS Push Button On/Off w/Speed Select..... \$84
- 18 OR On/Off Right..... No Charge
- 18 OL On/Off Left..... No Charge
- 19 ASLRDYM6 MK6 Display Only (**No Driver Control**)..... \$226

FOR ALL YOUR ASL DRIVE SYSTEM NEEDS PLEASE REFER TO THE FOLLOWING ASL PRICE LIST FORM 05-220C

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

- 1500M6 Proportional RIM Head Control..... \$2,263
- SMHD Small Headrest..... \$225
- LGHD Large Headrest..... \$225
- 1558MM6 Compact Joystick w/ Display..... \$1,024

ALTERNATIVE DIGITAL CONTROL SYSTEMS

- 20 SNPM6 Sip-N-Puff..... \$1,993
- 20 5018M6 Wafer Board w/SNP/DI..... \$2,375
- 20 5020M6 Mini Tash Joystick & Digital Interface Box..... \$2,506
- PKG32666 Therapin Sip-N-Puff Breath Tube Kit..... \$411

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount..... \$875
- 21 5018 Wafer Board..... \$382
- 21 5020 Mini Joystick w/Push..... \$513
- PACM6 Proportional Attendant Control..... \$956
- 21 1552M Attendant Control..... \$529

FOOTNOTES

- 16 TOD-SS: Toggle On/Off Drive Select, Speed, Speed Potentiometer Standard.
- 17 Drive Select through mode switch, or mode switch and Joystick Commands.
- 18 PSRM6 Only.
- 19 Display Only, No Driver Control. Allows choice from Price Sheet
- 20 Must also use either a Multiple Drive Joystick or MK6 Display.
- 21 Requires SNPM6 Interface Box.

JOYSTICK MOUNTING HARDWARE

22	<input type="checkbox"/> QLAM6	Quad Link Retract Joystick Mnt MPJ+.....	\$334
	<input type="checkbox"/> QLA	Quad Link Retract. Joystick Mnt SPJ+.....	\$334
	<input type="checkbox"/> LEFT	Left Hand Mounted Joystick.....	No Charge
	<input type="checkbox"/> RIMHW	Rim Hardware For 1500.....	\$229
	<input type="checkbox"/> PKG32669	Complete Bib Assembly (Tash/1558).....	\$379
	<input type="checkbox"/> RCM	Rear Cane Mount for Compact Joystick.....	\$149

<input type="checkbox"/> CGAT400-427I	Gatlin Midline Mount Horizontal - Invacare MK6i for MPJM6	\$914
<input type="checkbox"/> CGAT400-426D	Gatlin Midline Mount Horizontal - Invacare for SPJPLUS	\$914
<input type="checkbox"/> CGAT400V-427I	Gatlin Midline Mount Vertical (chin mount) - Invacare MMK6i for MPJM6	\$914
<input type="checkbox"/> CGAT400V-426D	Gatlin Midline Vertical (chin mount) - Invacare for SPJPLUS	\$914
<input type="checkbox"/> CGAT400-480R	Gatlin Midline Mount Horizontal - Invacare 1812	\$914
<input type="checkbox"/> CGAT400V-480R	Gatlin Midline Mount Vertical (chin mount) - Invacare 1812	\$914

ELECTRONIC ACCESSORIES

<input type="checkbox"/> 1813M6	MK6i Programmer w/Pro Memory Card.....	\$575	* UPDATED - 06/18/2007- Increase in Price
<input type="checkbox"/> 1813SD	Pro Memory Card.....	\$105	
<input type="checkbox"/> 1560	"T" Handle Flexible Joystick Extension.....	\$86	
<input type="checkbox"/> 1561	Straight Handle Flexible Joystick Extension.....	\$64	
<input type="checkbox"/> 1826	Chin Cup.....	\$36	
<input type="checkbox"/> A24VPS	24V MK6i Auxiliary Power Source.....	\$235	
<input type="checkbox"/> AUX12M6	Auxillary Module For 1 & 2.....	\$659	
<input type="checkbox"/> AUX34M6	Auxillary Module For 3 & 4.....	\$659	
23 <input type="checkbox"/> MPY	Mono Port "Y" Cable.....	\$20	

POWER SEATING ACCESSORIES

24 <input type="checkbox"/> SANODE	Single Actuator Operation through Driver Control.....	\$994
------------------------------------	---	-------

MISCELLANEOUS

BATTERY BOXES

25 <input type="checkbox"/> GP24	Group 24 Battery Boxes.....	\$351
----------------------------------	-----------------------------	-------

BATTERIES

27 <input type="checkbox"/> 22NFBATTERY	22 NF Battery (On Chair)	\$534
<input type="checkbox"/> 24BATTERY	24 Gel Battery (On Chair)	\$654

MISCELLANEOUS ACCESSORIES

28 <input type="checkbox"/> 1303	Wheel Locks.....	\$81
29 <input type="checkbox"/> 9205-3	Suspension Shocks.....	\$420
29 <input type="checkbox"/> 9505H-3	Suspension Shocks (Heavy Duty).....	\$420
29 <input type="checkbox"/> FDR	Fender.....	\$61
30 <input type="checkbox"/> EXTP	Extended Anti-Tippers.....	\$295

FRONT RIGGINGS

CENTRE PIVOT STYLE (ASB) (PAIR)

<input type="checkbox"/> 60HD	Centre Pivot Style Bracket.....	No Charge
<input type="checkbox"/> 70NDH	70° Swingaway Footrest.....	\$225
31 <input type="checkbox"/> 70STAPER	70° Tapered Footrest.....	\$262
<input type="checkbox"/> AT5544	Elevating Legrests.....	\$529

PIN STYLE (PSB) (PAIR)

<input type="checkbox"/> PW93	Swingaway Footrests 13 1/2" - 19".....	\$163
<input type="checkbox"/> P904A	Swingaway Elevating Leg 16-19 3/4".....	\$333
<input type="checkbox"/> PAL4A	Smart Artic Leg 17.25-21.5 (16" wide plus).....	\$366

FOOTPLATE OPTIONS (PAIR)

<input type="checkbox"/> 1651	Flip-Up Composite Footplates.....	No Charge
<input type="checkbox"/> 1350	Extra Large Aluminum Footplates.....	\$99
<input type="checkbox"/> AT5543	Adjustable Angle Flip Up Footplates.....	\$122

FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/> 1337	Calf Strap	\$59
<input type="checkbox"/> 1600BK	Heel Loops w/Ankle Straps	\$79
<input type="checkbox"/> LPT2	Longer Pivot Slide Tube For 1350.....	\$123
<input type="checkbox"/> IMPACTG	Impact Guards.....	No Charge
<input type="checkbox"/> ALPT	Longer Pivot And Slide Tube.....	\$123
<input type="checkbox"/> ALPT4	Longer Pivot And Slide Tube 4".....	\$123

FOOTNOTES

- 22 QLA includes ARM250 standard.
- 23 For use with MPJ+ or MK6i Dosplay to provide additional switch site.
- 24 Required for Elevate Operation through Driver Control - Multiple Drive Systems. Not compatible with SPJ+ Joysticks.
- 25 Only available on Torque SP (Standard on Arrow and Ranger X). Must select for TRBKTS.
- 26 Batteries will be installed by Invacare Corporation at additional charge. Chair will be shipped complete with batteries.
- 27 Not available with Arrow GB, Ranger X, or TrueTrack systems.
- 28 Not available with pneumatic tires.
- 29 Standard with Arrow GB or True Track Systems.
- 30 Comes standard with vent tray. Adds 2.5" to overall length.
- 31 Not available with 19" or 20" wide and 9" casters. Standard with impact guards on the footrests.

Battery Boxes	
*Standard Battery Tray	
Arrow	Group 24
Torque SP	22NF Standard Group 24 (optional)
*Sealed Gel Cell Batteries Recommended	

Batteries
22NF batteries with terminal configuration (positive on the left and negative on the right) **MUST** be used. 22NF batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection.

- Memory Card - BASIC**
- Standard with ALL Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).
 - Used to Back Up/Restore Programmed Settings for Only One Chair.
 - May be given to user or kept with provider's files for safe keeping.
 - Does not Contain Advanced Diagnostics/Help Screens/File Structure.
 - Not compatible with SPJ+ Joysticks.
- Memory Card - PROFESSIONAL**
- Standard with ALL MK6i Programmers.
 - Available also with USB Card Reader on Order Form.
 - Contains Advanced Diagnostics, Help Screens, Software Updates, File Storage/Retrieval.
 - Used to back up multiple chairs and programming settings.
 - Intended for use by qualified providers only.
 - Not compatible with SPJ+ Joysticks.

Curved ABS or Adjustable ASBA	
Seat Widths	16" - 22"
Seat Depths	16" - 22"
Back Heights	16" - 24"

**FOR YOUR SEATING NEEDS
 PLEASE REFER TO THE FOLLOWING
 REHAB SEAT AND BACK CUSHIONS
 SEE PRICE LIST FORM # 05-338C.**

Storm for Formula CG Elevating Seat Specifications		
	Arrow	Torque SP
Maximum Speed MPH		
TrueTrack	7	7 (option)
4-pole	N/A	N/A
4-pole (Performance)	N/A	6.5
4-pole HD	N/A	N/A
Seat Widths	17.75"	17.75"
Seat Depths	19.75"	19.75"
Weight Limits (Less 50 lb. for Vent)		
<300 lb.		X
<300 lb. True Track	X	X

- Formula CG Elevating Seat Features and Specifications**
- Up to 300 lb. weight capacity.
 - Traditional Rehab Seating (ASBA).
 - Seat-to-Floor height of 17.75" at 0° tilt.
 - Speed inhibit feature of 20% of maximum speed while elevated.

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X

*Not available with Contoura Seat Cushion
 Back Heights 25" or 29.5"

Extended Warranty Information:
 Electrical Components: **2 Years**
 Motors (2 and 4-pole): **2 Years**
 GB Motors: **5 Years**

* Warranty valid to original purchaser who purchases this product within Canada
 * For more product warranty information, please refer to the product owner's manual

*** UPDATED - 06/18/2007 - New Warranty Box**

Invacare Canada reserves the right to change prices and specifications without notice.
 HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4G 4Z4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)
 Regional Centre - 16769 Boul. Hymus, Kirkland, P.Q. H9H 3L4 Phones: (514) 630-6080 (800) 363-0137(P.Q.) Fax: (514) 630-6120
 Invacare®, and Yes, you can., are trademarks of Invacare Corporation.
 WWW.INVACARE.CA