



NAME \_\_\_\_\_ PHONE NO \_\_\_\_\_

ACCOUNT# \_\_\_\_\_ P.O. \_\_\_\_\_ DATE \_\_\_\_\_

ADDRESS \_\_\_\_\_

**BASE**

**CHAIR TYPE/POWER SEATING BASE**

- 3GAR-CG Arrow for Formula CG with MK6i..... \$8,879
- 3GRX-CG Ranger X for Formula CG with MK6i..... \$8,182
- 3GTQ-CG Torque SP for Formula CG with MK6i..... \$6,978

**\*Transport Tie Down**

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

\*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

**TRANSPORT TIE DOWN**

- 1  TRBKTS Wheelchair Transport Brackets..... \$325

**USER WEIGHT LIMITS**

**ARROW**

- U300 Under 300 lbs..... No Charge
- U400 Under 400 lbs..... No Charge

**RANGER X**

- U400 Under 400 lbs..... No Charge

**TORQUE SP**

- U300 Under 300 lbs..... No Charge

**DRIVE TYPE OPTION**

- 2  TRUTK True Track Upgrade..... \$1,465

**3 VENT TRAYS**

- VENT Articulating Vent Tray..... \$1,141

**FRAME TYPE**

- 4  SFS Short Base..... No Charge
- LFS Long Base..... \$403

**CASTER OPTIONS**

- 1259 6" X 2" Semi-Pneumatic Casters..... \$122
- 1228-3 8" X 2 1/4" Semi Pneumatic Casters Light Grey..... **Standard**
- 1251-3 8" X 2" Pneumatic Casters..... No Charge
- 5  1220-3 9" X 2 3/4" Pneumatic Casters..... \$151
- 5  1465-3 9" X 2 3/4" Foam Filled Casters..... \$276

**CASTER SUSPENSION**

- 6  9213-6 Shock Fork For 6" Casters..... \$195
- 9213-8 Shock Fork For 8" Casters..... \$195

**WHEEL OPTIONS**

- 1430-3 14" X 3" Pneumatic Tires..... **Standard**
- 1431-3 14" X 3" Tire w/Foam Filled Inserts..... \$184
- 7  1440 14" X 4" Pneumatic Tires..... \$231
- 7  1441 14" X 4" Tires w/Foam Filled Inserts..... \$402

**FREE WHEEL HUB**

- 8  FHUB Free Wheel Hub..... \$172

- FOOTNOTES**
- 1 Only Available with MED or TALL Seat to Floor Height.
  - 2 TrueTrack includes Gearless Brushless GB™ motors which reaches 7 mph and has a 400 lb. weight capacity and Monroe® spring shock suspension. Must upgrade to MK6TT Controller. This is Standard on 3GAR-CG and a pick for Ranger X or Torque SP.
  - 3 Vent tray on chair reduces chair weight capacity by approximately 50 lbs. Not available with TRBKTS.
  - 4 Not available with seat depths 19"-22". Must pick long base.
  - 5 Not available with short base.
  - 6 Not available with 9" casters or SFS.
  - 7 14 x 4" wheels increase the overall width of the base by 2.5". Does not include Storm Series design wheel. Not available with True Track.
  - 8 Will add .75" overall width. Not available with Arrow or TrueTrack system. Does not include 3G designed wheel or fender.

**FRAME FINISH**

**CONTEMPORARY**

- 24P Wet Black..... No Charge
- 60P Silver Metallic..... No Charge
- 71P Silver Vein..... No Charge
- 105P Midnight Blue..... No Charge
- 115P Black Prism..... No Charge

**TOO COOL**

- 30P Sunny Yellow..... No Charge
- 61P Electric Red..... No Charge
- 119P Bubblegum Pink..... No Charge
- 120P Grape Madness..... No Charge

**WAY COOL**

- 62P Electric Blue..... No Charge
- 104P Electric Teal..... No Charge
- 121P Lolly Pop Blue..... No Charge
- 122P Cosmic Blue..... No Charge
- 125P Emerald Green..... No Charge
- 127P Grasshopper Green..... No Charge

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X

\*Not available with Contoura Seat Cushion  
 Back Heights 25" or 29.5"

Conventional Seat Sizes	
Seat Widths	16" - 24"
Seat Depths	16" - 22"
Back Heights	20" - 26"

**SEAT**

**SYSTEM TYPE**

- CGR CG Recline Only System..... \$9,660

**SEAT-TO-FLOOR HEIGHT**

- 9 LOW 17.75"..... No Charge
- MED 18.75"..... No Charge
- TALL 19.75"..... \$252

**SEAT SELECTION**

- 10 CTS Contoura Seating..... \$553
- CNB Conventional Back..... Standard
- OMCU Omit Cushion..... Less Net (\$25)

**PUSH HANDLES**

- 11 PCNB Push Handles for Conventional and BPO..... No Charge

**SEAT WIDTH RANGE**

- 1620W 16" to 20" Width (Adjustable)..... No Charge
- 2024W 20" to 24" Width (Adjustable)..... \$429

**SEAT WIDTH SETTING**

- W16 16" Seat Width..... No Charge
- W17 17" Seat Width (Not Available w/Contoura)..... No Charge
- W18 18" Seat Width..... No Charge
- W19 19" Seat Width (Not Available w/Contoura)..... No Charge
- W20 20" Seat Width..... No Charge
- W21 21" Seat Width (Not Available w/Contoura)..... No Charge
- W22 22" Seat Width..... No Charge
- 12 W23 23" Seat Width (Not Available w/Contoura)..... No Charge
- 12 W24 24" Seat Width..... No Charge

**ADJUSTABLE SEAT DEPTH RANGE**

**CONVENTIONAL/CONTOURA**

- 1619D 16" to 19" Deep (Adjustable)..... No Charge
- 1922D 19" to 22" Deep (Adjustable)..... \$355

**CONVENTIONAL & CONTOURA**

**SEAT DEPTH SETTING**

- D16 16" Seat Depth..... No Charge
- D17 17" Seat Depth..... No Charge
- D18 18" Seat Depth..... No Charge
- D19 19" Seat Depth..... No Charge
- D20 20" Seat Depth..... No Charge
- D21 21" Seat Depth..... No Charge
- D22 22" Seat Depth..... No Charge

**FOOTNOTES**

- 9 Not available with Transport Tie Down option.
- 10 Selected Seat Depth on Contoura assumes 1" of back cushion compression. For example: When choosing an 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick includes leatherette seat cushion.
- 11 Push handles will add an additional 3.5" to back-cane height.
- 12 Available only with 400 lb. package.

**SEAT TILT**

- OSA 0°..... Standard
- 5SA 5°..... No Charge

**BACK HEIGHT**

**CONVENTIONAL BACK**

- TSBH20 Back Height 20"..... No Charge
- TSBH21 Back Height 21"..... No Charge
- TSBH22 Back Height 22"..... No Charge
- TSBH23 Back Height 23"..... No Charge
- TSBH24 Back Height 24"..... No Charge
- TSBH25 Back Height 25"..... \$311
- TSBH26 Back Height 26"..... \$311

**CONTOURA**

- TSBH25 Back Height 25"..... No Charge
- TSBH295 Back Height 29.5"..... No Charge

**HEADREST TYPE (OPTIONAL)**

- TAR127700 Curved Headrest ..... Standard
- TAR127700HW Curved Headrest Removable & Adjustable Hardware..... No Charge
- TAR128000 Two-Step Headrest ..... No Charge
- TAR128000HW Two-Step Headrest Removable & Adjustable Hardware..... No Charge
- T0022400 Large Headrest ..... \$97
- TAR22400HW Large Headrest Removable & Adjustable Hardware..... No Charge
- T0022600 Small Headrest ..... \$97
- TAR22600HW Small Headrest Removable & Adjustable Hardware..... No Charge

**FOR ALL YOUR STEALTH NEEDS PLEASE REFER TO THE FOLLOWING  
STEALTH PRICE LIST FORM 05-329C**

**ARM TYPE (PAIR)**

- RA19FL Reclining Full Length Adjustable Height - Left..... No Charge
- RA19FR Reclining Full Length Adjustable Height - Right..... No Charge
- RA29FL Reclining Desk Length Adjustable Height - Left..... No Charge
- RA29FR Reclining Desk Length Adjustable Height - Right..... No Charge

**13 OTTO BOCK ARM PADS (PAIR)**

- OB20400 One Piece Channel Armpads ..... \$165
- OB20500 Large Forearm Pads..... \$156
- OB20600 Medium Forearm Pads..... \$156
- OB20700 Small Forearm Pads..... \$156

**13 OTTO BOCK HAND PADS (PAIR)**

- OB23000 Large Flat Hand Pads..... \$128
- OB23300 Medium Flat Hand Pads..... \$128

**LATERALS (PAIR)**

**14 CONVENTIONAL BACK**

- TMSLM Swingaway Laterals Medium ..... \$315
- TMSLMHW S/A Laterals Medium Removable & Adjustable Hardware..... \$315
- TMSLL Swingaway Laterals Large ..... \$315
- TMSLLHW S/A Laterals Large Removable & Adjustable Hardware..... \$315

**15 CONTOURA SEATING**

- CSLM Swingaway Contoura Laterals Medium ..... \$315
- CSLMHW S/A Contoura Laterals Medium Removable & Adjustable Hardware..... \$315
- CSLL Swingaway Contoura Laterals Large ..... \$315
- CSLLHW S/A Contoura Laterals Large Removable & Adjustable Hardware..... \$315

**SEAT POSITIONING STRAPS**

- 1311BK Airline Buckle Seat Pos. Strap ..... \$79
- 1515 Push Button Style Seat Pos. Strap..... Standard

**FOOTNOTES**

- 13 Not compatible with quadlink retractable joystick mount.
- 14 For use with channel canes only.
- 15 For use with Contoura seating only.

**CHEST POSITIONING STRAPS**

<input type="checkbox"/> 7311BK	Airline Buckle Chest Pos. Strap .....	\$79
<input type="checkbox"/> 9515	Push Button Style Chest Pos. Strap .....	\$79
<input type="checkbox"/> 7321BK	Hook & Loop Chest Pos. Strap .....	\$61

**ELECTRONICS**

**CONTROLLER**

**TORQUE SP, RANGER X**

<input type="checkbox"/> MK690ACC	Controller w/Actuator Control .....	Standard
<input type="checkbox"/> MK6TT	True Track Controller.....	No Charge

**ARROW**

<input type="checkbox"/> MK6TT	True Track Controller.....	Standard
--------------------------------	----------------------------	----------

**STANDARD DRIVER CONTROLS**

<input type="checkbox"/> MPJM6	MPJ+ Multiple Drive Joystick.....	Standard
--------------------------------	-----------------------------------	----------

**ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS**

16 <input type="checkbox"/> PSFM6	Personal Joystick, Inductive on Front .....	\$956
16 <input type="checkbox"/> PSRM6	Personal Joystick, Inductive on Rear.....	\$956
18 <input type="checkbox"/> PBOD	Push Button On/Off w/PB Drive Select.....	\$84
17,18 <input type="checkbox"/> PBSS	Push Button On/Off w/Speed Select.....	\$84
18 <input type="checkbox"/> OR	On/Off Right.....	No Charge
18 <input type="checkbox"/> OL	On/Off Left.....	No Charge
19 ASLRDYM6	MK6 Display Only ( <b>No Driver Control</b> ).....	\$226

**FOR ALL YOUR ASL DRIVE SYSTEM NEEDS PLEASE REFER TO THE FOLLOWING ASL PRICE LIST FORM 05-220C**

<input type="checkbox"/> 1500M6	Proportional RIM Head Control.....	\$2,263
<input type="checkbox"/> SMHD	Small Headrest.....	\$225
<input type="checkbox"/> LGHD	Large Headrest.....	\$225
<input type="checkbox"/> 1558MM6	Compact Joystick w/ Display.....	\$1,024

**ALTERNATIVE DIGITAL CONTROL SYSTEMS**

20 <input type="checkbox"/> SNPM6	Sip-N-Puff.....	\$1,993
20 <input type="checkbox"/> 5018M6	Wafer Board w/SNP/DI.....	\$2,375
20 <input type="checkbox"/> 5020M6	Mini Tash Joystick & Digital Interface Box.....	\$2,506
<input type="checkbox"/> PKG32666	Therapin Sip-N-Puff Breath Tube Kit.....	\$411

**ADDITIONAL CONTROL CHOICES**

<input type="checkbox"/> 1558M6	Compact Joystick - Less Std Mount.....	\$875
21 <input type="checkbox"/> 5018	Wafer Board.....	\$382
21 <input type="checkbox"/> 5020	Mini Joystick w/Push.....	\$513
<input type="checkbox"/> PACM6	Proportional Attendant Control.....	\$956
21 <input type="checkbox"/> 1552M	Attendant Control.....	\$529

**JOYSTICK MOUNTING HARDWARE**

22 <input type="checkbox"/> QLAM6	Quad Link Retract Joystick Mount MPJ+.....	\$334
<input type="checkbox"/> LEFT	Left Hand Mounted Joystick.....	No Charge
<input type="checkbox"/> RIMHW	Rim Hardware For 1500.....	\$229
<input type="checkbox"/> PKG32669	Complete Bib Assembly (Tash/1558).....	\$379
<input type="checkbox"/> RCM	Rear Cane Mount for Compact Joystick.....	\$149

<input type="checkbox"/> CGAT400-427I	Gatlin Midline Mount Horizontal - Invacare MK6i for MPJM6	\$914
<input type="checkbox"/> CGAT400-426D	Gatlin Midline Mount Horizontal - Invacare for SPJPLUS	\$914
<input type="checkbox"/> CGAT400V-427I	Gatlin Midline Mount Vertical (chin mount) - Invacare MMK6i for MPJM6	\$914
<input type="checkbox"/> CGAT400V-426D	Gatlin Midline Mount Vertical (chin mount) - Invacare for SPJPLUS	\$914
<input type="checkbox"/> CGAT400-480R	Gatlin Midline Mount Horizontal - Invacare 1812	\$914
<input type="checkbox"/> CGAT400V-480R	Gatlin Midline Mount Vertical (chin mount) - Invacare 1812	\$914

**FOOTNOTES**

16 TOD-SS: Toggle On/Off Drive Select, Speed, Speed Potentiometer Standard.

**ELECTRONIC ACCESSORIES**

<input type="checkbox"/> 1813M6	MK6i Programmer w/Pro Memory Card.....	\$575	<b>* UPDATED - 06/18/2007- Increase in Price</b>
<input type="checkbox"/> 1813SD	Pro Memory Card.....	\$105	
<input type="checkbox"/> 1560	"T" Handle Flexible Joystick Extension.....	\$86	17 Drive Select through mode switch, or mode switch and Joystick Commands.
<input type="checkbox"/> 1561	Straight Handle Flexible Joystick Extension.....	\$64	18 PSRM6 Only.
<input type="checkbox"/> 1826	Chin Cup.....	\$36	19 Display Only, No Driver Control. Allows choice from Price Sheet
<input type="checkbox"/> A24VPS	24V MK6i Auxiliary Power Source.....	\$235	20 Must also use either a Multiple Drive Joystick or MK6 Display.
<input type="checkbox"/> AUX12M6	Auxiliary Module For 1 & 2.....	\$659	21 Requires SNPM6 Interface Box.
<input type="checkbox"/> AUX34M6	Auxiliary Module For 3 & 4.....	\$659	22 QLAM6 includes ARM250 standard.

NOTE: All specifications and dimensions are approximate.  
For the most current pricing information visit [www.invacare.ca](http://www.invacare.ca)

**POWER SEATING ACCESSORIES**

23	<input type="checkbox"/> SANODE	Single Actuator Operation through Driver Control.....	\$994
24	<input type="checkbox"/> 4WSB	4 Way Switch Box for Power Legs.....	No Charge
25	<input type="checkbox"/> S4WSB	Multiple Actuator Interface Box.....	\$1,973
	<input type="checkbox"/> FWT	4 Way Toggle.....	No Charge
	<input type="checkbox"/> QPB	4 Push Buttons.....	No Charge

**MISCELLANEOUS**

**BATTERY BOXES**

26	<input type="checkbox"/> GP24	Group 24 Battery Boxes.....	\$351
	<input type="checkbox"/> 22NF	22NF Battery.....	Standard

**27 BATTERIES**

28	<input type="checkbox"/> 22NFBATTERY	22 NF Battery (On Chair) .....	\$534
	<input type="checkbox"/> 24BATTERY	24 Gel Battery (On Chair) .....	\$654

**MISCELLANEOUS ACCESSORIES**

	<input type="checkbox"/> 1303	Wheel Locks.....	\$81
29	<input type="checkbox"/> 9205-3	Suspension Shocks.....	\$420
29	<input type="checkbox"/> 9205H-3	Suspension Shocks (Heavy Duty).....	\$420
29	<input type="checkbox"/> FDR	Fender.....	\$61
30	<input type="checkbox"/> EXTP	Extended Anti-Tippers.....	\$295

**FRONT RIGGINGS**

**HANGER BRACKET TYPE**

31	<input type="checkbox"/> ASB	Centre Pivot Style Bracket.....	No Charge
32	<input type="checkbox"/> PSB	Pin Style Bracket.....	No Charge

**CENTRE PIVOT STYLE (ASB) (PAIR)**

33	<input type="checkbox"/> 60HD	60° Swingaway Footrest.....	No Charge
33	<input type="checkbox"/> 70NDH	70° Non-Tapered Swingaway Footrest.....	\$225
33	<input type="checkbox"/> 70STAPER	70° Tapered Footrest.....	\$262
34	<input type="checkbox"/> AT5544	Elevating Legrests.....	\$529

**PIN STYLE (PSB) (PAIR)**

	<input type="checkbox"/> PW93	Swingaway Footrests 13 1/2" - 19".....	\$163
	<input type="checkbox"/> P904A	Swingaway Elevating Leg 16-19 3/4".....	\$333
	<input type="checkbox"/> PAL4A	Smart Artic Leg 17.25-21.5 (16" wide plus).....	\$366

**POWER ELEVATING LEGRESTS-ARTICULATING**

35,36	<input type="checkbox"/> ELRPW	Power ELR with Articulation (Aluminum).....	\$2,856
-------	--------------------------------	---------------------------------------------	---------

**FOOTPLATE OPTIONS (PAIR)**

	<input type="checkbox"/> 1651	Flip-Up Composite Footplates.....	No Charge
	<input type="checkbox"/> 1350	Extra Large Aluminum Footplates.....	\$99
	<input type="checkbox"/> AT5543	Adjustable Angle Flip Up Footplates.....	\$122

**FOOTREST/LEGREST ACCESSORIES**

	<input type="checkbox"/> 1337	Calf Strap .....	\$59
	<input type="checkbox"/> 1600BK	Heel Loops w/Ankle Straps .....	\$79
	<input type="checkbox"/> LPT2	Longer Pivot Slide Tube For 1350.....	\$123
	<input type="checkbox"/> IMPACTG	Impact Guards.....	No Charge
	<input type="checkbox"/> ALPT	Longer Pivot And Slide Tube.....	\$123
	<input type="checkbox"/> ALPT4	Longer Pivot And Slide Tube 4".....	\$123

**SEATING**

**FOR YOUR SEATING NEEDS  
PLEASE REFER TO THE FOLLOWING**

**REHAB SEAT AND BACK CUSHIONS**  
SEE PRICE LIST FORM # 05-338C.

**FOOTNOTES**

23 Required for Recline Operation through Driver Control - Single Actuator Systems only. The Egg Switch with 611 Mount does not come standard when ordering a SANODE.

24 Required with all multiple actuator systems. Provides a D-9 port for the 4 Way Toggle or 4 Push Buttons.

25 Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control.

26 Only available on Torque SP (Standard on Arrow and RangerX). Must select for TRBKTS.

27 Batteries will be installed by Invacare Corporation at additional charge. Chair will be shipped complete with batteries.

28 Not available with Arrow GB, Ranger X, or TrueTrack systems.

29 Standard with Arrow GB or True Track Systems.

30 Comes standard with vent tray. Adds 2.5" to overall length.

31 For Center Pivot Style Front Rigging Only.

32 For Pin & Saddle Front Rigging Only.

33 Not available with 19" or 20" wide and 9" casters. Standard with impact guards on the footrests.

34 Manually elevating legrests.

35 Independent powered elevating legrests include elevating and articulating legrests with adjustable angle footplates. Extension range is 13"-20". Sold as pair.

36 Must choose 4WSB or S4WSB to operate power legs.

Formula CG Recline Features and Specifications	
•	Weight capacity of 300 lb. or 400 lb. (with heavy-duty motors).
•	Recline range of 170°.
•	Mechanical Shear reduction of 3".
•	Conventional or Contoura style seating.
•	Adjustable seat frame with Conventional style.
•	Seat-to-Floor height of 17.75" at 0° tilt.

Tech Tips	
Operating Powered Seating through the Driver's Control. All chairs require adding or replacing one component to allow actuator operation through the driver control.	
<b>SPJ+ with PSS:</b>	Replace joystick with SPJ+ with ACC (Actuator Control) - E2310
<b>Single Actuator Systems with MPJ+, PSR+, PSF+ MK6i Display:</b>	Add a SANODE (Single Actuator Node) to the system - E2310
<b>Multiple Actuator Systems:</b>	Replace the standard 4Way Switch box with a Multiple Actuator Interface Box - E2311

Batteries	
<b>22NF</b> batteries with terminal configuration (positive on the left and negative on the right) MUST be used. 22NF batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection.	
<b>GP24</b> batteries with terminal configuration (negative on the left and positive on the right) MUST be used. GP24 batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection. See Owner's Manual, part number 1114809. These recommendations MUST be followed otherwise injury and damage may occur.	

Battery Boxes	
	*Standard Battery Tray
Arrow	Group 24
Ranger X	Group 24
Torque SP	22 NF Group 24 (Optional)
*Sealed Gel Cell Batteries Recommended	

Storm for Formula CG Recline Specifications			
	Arrow	Ranger X	Torque SP
<b>Maximum Speed MPH</b>			
TrueTrack	7	7 (option)	7 (option)
4-pole	N/A	N/A	NA/SP
4-pole Perf	N/A	N/A	6.5 /SP
4-pole HD	N/A	5	N/A
<b>Seat-To-Floor at 0°</b>			
Low	17.75"	17.75"	17.75"
High	19.75"	19.75"	19.75"
<b>Weight Limits</b>			
<300 lb.			X
<400 lb.		X	
<300 lb. True Track	X		X
<400 lb. True Track	X	X	

Memory Card - BASIC	
* Standard with ALL Expandable Driver Controls. (MPJ+, PSR+, PSF+, MK6i Display)	
* Used to Back Up/Restore Programmed Settings for Only One Chair.	
* May be given to user or kept with provider's files for safe keeping.	
* Does not Contain Advanced Diagnostics/Help Screens/File Structure.	
* Not compatible with SPJ+ Joysticks.	
<b>Memory Card - PROFESSIONAL</b>	
* Standard with ALL MK6i Programmers.	
* Available also with USB Card Reader on Order form.	
* Contains Advanced Diagnostics, Help Screens, Software Updates, File Storage/Retrieval.	
* Used to back up multiple chairs and programming settings.	
* Intended for use by qualified providers only.	
* Not compatible with SPJ+ Joysticks	

Extended Warranty Information:	
Electrical Components:	<b>2 Years</b>
Motors (2 and 4-pole):	<b>2 Years</b>
GB Motors:	<b>5 Years</b>
* Warranty valid to original purchaser who purchases this product within Canada	
* For more product warranty information, please refer to the product owner's manual	

**\* UPDATED - 06/18/2007 - New Warranty Box**

Invacare Canada reserves the right to change prices and specifications without notice.  
 HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4G 4Z4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)  
 Regional Centre - 16769 Boul. Hymus, Kirkland, P.Q. H9H 3L4 Phones: (514) 630-6080 (800) 363-0137(P.Q.) Fax: (514) 630-6120  
 Invacare®, and Yes, you can., are trademarks of Invacare Corporation.  
 WWW.INVACARE.CA