



NAME _____ PHONE NO _____
 ACCOUNT# _____ P.O. _____ DATE _____
 ADDRESS _____

BASE

CHAIR TYPE/POWER SEATING BASE

- 3GAR-CG Arrow for Formula CG with MK6i..... \$8,879
- 3GRX-CG Ranger X for Formula CG with MK6i..... \$8,182
- 3GTQ-CG Torque SP for Formula CG with MK6i..... \$6,978

TRANSPORT TIE DOWN

- 1 TRBKTS Wheelchair Transport Brackets..... \$325

USER WEIGHT LIMITS

ARROW

- 2 U300 Under 300 lbs..... No Charge
- 3 U400 Under 400 lbs..... No Charge

RANGER X

- 3 U400 Under 400 lbs..... No Charge

TORQUE SP

- 2 U300 Under 300 lbs..... No Charge

DRIVE TYPE OTION

- 4 TRUTK True Track upgrade..... \$1,465

5 VENT TRAYS

- VENT Articulating Vent Tray..... \$1,141

FRAME TYPE

- 6 SFS Short Base..... No Charge
- LFS Long Base..... \$403

CASTER OPTIONS

- 1259 6" X 2" Semi-Pneumatic Casters..... \$122
- 1228-3 8" X 2 1/4" Semi Pneu. Casters Lt Grey..... **Standard**
- 1251-3 8" X 2" Pneumatic Casters..... **No Charge**
- 7 1220-3 9" X 2 3/4" Pneumatic Casters..... \$151
- 7 1465-3 9" X 2 3/4" Foam Filled Casters..... \$276

8 CASTER SUSPENSION

- 9213-6 Shock Fork For 6" Casters..... \$195
- 9213-8 Shock Fork For 8" Casters..... \$195

WHEEL OPTIONS

- 1430-3 14" X 3" Pneumatic Tires..... **Standard**
- 1431-3 14" X 3" Tire w/Foam Filled Inserts..... \$184
- 9 1440 14" X 4" Pneumatic Tires..... \$231
- 9 1441 14" X 4" Tires w/Foam Filled Inserts..... \$402

10 FREE WHEEL HUB

- FHUB Free Wheel Hub..... \$172

***Transport Tie Down**

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

- FOOTNOTES**
- 1 Only Available with MED or TALL Seat to Floor Height.
 - 2 With 300 lb. weight capacity, maximum range of tilt is 55 degrees.
 - 3 With 400 lb. weight capacity, maximum range of tilt is 45 degrees.
 - 4 TrueTrack includes Gearless Brushless GB™ motors which reaches 7 mph and has a 400 lb. weight capacity and Monroe® spring shock suspension. Must upgrade to MK6TT Controller.
 - 5 Vent tray on chair reduces chair weight capacity by approximately 50 lbs.
 - 6 Not available with 9213-6 Shock Fork for 6" Casters or Seat depths 19"-22". Must pick long base.
 - 7 Not available with short base.
 - 8 Not available with 9" casters.
 - 9 14 x 4" wheels increase the overall width of the base by 2.5". Does not include 3G design wheel. Not available with True Track.
 - 10 Will add .75" overall width. Not available with Arrow GB or TrueTrack system. Does not include 3G designed wheel or fender.

FRAME FINISH

CONTEMPORARY

- 24P Wet Black..... No Charge
- 60P Silver Metallic..... No Charge
- 71P Silver Vein..... No Charge
- 105P Midnight Blue..... No Charge
- 115P Black Prism..... No Charge

TOO COOL

- 30P Sunny Yellow..... No Charge
- 61P Electric Red..... No Charge
- 119P Bubblegum PinkV..... No Charge
- 120P Grape Madness..... No Charge

WAY COOL

- 62P Electric Blue..... No Charge
- 104P Electric Teal..... No Charge
- 121P Lolly Pop Blue..... No Charge
- 122P Cosmic Blue..... No Charge
- 125P Emerald Green..... No Charge
- 127P Grasshopper Green..... No Charge

SEAT

SYSTEM TYPE

- CGTR Tilt & Recl. w/Mechanical Shear Reduction..... \$11,541

SEAT-TO-FLOOR HEIGHT

- ¹¹ LOW 17.75"..... No Charge
- MED 18.75"..... No Charge
- TALL 19.75"..... \$252

SEAT SELECTION

- ¹² CTS Contoura Seating..... \$553
- CNB Conventional Back..... Standard
- OMCU Omit Cushion..... Less Net (\$25)

PUSH HANDLES

- ¹³ PCNB Push Handles for Conventional and BPO..... No Charge

SEAT WIDTH RANGE

- 1620W 16" to 20" Width (Adjustable)..... No Charge
- 2024W 20" to 24" Width (Adjustable)..... \$429

SEAT WIDTH SETTING

- W16 16" Seat Width..... No Charge
- W17 17" Seat Width (Not Available w/Contoura)..... No Charge
- W18 18" Seat Width..... No Charge
- W19 19" Seat Width (Not Available w/Contoura)..... No Charge
- W20 20" Seat Width..... No Charge
- W21 21" Seat Width ((Not Available w/Contoura)..... No Charge
- W22 22" Seat Width..... No Charge
- W23 23" Seat Width (Not Available w/Contoura)..... No Charge
- W24 24" Seat Width..... No Charge

ADJUSTABLE SEAT DEPTH RANGE

CONVENTIONAL/CONTOURA

- 1619D 16" to 19" Deep (Adjustable)..... No Charge
- 1922D 19" to 22" Deep (Adjustable)..... \$355

CONVENTIONAL & CONTOURA

SEAT DEPTH SETTING

- D16 16" Seat Depth..... No Charge
- D17 17" Seat Depth..... No Charge
- D18 18" Seat Depth..... No Charge
- D19 19" Seat Depth..... \$278
- D20 20" Seat Depth..... \$278
- D21 21" Seat Depth..... \$454
- D22 22" Seat Depth..... \$454

Conventional Seat Sizes	
Seat Widths	16" - 24"
Seat Depths	16" - 22"
Back Heights	20" - 26"

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X

*Not available with Contoura Seat Cushion	
Back Heights	25" or 29.5"

FOOTNOTES

11 Not available with Transport Tie Down option.

12 Selected Seat Depth on Contoura assumes 1" of back cushion compression. For example: When choosing an 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick includes leatherette seat cushion.

13 Push handles will add an additional 3.5" to back-cane height.

SEAT TILT

- OSA 0°..... Standard
- 5SA 5°..... No Charge

BACK HEIGHT

CONVENTIONAL BACK

- TSBH20 Back Height 20"..... No Charge
- TSBH21 Back Height 21"..... No Charge
- TSBH22 Back Height 22"..... No Charge
- TSBH23 Back Height 23"..... No Charge
- TSBH24 Back Height 24"..... No Charge
- TSBH25 Back Height 25"..... \$311
- TSBH26 Back Height 26"..... \$311

CONTOURA

- TSBH25 Back Height 25"..... No Charge
- TSBH295 Back Height 29.5"..... No Charge

HEADREST TYPE (OPTIONAL)

- TAR127700 Curved Headrest Standard
- TAR127700HW Curved Headrest Removable & Adjustable Hardware..... No Charge
- TAR128000 Two-Step Headrest No Charge
- TAR128000HW Two-Step Headrest Removable & Adjustable Hardware..... No Charge
- T0022400 Large Headrest \$97
- TAR22400HW Large Headrest Removable & Adjustable Hardware..... No Charge
- T0022600 Small Headrest \$97
- TAR22600HW Small Headrest Removable & Adjustable Hardware..... No Charge

**FOR ALL YOUR STEALTH NEEDS PLEASE REFER TO THE FOLLOWING
STEALTH PRICE LIST FORM 05-329C**

ARM TYPE (PAIR)

- RA19FL Reclining Full Length Adj. Hgt. - Left..... No Charge
- RA19FR Reclining Full Length Adj. Hgt. - Right..... No Charge
- RA29FL Reclining Desk Length Adj. Hgt. - Left..... No Charge
- RA29FR Reclining Desk Length Adj. Hgt. - Right..... No Charge

14 OTTO BOCK ARM PADS (PAIR)

- OB20400 One Piece Channel Armpads \$165
- OB20500 Large Forearm Pads..... \$156
- OB20600 Medium Forearm Pads..... \$156
- OB20700 Small Forearm Pads..... \$156

14 OTTO BOCK HAND PADS (PAIR)

- OB23000 Large Flat Hand Pads..... \$128
- OB23300 Medium Flat Hand Pads..... \$128

LATERALS (PAIR)

15 CONVENTIONAL BACK

- TMSLM Swingaway Laterals Medium \$315
- TMSLMHW S/A Laterals Medium Removable & Adjustable Hardware..... \$315
- TMSLL Swingaway Laterals Large \$315
- TMSLLHW S/A Laterals Large Removable & Adjustable Hardware..... \$315

16 CONTOURA SEATING

- CSLM Swingaway Contoura Laterals Medium \$315
- CSLMHW S/A Contoura Laterals Medium Removable & Adjustable Hardware..... \$315
- CSLL Swingaway Contoura Laterals Large \$315
- CSLLHW S/A Contoura Laterals Large Removable & Adjustable Hardware..... \$315

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Pos. Strap \$79
- 1515 Push Button Style Seat Pos. Strap..... Standard

FOOTNOTES

- 14 Not compatible with quadlink retractable joystick mount.
- 15 For use with channel canes only.
- 16 For use with Contoura seating only.

CHEST POSITIONING STRAPS

<input type="checkbox"/> 7311BK	Airline Buckle Chest Pos. Strap	\$79
<input type="checkbox"/> 9515	Push Button Style Chest Pos. Strap	\$79
<input type="checkbox"/> 7321BK	Hook & Loop Chest Pos. Strap	\$61

ELECTRONICS

CONTROLLER

TORQUE SP, RANGER X

<input type="checkbox"/> MK690ACC	Controller w/Actuator Control	Standard
<input type="checkbox"/> MK6TT	True Track Controller.....	No Charge

ARROW

<input type="checkbox"/> MK6TT	True Track Controller.....	Standard
--------------------------------	----------------------------	----------

STANDARD DRIVER CONTROL

<input type="checkbox"/> MPJM6	MPJ+ Multiple Drive Joystick.....	Standard
--------------------------------	-----------------------------------	----------

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

17 <input type="checkbox"/> PSFM6	Personal Joystick, Inductive on Front.....	No Charge
17 <input type="checkbox"/> PSRM6	Personal Joystick, Inductive on Rear.....	No Charge
19 <input type="checkbox"/> PBOD	Push Button On/Off w/PB Drive Select.....	\$84
18,19 <input type="checkbox"/> PBSS	Push Button On/Off w/Speed Select.....	\$84
19 <input type="checkbox"/> OR	On/Off Right.....	No Charge
19 <input type="checkbox"/> OL	On/Off Left.....	No Charge
20 <input type="checkbox"/> ASLRDYM6	MK6 Display Only (No Driver Control).....	\$226

FOR ALL YOUR ASL DRIVE SYSTEM NEEDS PLEASE REFER TO THE FOLLOWING ASL PRICE LIST FORM 05-220C

<input type="checkbox"/> 1500M6	Proportional RIM Head Control.....	\$2,263
<input type="checkbox"/> SMHD	Small Headrest.....	\$225
<input type="checkbox"/> LGHD	Large Headrest.....	\$225
<input type="checkbox"/> 1558MM6	Compact Joystick w/ Display.....	\$1,024

ALTERNATIVE DIGITAL CONTROL SYSTEMS

21 <input type="checkbox"/> SNPM6	Sip-N-Puff.....	\$1,993
21 <input type="checkbox"/> 5018M6	Wafer Board w/SNP/DI.....	\$2,375
21 <input type="checkbox"/> 5020M6	Mini Tash Joystick & Digital Interface Box.....	\$2,506
<input type="checkbox"/> PKG32666	Therapin Sip-N-Puff Breath Tube Kit.....	\$411

ADDITIONAL CONTROL CHOICES

<input type="checkbox"/> 1558M6	Compact Joystick - Less Std Mount.....	\$875
22 <input type="checkbox"/> 5018	Wafer Board.....	\$382
22 <input type="checkbox"/> 5020	Mini Joystick w/Push.....	\$513
<input type="checkbox"/> PACM6	Proportional Attendant Control.....	\$956
22 <input type="checkbox"/> 1552M	Attendant Control.....	\$529

JOYSTICK MOUNTING HARDWARE

23 <input type="checkbox"/> QLAM6	Quad Link Retract Joystick Mnt MPJ+.....	\$334
<input type="checkbox"/> LEFT	Left Hand Mounted Joystick.....	No Charge
<input type="checkbox"/> RIMHW	Rim Hardware For 1500.....	\$229
<input type="checkbox"/> PKG32669	Complete Bib Assembly (Tash/1558).....	\$379
<input type="checkbox"/> RCM	Rear Cane Mount for Compact Joystick.....	\$149

FOOTNOTES

<input type="checkbox"/> CGAT400-427I	Gatlin Midline Mount Horizontal - Invacare MK6i for MPJM6	\$914	17 TOD-SS: Toggle On/Off Drive Select, Speed, Speed Potentiometer Standard.
<input type="checkbox"/> CGAT400-426D	Gatlin Midline Mount Horizontal - Invacare for SPJPLUS	\$914	18 Drive Select through mode switch, or mode switch and Joystick Commands.
<input type="checkbox"/> CGAT400V-427I	Gatlin Midline Mount Vertical (chin mount) - Invacare MMK6i for MPJM6	\$914	19 PSRM6 Only.
<input type="checkbox"/> CGAT400V-426D	Gatlin Midline Mount Vertical (chin mount) - Invacare for SPJPLUS	\$914	20 Display Only, No Driver Control. Allows choice from Price Sheet
<input type="checkbox"/> CGAT400-480R	Gatlin Midline Mount Horizontal - Invacare 1812	\$914	21 Must also use either a Multiple Drive Joystick or MK6 Display.
<input type="checkbox"/> CGAT400V-480R	Gatlin Midline Mount Vertical (chin mount) - Invacare 1812	\$914	22 Requires SNPM6 Interface Box.

ELECTRONIC ACCESSORIES

<input type="checkbox"/>	1813M6	MK6i Programmer w/Pro Memory Card.....	\$575
<input type="checkbox"/>	1813SD	Pro Memory Card.....	\$105
<input type="checkbox"/>	1560	"T" Handle Flexible Joystick Extension.....	\$86
<input type="checkbox"/>	1561	Straight Handle Flexible Joystick Extension.....	\$64
<input type="checkbox"/>	1826	Chin Cup.....	\$36
<input type="checkbox"/>	A24VPS	24V MK6i Auxiliary Power Source.....	\$235
<input type="checkbox"/>	AUX12M6	Auxiliary Module For 1 & 2.....	\$659
<input type="checkbox"/>	AUX34M6	Auxiliary Module For 3 & 4.....	\$659
24 <input type="checkbox"/>	MPY	Mono Port "Y" Cable.....	\$20

* **UPDATED - 06/18/2007- Increase in Price**

Battery Boxes	
*Standard Battery Tray	
TDX5	Group 24
TDX4	Group 24 22 NF (optional)
TDX3	22 NF Group 24 (Optional)
*Sealed Gel Cell Batteries Recommended	

POWER SEATING ACCESSORIES

25 <input type="checkbox"/>	4WSB	4 Way Switch Box for Power Legs.....	No Charge
26 <input type="checkbox"/>	S4WSB	Multiple Actuator Interface Box.....	\$1,973
<input type="checkbox"/>	FWT	4 Way Toggle.....	No Charge
<input type="checkbox"/>	QPB	4 Push Buttons.....	No Charge

Batteries
22NF batteries with terminal configuration (positive on the left and negative on the right) **MUST** be used. 22NF batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection.

GP24 batteries with terminal configuration (negative on the left and positive on the right) **MUST** be used. GP24 batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection. See Owner's Manual, part number 1114809. These recommendations **MUST** be followed otherwise injury and damage may occur.

MISCELLANEOUS

BATTERY BOXES

27 <input type="checkbox"/>	GP24	Group 24 Battery Boxes.....	\$351
-----------------------------	------	-----------------------------	-------

BATTERIES

29 <input type="checkbox"/>	22NFBATTERY	22 NF Battery (On Chair)	\$534
<input type="checkbox"/>	24BATTERY	24 Gel Battery (On Chair)	\$654

MISCELLANEOUS ACCESSORIES

30 <input type="checkbox"/>	1303	Wheel Locks.....	\$81
31 <input type="checkbox"/>	9205-3	Suspension Shocks.....	\$420
31 <input type="checkbox"/>	9505H-3	Suspension Shocks (Heavy Duty).....	\$420
31 <input type="checkbox"/>	FDR	Fender.....	\$61
32 <input type="checkbox"/>	EXTP	Extended Anti-Tippers.....	\$295

FRONT RIGGINGS

HANGER BRACKET TYPE

33 <input type="checkbox"/>	ASB	Centre Pivot Style Bracket.....	No Charge
34 <input type="checkbox"/>	PSB	Pin Style Bracket.....	No Charge

CENTRE PIVOT STYLE (ASB) (PAIR)

35 <input type="checkbox"/>	60HD	Centre Pivot Style Bracket.....	No Charge
35 <input type="checkbox"/>	70NDH	70° Swingaway Footrest.....	\$225
35 <input type="checkbox"/>	70STAPER	70° Tapered Footrest.....	\$262
36 <input type="checkbox"/>	AT5544	Elevating Legrests.....	\$529

PIN STYLE (PSB) (PAIR)

<input type="checkbox"/>	PW93	Swingaway Footrests 13 1/2" - 19".....	\$163
36 <input type="checkbox"/>	P904A	Swingaway Elevating Leg 16-19 3/4".....	\$333
36 <input type="checkbox"/>	PAL4A	Smart Artic Leg 17.25-21.5 (16" wide plus).....	\$366

FOOTNOTES

- 23 QLAM6 includes ARM250 standard.
- 24 For use with MPJ+ or MK6i Display to provide additional switch site.
- 25 Required with all multiple actuator systems. Provides a D-9 port for the 4 Way Toggle or 4 Push Buttons.
- 26 Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control.
- 27 Only available on Torque SP (Standard on Arrow and Ranger X). Must select for TRBKTS.
- 28 Batteries will be installed by Invacare Corporation at additional charge. Chair will be shipped complete with batteries.
- 29 Not available with Arrow GB, Ranger X, or TrueTrack systems.
- 30 Not available with pneumatic tires.
- 31 Standard with Arrow GB or True Track Systems.
- 32 Comes standard with vent tray. Adds 2.5" to overall length.
- 33 For Center Pivot Style Front Rigging Only.
- 34 For Pin & Saddle Front Rigging Only.
- 35 Not available with 19" or 20" wide and 9" casters. Standard with impact guards on the footrests.
- 36 Manually elevating legrests.

POWER LEGS

<input type="checkbox"/>	ELRPW	Power ELR with Articulation (Aluminum).....	\$2,856
--------------------------	-------	---	---------

FOOTPLATE OPTIONS (PAIR)

<input type="checkbox"/>	1651	Flip-Up Composite Footplates.....	No Charge
<input type="checkbox"/>	1350	Extra Large Aluminum Footplates.....	\$99
<input type="checkbox"/>	AT5543	Adjustable Angle Flip Up Footplates.....	\$122

FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/>	1337	Calf Strap	\$59
<input type="checkbox"/>	1600BK	Heel Loops w/Ankle Straps	\$79
<input type="checkbox"/>	LPT2	Longer Pivot Slide Tube For 1350.....	\$123
<input type="checkbox"/>	IMPACTG	Impact Guards.....	No Charge
<input type="checkbox"/>	ALPT	Longer Pivot And Slide Tube.....	\$123
<input type="checkbox"/>	ALPT4	Longer Pivot And Slide Tube 4".....	\$123

Storm for Formula CG Tilt/Recline Specifications			
	Arrow	Ranger X	Torque SP
Maximum Speed MPH			
TrueTrack	7	7	7
4-pole	N/A	N/A	4.5
4-pole HD	N/A	N/A	4.5
Perf 4-pole	N/A	6	6
Seat-To-Floor at 0°			
Low	17.75"	17.75"	17.75"
High	19.75"	19.75"	19.75"
Weight Limits			
<300 lb. (250 lb. w/vent)	X	X	X
<400 lb. (350 lb. w/vent)	N/A	N/A	X
<400 lb. True Track (350 lb. w/vent)	X	X	X
Warranty	Lifetime Limited on Base Frame Three (3) Year Limited on Seat Frame One (1) Year Limited on Electronics One (1) Year Limited on Actuators Five (5) Year Limited Warranty on GB™ Motors Eighteen (18) Months Limited Warranty on Gear Motor		

- Formula CG Tilt/Recline Features and Specifications**
- Weight capacity of 300 lb. or 400 lb. (with heavy-duty motors).
 - Recline range of 170°
 - Mechanical Shear reduction of 3".
 - Conventional or Contoura style seating.
 - Adjustable seat frame with Conventional style.
 - Seat-to-Floor height of 17.75" at 0° tilt.

**FOR YOUR SEATING NEEDS
 PLEASE REFER TO THE FOLLOWING
 REHAB SEAT AND BACK CUSHIONS
 SEE PRICE LIST FORM # 05-338C.**

Extended Warranty Information:
 Electrical Components: **2 Years**
 Motors (2 and 4-pole): **2 Years**
 GB Motors: **5 Years**

* Warranty valid to original purchaser who purchases this product within Canada
 * For more product warranty information, please refer to the product owner's manual

*** UPDATED - 06/18/2007 - New Warranty Box**

Invacare Canada reserves the right to change prices and specifications without notice.
 HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4G 4Z4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)
 Regional Centre - 16769 Boul. Hymus, Kirkland, P.Q. H9H 3L4 Phones: (514) 630-6080 (800) 363-0137(P.Q.) Fax: (514) 630-6120
 Invacare®, and Yes, you can., are trademarks of Invacare Corporation.
 WWW.INVACARE.CA