



NAME \_\_\_\_\_ PHONE NO \_\_\_\_\_  
 ACCOUNT# \_\_\_\_\_ P.O. \_\_\_\_\_ DATE \_\_\_\_\_  
 ADDRESS \_\_\_\_\_

**FRAME**

**CHAIR TYPE/POWER SEATING BASE**

- TDX5-PS TDX5 for Formula w/MK6i..... \$8,974
- TDX3-PS TDX3 for Formula w/MK6i...(includes Std MPJM6)..... \$8,164

**TRANSPORT TIE DOWN**

- TRBKTS Wheelchair Transport Brackets..... \$325

**SYSTEM TYPE**

- E-TDX Formula Elevate for TDX..... \$2,716

**USER WEIGHT LIMITS**

**TDX5**

- U300 User Weight <300lbs. (<250 lbs. w/Vent)..... No Charge

**TDX3**

- U300 User Weight <300lbs. (No Vent)..... No Charge

**2 TRUETRACK PACKAGE FOR TDX3 AND TDX4**

- TRUTK TrueTrack Upgrade..... \$1,465

**WIDTH AND DEPTH**

**4 ADJUSTABLE ASBA SEAT WIDTH RANGE**

- 1620W 16" to 20" Width (Adjustable)..... No Charge
- 2024W 20" to 24" Width (Adjustable)..... \$429

**ADJASBA AND CONTOURA SEAT WIDTH SETTING**

- W16 16" Seat Width..... No Charge
- W17 17" Seat Width..... No Charge
- W18 18" Seat Width..... No Charge
- W19 19" Seat Width..... No Charge
- W20 20" Seat Width..... No Charge
- W21 21" Seat Width..... No Charge
- W22 22" Seat Width..... No Charge

**ADJUSTABLE SEAT DEPTH RANGE**

**ADJASBA**

- 1619D 16" to 19" Deep (Adjustable)..... No Charge
- 1922D 19" to 22" Deep (Adjustable)..... \$355

**CONTOURA**

- 1618D 16" to 18" Deep (Adjustable)..... No Charge
- 1921D 19" to 21" Deep (Adjustable)..... No Charge
- 22D 22" Deep..... No Charge

**\*Transport Tie Down**

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

\*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

**\*NOTE\***

**If you need TDX4-PS Formula Elevate with 4-Pole High Performance Motors, please use special quote number C-070118031949**  
**No Charge**

**FOOTNOTES**

- 1 Package includes GB Motors and TrueTrack driving technology, which enables driver to track straighter regardless of the terrain. Top Speed is 7.5 MPH. Must select MK6TT controller.
- 2 It is recommended to upgrade to Group 24 batteries and/or TrueTrack for improved battery charge life.
- 3 Package includes 4-pole motors and 6.5 MPH top speed unless TRUTK upgrade is selected.
- 4 Contoura seating available in even widths 16" - 22" and depths between 16" - 22" in 1" increments.

**4 ADJASBA SEAT DEPTH SETTING**

<input type="checkbox"/> D16	16" Seat Depth.....	No Charge
<input type="checkbox"/> D17	17" Seat Depth.....	No Charge
<input type="checkbox"/> D18	18" Seat Depth.....	No Charge
<input type="checkbox"/> D19	19" Seat Depth.....	No Charge
<input type="checkbox"/> D20	20" Seat Depth.....	No Charge
<input type="checkbox"/> D21	21" Seat Depth.....	No Charge
<input type="checkbox"/> D22	22" Seat Depth.....	No Charge

**SEAT OPTIONS**

<input type="checkbox"/> ADJASBA	Seat w/Adjustable Seat Tilt & Back Angle.....	\$774
5 <input type="checkbox"/> CTS	Contoura Seating.....	\$553

**ARMS**

**6 FIXED ARM TYPE (PAIR)**

<input type="checkbox"/> 18F	Full Length Fixed Height Conventional.....	\$244
<input type="checkbox"/> 28F	Desk Length Fixed Height Conventional.....	\$244

**6 ADJUSTABLE ARM TYPE (PAIR)**

<input type="checkbox"/> 19F	Full Length Adjustable Height Arms.....	\$348
<input type="checkbox"/> 29F	Desk Length Adjustable Height Arms.....	\$348

**6 MIXED ADJUSTABLE ARM TYPE (EACH)**

<input type="checkbox"/> 19FR	Full Length Adjustable Height Arm Right.....	\$174
<input type="checkbox"/> 29FR	Desk Length Adjustable Height Arm Right.....	\$174
<input type="checkbox"/> 19FL	Full Length Adjustable Height Arm Left.....	\$174
<input type="checkbox"/> 29FL	Desk Length Adjustable Height Arm Left.....	\$174

**BACKS**

**ASBA BACK TYPE**

<input type="checkbox"/> SECT	Sectional Rear Uprights.....	No Charge
<input type="checkbox"/> SRUP	Straight Rear Uprights.....	\$283
<input type="checkbox"/> 9720P	No Push Handles (Straight Canes).....	\$111
<input type="checkbox"/> L8661	Omit Rear Uprights.....	Less Net (\$25)

**ASBA BACK ANGLE DEGREE**

<input type="checkbox"/> 90BA	90°.....	<b>Standard</b>
<input type="checkbox"/> 95BA	95°.....	\$66

**BACK HEIGHT**

**ADJASBA**

<input type="checkbox"/> BH12	Back Height 12".....	No Charge
<input type="checkbox"/> BH13	Back Height 13".....	No Charge
<input type="checkbox"/> BH14	Back Height 14".....	No Charge
<input type="checkbox"/> BH15	Back Height 15".....	No Charge
<input type="checkbox"/> BH16	Back Height 16".....	No Charge
<input type="checkbox"/> BH17	Back Height 17".....	No Charge
<input type="checkbox"/> BH18	Back Height 18".....	No Charge
<input type="checkbox"/> BH19	Back Height 19".....	No Charge
<input type="checkbox"/> BH20	Back Height 20".....	No Charge
<input type="checkbox"/> BH21	Back Height 21".....	No Charge
<input type="checkbox"/> BH22	Back Height 22".....	No Charge
<input type="checkbox"/> BH23	Back Height 23".....	No Charge
<input type="checkbox"/> BH24	Back Height 24".....	No Charge

**CONTOURA**

<input type="checkbox"/> TSBH25	Back Height 25".....	No Charge
<input type="checkbox"/> TSBH295	Back Height 29.5".....	No Charge

**7 SEAT TILT**

<input type="checkbox"/> 0SA	0°.....	<b>Standard</b>
<input type="checkbox"/> 5SA	5°.....	No Charge

**FOOTNOTES**

- 5 Please refer to table on page 6 for seat style combinations.
- 6 Contoura available in even widths 16" - 22" and depths between 16" - 22" in 1" increments.
- 7 All arms are flip-up and detachable.

**WHEELS AND COLOURS**

**WHEEL OPTION**

**8 TDX3, TDX5**

- 1430-3T 14" X 3" Pneumatic Tires..... No Charge
- 1431-3T 14" X 3" Tire w/Foam Filled Inserts..... \$184

**FRAME FINISH**

**CONTEMPORARY**

- 24P Wet Black..... No Charge
- 60P Silver Metallic..... No Charge
- 71P Silver Vein..... No Charge
- 105P Midnight Blue..... No Charge
- 115P Black Prism..... No Charge

**TOO COOL**

- 30P Sunny Yellow..... No Charge
- 61P Electric Red..... No Charge
- 119P Bubblegum Pink..... No Charge
- 120P Grape Madness..... No Charge

**WAY COOL**

- 62P Electric Blue..... No Charge
- 104P Electric Teal..... No Charge
- 121P Lolly Pop Blue..... No Charge
- 122P Cosmic Blue..... No Charge
- 125P Emerald Green..... No Charge
- 127P Grasshopper Green..... No Charge

**MK6i SELECTION**

**MK6i CONTROL**

**TDX3**

- MK690ACC Controller w/Actuator Control ..... Standard
- MK6TT True Track Controller..... No Charge

**TDX5**

- MK6TT True Track Controller..... No Charge

**13 STANDARD DRIVER CONTROL**

- MPJM6 MPJ+ Multiple Joystick Drive Joystick..... **Standard**

**ALTERNATIVE PROPORTIONAL CONTROL OPTIONS**

- 9  PSFM6 Personal Joystick, Inductive on Front..... No Charge
- 9  PSRM6 Personal Joystick, Inductive on Rear..... No Charge
- 11  PBOD Push Button On/Off w/PB Drive Select..... No Charge
- 10,11  PBSS Push Button On/Off w/Speed Select..... No Charge
- 11  OR On/Off Right..... \$84
- 11  OL On/Off Left..... \$84
- 12  ASLRDYM6 MK6 Display Only (No Driver Control)..... \$226

**FOR ALL YOUR ASL DRIVE SYSTEM NEEDS PLEASE REFER TO THE FOLLOWING ASL PRICE LIST FORM 05-220C**

**FOOTNOTES**

- 8 The 12" solid drive wheel is standard on TDX3 SE and 14" drive wheels are optional. 14" flat free drive wheel standard on TDX3,TDX5 and anytime TrueTrack upgrade is selected.
- 9 TOD-SS: Toggle On/Off Drive Select, Speed, Speed Potentiometer Standard.
- 10 Drive Select through mode switch, or mode switch and Joystick Commands.
- 11 PSRM6 Only.
- 12 Display Only, No Driver Control. Allows choice from Order Form.

**ALTERNATIVE PROPORTIONAL CONTROL OPTIONS CONTINUED**

- 1500M6 Proportional RIM Head Control..... \$2,263
- SMHD Small Headrest..... \$225
- LGHD Large Headrest..... \$225
- 1558MM6 Compact Joystick w/ Display..... \$1,024

**ALTERNATIVE DIGITAL CONTROL SYSTEMS**

13	<input type="checkbox"/>	SNPM6	Sip-N-Puff.....	\$1,993
13	<input type="checkbox"/>	5018M6	Wafer Board w/SNP/DI.....	\$2,375
13	<input type="checkbox"/>	5020M6	Mini Tash Joystick & Digital Interface Box.....	\$2,506

**SIP N PUFF KIT**

	<input type="checkbox"/>	PKG32666	Therapin Sip-N-Puff Breath Tube Kit.....	\$411
--	--------------------------	----------	--	-------

**ADDITIONAL CONTROL CHOICES**

	<input type="checkbox"/>	1558M6	Compact Joystick - Less Std Mount.....	\$875
14	<input type="checkbox"/>	5018	Wafer Board.....	\$382
14	<input type="checkbox"/>	5020	Mini Joystick w/Push.....	\$513
	<input type="checkbox"/>	PACM6	Proportional Attendant Control.....	\$965
14	<input type="checkbox"/>	1552M	Attendant Control.....	\$529

**JOYSTICK MOUNTING HARDWARE**

15	<input type="checkbox"/>	QLAM6	Quad Link Retract Joystick Mount MPJ+.....	\$334
	<input type="checkbox"/>	LEFT	Left Hand Mounted Joystick.....	No Charge
	<input type="checkbox"/>	RIMHW	Rim Hardware For 1500.....	\$229
	<input type="checkbox"/>	PKG32669	Complete Bib Assembly (Tash/1558).....	\$379
	<input type="checkbox"/>	RCM	Rear Cane Mount for Compact Joystick.....	\$149

<input type="checkbox"/>	CGAT400-427I	Gatlin Midline Mount Horizontal - Invacare MK6i for MPJM6	\$914
<input type="checkbox"/>	CGAT400-426D	Gatlin Midline Mount Horizontal - Invacare for SPJPLUS	\$914
<input type="checkbox"/>	CGAT400V-427I	Gatlin Midline Mount Vertical (chin mount) - Invacare MMK6i for MPJM6	\$914
<input type="checkbox"/>	CGAT400V-426D	Gatlin Midline Mount Vertical (chin mount) - Invacare for SPJPLUS	\$914
<input type="checkbox"/>	CGAT400-480R	Gatlin Midline Mount Horizontal - Invacare 1812	\$914
<input type="checkbox"/>	CGAT400V-480R	Gatlin Midline Mount Vertical (chin mount) - Invacare 1812	\$914

**ELECTRONIC ACCESSORIES**

	<input type="checkbox"/>	1813M6	MK6i Programmer w/Pro Memory Card.....	\$575	<b>* UPDATED - 06/18/2007- Increase in Price</b>
	<input type="checkbox"/>	1813SD	Pro Memory Card.....	\$105	
	<input type="checkbox"/>	1560	"T" Handle Flexible Joystick Extension.....	\$86	
	<input type="checkbox"/>	1561	Straight Handle Flexible Joystick Extension.....	\$64	
	<input type="checkbox"/>	1826	Chin Cup.....	\$36	
	<input type="checkbox"/>	A24VPS	24V MK6i Auxiliary Power Source.....	\$235	
	<input type="checkbox"/>	AUX12M6	Auxiliary Module For 1 & 2.....	\$659	
	<input type="checkbox"/>	AUX34M6	Auxiliary Module For 3 & 4.....	\$659	
16	<input type="checkbox"/>	MPY	Mono Port "Y" Cable.....	\$20	

**POWER SEATING ACCESSORIES**

17	<input type="checkbox"/>	SANODE	Single Actuator Operation Through Driver Control.....	\$994
18	<input type="checkbox"/>	4WSB	4 Way Switch Box for Power Legs.....	No Charge
19	<input type="checkbox"/>	S4WSB	Multiple Actuator Interface Box.....	\$1,973
	<input type="checkbox"/>	FWT	4 Way Toggle.....	No Charge
	<input type="checkbox"/>	QPB	4 Push Buttons.....	No Charge

**POSITIONING STRAPS**

**SEAT POSITIONING STRAPS**

	<input type="checkbox"/>	1311BK	Airline Buckle Seat Positioning Strap.....	\$79
	<input type="checkbox"/>	1515	Auto Style Buckle Seat Positioning Strap.....	<b>Standard</b>

**CHEST POSITIONING STRAPS**

	<input type="checkbox"/>	7311BK	Airline Buckle Chest Positioning Strap.....	\$79
	<input type="checkbox"/>	9515	Auto Style Chest Positioning Strap.....	\$79
	<input type="checkbox"/>	7321BK	Hook & Loop Chest Pos. Strap.....	\$61

**BATTERY**

**BATTERY TRAY TYPE**

	<input type="checkbox"/>	22TRY	22NF Style Tray.....	<b>Standard</b>
--	--------------------------	-------	----------------------	-----------------

**BATTERIES**

	<input type="checkbox"/>	22NFBATTERY	22 NF Battery (On Chair).....	\$534
--	--------------------------	-------------	-------------------------------	-------

**FOOTNOTES**

- 13 Must also use either a Multiple Drive Joystick or MK6 Display.
- 14 Requires SNPM6 Sip-N-Puff Interface Box.
- 15 QLA includes ARM250 standard.
- 16 For use with MPJ+ or MK6i Display to provide additional switch site.
- 17 Required for Elevate Only Systems Operation through Driver Control - Multiple Drive Systems.
- 18 Required with all multiple actuator systems. Provides a D-9 port for the 4-Way Toggle or 4 Push Buttons.
- 19 Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control.

**ACCESSORIES**

**BACK ACCESSORIES**

- L8601 Omit Back Upholstery..... Less Net (\$20)
- RBP Reinforced Back Upholstery..... \$84

**MISCELLANEOUS ACCESSORIES**

- 1303 Wheel Locks..... \$81
- OMCU Omit Cushion..... Less Net (\$25)

**FRONT RIGGING OPTIONS**

<sup>20</sup> **FRONT RIGGING**

<sup>21</sup> **CENTER MOUNT**

- CTMT-F Centre Mount Adjustable Knee Angle 12 1/4"-15"..... \$262
- CTMT-FJ Junior Centre Mount Adj. Knee Angle 9 1/4"-12"..... \$262
- CTMT-FL Center Mount Rigging Extension 16"-18 3/4"..... \$416

**CENTRE PIVOT STYLE**

- 70NHD 70° Swingaway Footrests (Steel)..... \$225
- <sup>22</sup>  70STAPER 70° Tapered Footrests (Aluminum)..... \$262
- AT5544 Manual Elevating Legrests (Aluminum)..... \$529

**PIN STYLE-HEMI (PAIR)**

- PHW93 70° Swingaway Footplates (Steel)..... \$163
- PHAL4A Manual Articulating/Elevating Legrests (Steel)..... \$366

**POWER ELEVATING LEGRESTS - ARTICULATING**

- ELRPW Power Elev. Legrests w/Articulation 20" Ext..... \$2,856
- ELRPW22 Power Elev. Legrests w/Articulation 22" Ext..... \$2,856

**FOOTPLATE OPTIONS (PAIR)**

- 1651 Flip-Up Composite Footplates..... No Charge
- 1350 Extra Large Aluminum Footplates..... \$99
- AT5543 Adjustable Angle Flip Up Footplates..... \$122
- AT55430 Adjustable Angle Flip Up Footplates Not Selected..... No Charge

**FOOTREST/LEGREST ACCESSORIES**

- 1337 Calf Strap..... \$59
- 1600BK Heel Loops w/Ankle Straps..... \$79
- LPT2 Longer Pivot Slide Tube For 1350..... \$123
- IMPACTG Impact Guards..... No Charge
- ALPT Longer Pivot And Slide Tube..... \$123
- ALPT4 Longer Pivot And Slide Tube 4"..... \$123

**FOOTNOTES**

<sup>20</sup> All riggings include a standard footplate with the exception of the 70NHD.  
<sup>21</sup> Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring.  
<sup>22</sup> Standard with impact guards.

Storm TDX for Formula Elevate Specifications			
	TDX5	TDX4	TDX3
<b>Maximum Speed MPH</b>			
TrueTrack	7.5 (std)	7.5	7.5
4-Pole	N/A	6.5	6.0
<b>Seat Widths</b>	16-22"	16-22"	16-22"
<b>Seat Depths</b>	16-22"	16-22"	16-22"
<b>Weight Limits</b>			
<300 lb.	N/A	X	X
<300 lb. TrueTrack	X	X	X
<b>Warranty</b>	Lifetime Limited on Base Frame Three (3) Year Limited on Seat Frame One (1) Year Limited on Electronics One (1) Year Limited on Actuators Five (5) Year Limited Warranty on GB™ Motors Eighteen (18) Months Limited Warranty on Gear Motor		

Formula Elevate Features and Specifications
<ul style="list-style-type: none"> <li>Standard with 300-lb. weight capacity.</li> <li>Traditional Rehab Seating (ASBA)</li> <li>Seat-to-floor height of 18" at 0° tilt.</li> <li>Speed Inhibit feature of 20% of maximum speed while elevated.</li> </ul>

**Extended Warranty Information:**  
**Electrical Components: 2 Years**  
**Motors (2 and 4-pole): 2 Years**  
**GB Motors: 5 Years**

\* Warranty valid to original purchaser who purchases this product within Canada  
 \* For more product warranty information, please refer to the product owner's manual

**\* UPDATED - 06/18/2007 - New Warranty Box**

Battery Boxes	
*Standard Battery Tray	
TDX5	Group 24
TDX4	Group 24 22 NF (optional)
TDX3	22 NF Group 24 (Optional)
*Sealed Gel Cell Batteries Recommended	

**Batteries (not included)**  
**22NF** batteries with terminal configuration (positive on the left and negative on the right) **MUST** be used. 22NF batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection.

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X
*Not available with Contoura Seat Cushion							
Back Heights	25" or 29.5"						

**Invacare Canada reserves the right to change prices and specifications without notice.**  
 HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4G 4Z4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)  
 Regional Centre - 16769 Boul. Hymus, Kirkland, P.Q. H9H 3L4 Phones: (514) 630-6080 (800) 363-0137(P.Q.) Fax: (514) 630-6120  
 Invacare®, and Yes, you can., are trademarks of Invacare Corporation.  
 WWW.INVACARE.CA