



NAME _____ PHONE NO _____

ACCOUNT# _____ P.O. _____ DATE _____

ADDRESS _____

FRAME

CHAIR TYPE/POWER SEATING BASE

- TDX5-PS TDX5 for Formula w/MK6i..... \$8,974
- TDX4-PS TDX4 for Formula w/MK6i..... \$7,791
- TDX3-PS TDX3 for Formula w/MK6i...(includes Std MPJM6)..... \$8,164

TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets..... \$325

SYSTEM TYPE

- T-TDX Formula Tilt for TDX..... \$6,711

USER WEIGHT LIMITS

TDX5

- U300 User Weight <300lbs. (<250 lbs. w/Vent)..... No Charge
- U400 User Weight <400lbs. (<350 lbs. w/Vent)..... No Charge

2 TDX4

- U400 User Weight <400lbs. (<350 lbs. w/Vent)..... No Charge

2 TDX3

- U300 User Weight <300lbs. (<250 lbs. w/Vent)..... No Charge

2 TRUETRACK PACKAGE FOR TDX3 AND TDX4

- TRUTK TrueTrack Upgrade..... \$1,465

5 PUSH HANDLES

- PCNB Push Handles for Conventional and BPO..... No Charge

WIDTH AND DEPTH

SEAT WIDTH RANGE

- 1620W 16" to 20" Width (Adjustable)..... No Charge
- 2024W 20" to 24" Width (Adjustable)..... \$429

6 CONVENTIONAL & CONTOURA

SEAT WIDTH SETTING

- W16 16" Seat Width..... No Charge
- W17 17" Seat Width (Not Available w/Contoura)..... No Charge
- W18 18" Seat Width..... No Charge
- W19 19" Seat Width (Not Available w/Contoura)..... No Charge
- W20 20" Seat Width..... No Charge
- W21 21" Seat Width ((Not Available w/Contoura)..... No Charge
- W22 22" Seat Width..... No Charge
- 7 W23 23" Seat Width (Not Available w/Contoura)..... No Charge
- 7 W24 24" Seat Width..... No Charge

ADJUSTABLE SEAT DEPTH RANGE

CONVENTIONAL

- 1619D 16" to 19" Deep (Adjustable)..... No Charge
- 2022D 20" to 22" Deep (Adjustable)..... \$355

***Transport Tie Down**

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair; in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

FOOTNOTES

- 1 Package includes GB motors and TrueTrack driving technology, which enables driver to track straighter regardless of the terrain. Top speed is 7.5 MPH. Must select MK6TT controller.
- 2 It is recommended to upgrade to Group 24 batteries and/or TrueTrack for improved battery charge life.
- 3 Package includes heavy-duty 4-pole motors and 5 MPH top speed unless TRUTK upgrade is selected.
- 4 Package includes 4-pole motors and 6.5 MPH top speed unless TRUTK upgrade is selected.
- 5 Push handles are standard on Contoura Seating and available as an option on the Conventional Seating and Back Posts Only. The push handle option adds approx. 3.5" to the height of the back cane on the Conventional Seating and Back Post Only.
- 6 Contoura only available in even widths - non adjustable in width. Please see chart on page 6.
- 7 Must pick 400 lb. user weight limit.

CONTOURA

- 1618D 16" to 18" Deep (Adjustable)..... No Charge
- 1921D 19" to 21" Deep (Adjustable)..... No Charge
- 22D 22" Deep..... No Charge

CONVENTIONAL & CONTOURA

SEAT DEPTH SETTING

- D16 16" Seat Depth..... No Charge
- D17 17" Seat Depth..... No Charge
- D18 18" Seat Depth..... No Charge
- D19 19" Seat Depth..... No Charge
- D20 20" Seat Depth..... No Charge
- D21 21" Seat Depth..... No Charge
- D22 22" Seat Depth..... No Charge

SEAT OPTIONS

- CNB Conventional Back..... Standard
- ⁸ BPO Back Posts Only..... No Charge
- ⁹ CTS Contoura Seating..... \$553
- OMCU Omit Contoura Seat Cushion..... Less Net (\$25)

ARMS

^{10,11} **ADJUSTABLE ARM TYPE (PAIR)**

- 19F Full Length Adjustable Height Arms..... \$348
- 29F Desk Length Adjustable Height Arms..... \$348

¹¹ **MIXED ADJUSTABLE ARM TYPE (EACH)**

- 19FR Full Length Adjustable Height Arm Right..... \$174
- 29FR Desk Length Adjustable Height Arm Right..... \$174
- 19FL Full Length Adjustable Height Arm Left..... \$174
- 29FL Desk Length Adjustable Height Arm Left..... \$174

BACK HEIGHT

CURVED ABS OR CONVENTIONAL BACK

- TSBH20 Back Height 20"..... No Charge
- TSBH21 Back Height 21"..... No Charge
- TSBH22 Back Height 22"..... No Charge
- TSBH23 Back Height 23"..... No Charge
- TSBH24 Back Height 24"..... No Charge
- TSBH25 Back Height 25"..... \$311
- TSBH26 Back Height 26"..... \$311

CONTOURA

- TSBH25 Back Height 25"..... No Charge
- TSBH295 Back Height 29.5"..... No Charge

HEADREST TYPE (OPTIONAL)

- TAR127700 Curved Headrest Standard
- TAR128000 Two-Step Headrest No Charge
- T0022400 Large Headrest \$97
- T0022600 Small Headrest \$97

**FOR ALL YOUR STEALTH NEEDS PLEASE REFER TO THE FOLLOWING
STEALTH PRICE LIST FORM 05-329C**

¹² **SEAT TILT**

- 0SA 0°..... Standard
- 5SA 5°..... No Charge

¹³ **OTTO BOCK ARM PADS (PAIR)**

- OB20400 One Piece Channel Armpads \$165
- OB20500 Large Forearm Pads..... \$156
- OB20600 Medium Forearm Pads..... \$156
- OB20700 Small Forearm Pads..... \$156

FOOTNOTES

- 8 Round canes are provided with BPO. BPO comes standard with a 20" back cane height. A 22", 23" or 24" back cane height is available through specials, quote # SM-050318052054 as a no charge.
- 9 Selected seat depth on Contoura assumes 1" of back cushion compression. For example: When choosing an 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick includes leatherette seat cushion.
- 10 Not available with Contoura seat cushion on 22" or 24" wide.
- 11 All arms are flip-up and detachable.
- 12 Indicates static tilt of seat. Impacts seat-to-floor height, every 5° raises the seat-to-floor height.
- 13 Not available with Contoura seat cushion on 16" or 18" wide.

OTTO BOCK HAND PADS (PAIR)

<input type="checkbox"/> OB23000	Large Flat Hand Pads.....	\$128
<input type="checkbox"/> OB23300	Medium Flat Hand Pads.....	\$128

LATERALS (PAIR)

CONVENTIONAL BACK

14 <input type="checkbox"/> TMSLM	Swingaway Laterals Medium	\$315
14 <input type="checkbox"/> TMSLL	Swingaway Laterals Large	\$315

CONTOURA SEATING

15 <input type="checkbox"/> CSLM	Swingaway Contoura Laterals Medium	\$315
15 <input type="checkbox"/> CSLL	Swingaway Contoura Laterals Large	\$315

WHEELS AND COLOURS

SEAT TO FLOOR HEIGHTS

<input type="checkbox"/> LOW	Seat To Floor 17 1/2".....	No Charge
<input type="checkbox"/> TALL	Seat To Floor 21 1/2".....	\$252

WHEEL OPTION

<input type="checkbox"/> 1430-3	14" X 3" Pneumatic Tires.....	Standard
<input type="checkbox"/> 1431-3	14" X 3" Tire w/Foam Filled Inserts.....	\$184

FRAME FINISH

CONTEMPORARY

<input type="checkbox"/> 24P	Wet Black.....	No Charge
<input type="checkbox"/> 60P	Silver Metallic.....	No Charge
<input type="checkbox"/> 71P	Silver Vein.....	No Charge
<input type="checkbox"/> 105P	Midnight Blue.....	No Charge
<input type="checkbox"/> 115P	Black Prism.....	No Charge

TOO COOL

<input type="checkbox"/> 30P	Sunny Yellow.....	No Charge
<input type="checkbox"/> 61P	Electric Red.....	No Charge
<input type="checkbox"/> 119P	Bubblegum PinkV.....	No Charge
<input type="checkbox"/> 120P	Grape Madness.....	No Charge

WAY COOL

<input type="checkbox"/> 62P	Electric Blue.....	No Charge
<input type="checkbox"/> 104P	Electric Teal.....	No Charge
<input type="checkbox"/> 121P	Lolly Pop Blue.....	No Charge
<input type="checkbox"/> 122P	Cosmic Blue.....	No Charge
<input type="checkbox"/> 125P	Emerald Green.....	No Charge
<input type="checkbox"/> 127P	Grasshopper Green.....	No Charge

MK6i SELECTION

MK6i CONTROL

TDX3, TDX4

<input type="checkbox"/> MK690	90 Amp Controller	Standard
<input type="checkbox"/> MK6TT	True Track Controller.....	No Charge

TDX5

<input type="checkbox"/> MK6TT	True Track Controller.....	No Charge
--------------------------------	----------------------------	-----------

EXPANDABLE CONTROLLER SYSTEMS

<input type="checkbox"/> MPJM6	MPJ+ Multiple Drive Joystick.....	Standard
15 <input type="checkbox"/> PSFM6	Personal Joystick, Inductive on Front	No Charge
15 <input type="checkbox"/> PSRM6	Personal Joystick, Inductive on Rear	No Charge
17 <input type="checkbox"/> PBOD	Push Button On/Off w/PB Drive Select.....	\$84
16,17 <input type="checkbox"/> PBSS	Push Button On/Off w/Speed Select.....	\$84
17 <input type="checkbox"/> OR	On/Off Right.....	No Charge
17 <input type="checkbox"/> OL	On/Off Left.....	No Charge
18 <input type="checkbox"/> ASLRDYM6	MK6 Display Only (No Driver Control).....	\$226

FOOTNOTES

- 14 For use with Contoura Seating only.
- 15 TOD-SS: Toggle On/Off Drive Select, Speed, Speed Potentiometer Standard.
- 16 Drive Select through mode switch, or mode switch and Joystick Commands.
- 17 PSRM6 Only.
- 18 Display Only, No Driver Control. Allows choice from Order Form.

FOR ALL YOUR ASL DRIVE SYSTEM NEEDS PLEASE REFER TO THE FOLLOWING ASL PRICE LIST FORM 05-220C

<input type="checkbox"/>	1500M6	Proportional RIM Head Control.....	\$2,263
<input type="checkbox"/>	SMHD	Small Headrest.....	\$225
<input type="checkbox"/>	LGHD	Large Headrest.....	\$225
<input type="checkbox"/>	1558MM6	Compact Joystick w/ Display.....	\$1,024

ALTERNATIVE DIGITAL CONTROL SYSTEMS

19	<input type="checkbox"/>	SNPM6	Sip-N-Puff.....	\$1,993
19	<input type="checkbox"/>	5018M6	Wafer Board w/SNP/DI.....	\$2,375
19	<input type="checkbox"/>	5020M6	Mini Tash Joystick & Digital Interface Box.....	\$2,506

SIP N PUFF KIT

<input type="checkbox"/>	PKG32666	Therapin Sip-N-Puff Breath Tube Kit.....	\$411
--------------------------	----------	--	-------

ADDITIONAL CONTROL CHOICES

	<input type="checkbox"/>	1558M6	Compact Joystick - Less Std Mount.....	\$875
20	<input type="checkbox"/>	5018	Wafer Board.....	\$382
20	<input type="checkbox"/>	5020	Mini Joystick w/Push.....	\$513
	<input type="checkbox"/>	PACM6	Proportional Attendant Control.....	\$965
20	<input type="checkbox"/>	1552M	Attendant Control.....	\$529

JOYSTICK MOUNTING HARDWARE

21	<input type="checkbox"/>	QLAM6	Quad Link Retract Joystick Mnt MPJ+.....	\$334
21	<input type="checkbox"/>	QLA	Quad Link Retract Joystick Mnt SPJ+.....	\$334
	<input type="checkbox"/>	LEFT	Left Hand Mounted Joystick.....	No Charge
	<input type="checkbox"/>	RIMHW	Rim Hardware For 1500.....	\$229
	<input type="checkbox"/>	PKG32669	Complete Bib Assembly (Tash/1558).....	\$379
	<input type="checkbox"/>	RCM	Rear Cane Mount for Compact Joystick.....	\$149

<input type="checkbox"/>	CGAT400-427I	Gatlin Midline Mount Horizontal - Invacare MK6i for MPJM6	\$914
<input type="checkbox"/>	CGAT400-426D	Gatlin Midline Mount Horizontal - Invacare for DPJ/SPJ & SPJPLUS	\$914
<input type="checkbox"/>	CGAT400V-427I	Gatlin Midline Mount Vertical (chin mount) - Invacare MMK6i for MPJM6	\$914
<input type="checkbox"/>	CGAT400V-426D	Gat.Mid. Mnt Vertical (chin mount) - Invacare for DPJ/SPJ & SPJPLUS	\$914
<input type="checkbox"/>	CGAT400-480R	Gatlin Midline Mount Horizontal - Invacare 1812	\$914
<input type="checkbox"/>	CGAT400V-480R	Gatlin Midline Mount Vertical (chin mount) - Invacare 1812	\$914

ELECTRONIC ACCESSORIES

<input type="checkbox"/>	1813M6	MK6i Programmer w/Pro Memory Card.....	\$575	* UPDATED - 06/18/2007- Increase in Price
<input type="checkbox"/>	1813SD	Pro Memory Card.....	\$105	
<input type="checkbox"/>	1560	"T" Handle Flexible Joystick Extension.....	\$86	
<input type="checkbox"/>	1561	Straight Handle Flexible Joystick Extension.....	\$64	
<input type="checkbox"/>	MPJLHS	MPJ+ Lateral Support Platform, Left	\$58	
<input type="checkbox"/>	MPJRHS	MPJ+ Lateral Support Platform, Right	\$58	
<input type="checkbox"/>	1826	Chin Cup.....	\$36	
<input type="checkbox"/>	A24VPS	24V MK6i Auxiliary Power Source.....	\$235	
<input type="checkbox"/>	AUX12M6	Auxillary Module For 1 & 2.....	\$659	
<input type="checkbox"/>	AUX34M6	Auxillary Module For 3 & 4.....	\$659	

POWER SEATING ACCESSORIES

22	<input type="checkbox"/>	SANODE	Single Actuator Operation Through Driver Control.....	\$994
23	<input type="checkbox"/>	4WSB	4 Way Switch Box for Power Legs.....	No Charge
24	<input type="checkbox"/>	S4WSB	Multiple Actuator Interface Box.....	\$1,973
	<input type="checkbox"/>	FWT	4 Way Toggle.....	No Charge
	<input type="checkbox"/>	QPB	4 Push Buttons.....	No Charge

POSITIONING STRAPS

SEAT POSITIONING STRAPS

<input type="checkbox"/>	1311BK	Airline Buckle Seat Positioning Strap.....	\$79
<input type="checkbox"/>	1515	Auto Style Buckle Seat Positioning Strap.....	Standard

CHEST POSITIONING STRAPS

<input type="checkbox"/>	7311BK	Airline Buckle Chest Positioning Strap.....	\$79
<input type="checkbox"/>	9515	Auto Style Chest Positioning Strap.....	\$79
<input type="checkbox"/>	7321BK	Hook & Loop Chest Pos. Strap.....	\$61

FOOTNOTES

- 19 Must also use either a Multiple Drive Joystick or MK6i Display.
- 20 Requires SNPM6 Sip-N-Puff Interface Box.
- 21 QLA includes ARM250 standard.
- 22 Allows for single actuator operation through the joystick.
- 23 Required with all multiple actuator systems. Provides a D-9 port for the 4-Way Toggle or 4 Push Buttons.
- 24 Replaces the Four Way Switch Box and allows multiple actuator operation

BATTERY

25,26 BATTERY TRAY TYPE

27	<input type="checkbox"/> 22TRY	22NF Style Tray.....	Standard
28,29	<input type="checkbox"/> 24TRY	24 Group Style Tray.....	No Charge
	<input type="checkbox"/> 22TRY3	22NF Style Tray for Vent	\$201

26,30 BATTERIES

	<input type="checkbox"/> 24BATTERY	24 Gel Battery (On Chair).....	\$654
29	<input type="checkbox"/> 22NFBATT3	22 NF Battery for Vent On Chair.....	\$800
31	<input type="checkbox"/> 22NFBATTERY	22 NF Battery (On Chair).....	\$534
32	<input type="checkbox"/> 22NFTRY3BAT2	22NF Battery, 2 Batteries (On Chair).....	\$534

33,34 VENT TRAYS

	<input type="checkbox"/> VENT	Articulating Vent Tray.....	\$1,141
--	-------------------------------	-----------------------------	----------------

ACCESSORIES

BACK ACCESSORIES

<input type="checkbox"/> L8601	Omit Back Upholstery.....	Less Net (\$20)
<input type="checkbox"/> RBP	Reinforced Back Upholstery.....	\$84

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/> 1303	Wheel Locks.....	\$81
<input type="checkbox"/> OMCU	Omit Cushion.....	Less Net (\$25)

FRONT RIGGING OPTIONS

35 FRONT RIGGING

36 CENTER MOUNT

<input type="checkbox"/> CTMT-F	Centre Mount Adj.Knee Angle 12 1/4"-15".....	\$262
<input type="checkbox"/> CTMT-FJ	Junior Centre Mount Adj. Knee Angle 9 1/4"-12".....	\$262
<input type="checkbox"/> CTMT-FL	Center Mount Rigging Extension 16"-18 3/4".....	\$416

37 CENTRE PIVOT STYLE

<input type="checkbox"/> 70NHD	70° Swingaway Footrests (Steel).....	\$225
38 <input type="checkbox"/> 70STAPER	70° Tapered Footrests (Aluminum).....	\$262
<input type="checkbox"/> AT5544	Manual Elevating Legrests (Aluminum).....	\$529

37 PIN STYLE-HEMI (PAIR)

<input type="checkbox"/> PHW93	70° Swingaway Footplates (Steel).....	\$163
<input type="checkbox"/> PHAL4A	Manual Articulating/Elevating Legrests (Steel).....	\$366

POWER ELEVATING LEGRESTS - ARTICULATING

39 <input type="checkbox"/> PCTMT	Power Elev. Center Mount.....	\$2,898
39 <input type="checkbox"/> PCTMT-L	Power Elev. Center Mount-Long.....	\$2,898
38,40 <input type="checkbox"/> ELRPW	Power Elev. Legrests w/Articulation 20" Ext.....	\$2,856
38,41 <input type="checkbox"/> ELRPW22	Power Elev. Legrests w/Articulation 22" Ext.....	\$2,856

FOOTPLATE OPTIONS (PAIR)

<input type="checkbox"/> 1651	Composite Flip-Up Footplates.....	No Charge
<input type="checkbox"/> 1350	Extra Large Aluminum Footplates.....	\$99
<input type="checkbox"/> AT5543	Adjustable Angle Flip Up Footplates.....	\$122
<input type="checkbox"/> AT5543O	Adjustable Angle Flip Up Footplates Not Selected.....	No Charge
42 <input type="checkbox"/> FBCTMT	Footboard for Centre Mount.....	\$262

FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/> 1337	Calf Strap.....	\$59
<input type="checkbox"/> 1600BK	Heel Loops w/Ankle Straps.....	\$79
<input type="checkbox"/> LPT2	Longer Pivot Slide Tube For 1350.....	\$123
<input type="checkbox"/> IMPACTG	Impact Guards.....	No Charge
<input type="checkbox"/> ALPT	Longer Pivot And Slide Tube.....	\$123
<input type="checkbox"/> ALPT4	Longer Pivot And Slide Tube 4".....	\$123

FOOTNOTES

- 25 Due to battery harness & tray, only batteries with terminal configuration & dimension per page 4 will work with the Invacare TDX chairs.
- 26 It is recommended to upgrade to Group 24 batteries and/or TrueTrack for improved battery charge life.
- 27 Standard with TDX3
- 28 Standard with TDX5 and TDX4, optional with TDX3
- 29 Not available with vent tray.
- 30 Sealed Gel Batteries highly recommended for use in Invacare power wheelchairs.
- 31 Available with Vent Option Only.
- 32 22NFTRY3BAT2 use when adding after Market Vent.
- 33 With vent use it is strongly recommended to upgrade to TrueTrack for improved battery charge life.
- 34 Not available with TRBKTS.
- 35 All riggings include a standard footplate with the exception of the 70NHD.
- 36 Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring.
- 37 Provides up to 1 1/2" of depth adjustment in 1/2" increments.
- 38 Standard with impact guards.
- 39 Provides 6.5" of articulation and 90°-25° of elevation Knee-to-heel is 10"-13" or 14"-17".
- 40 Independent powered elevating legrests include elevating and articulating legrests with angle adjustable footplates. Extension range 13" - 20". Sold as pair.
- 41 Independent powered elevating legrests include elevating and articulating legrests with angle adjustable footplates. Extension range 15" - 22". Sold as pair.
- 42 One piece footboard for Center Mount (CTMT-F/CTMT-FL and PCTMT/CTMT-L) Riggings.

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEAT AND BACK CUSHIONS
SEE PRICE LIST FORM # 05-338C.**

Conventional Seat Sizes	
Seat Widths	12" - 24"
Seat Depths	12" - 22"
Back Heights	16" - 26" (16"-20" only available with 12"-15" seat widths)

TDX for Formula Tilt Only Specifications				
	TDX5	TDX4	TDX3	TDX3SE
Maximum Speed MPH				
TrueTrack	7.5	7.5 (option)	7.5 (option)	7.5 (option)
4-pole	N/A	N/A	6.5	5
4-pole HD	N/A	5	N/A	N/A
Seat-To-Floor at 0°				
Low	18"	18"	18"	18"
High	19.25"	19.25"	19.25"	19.25"
at 5°				
Low	19.75"	19.75"	19.75"	19.75"
High	21"	21"	21"	21"
Weight Limits				
<300 lb. (250 lb. w/vent)			X	X (No Vent)
<400 lb. (350 lb. w/vent)		X		
<300 lb. TrueTrack (250 lb. w/vent)	X			
<400 lb. TrueTrack (350 lb. w/vent)	X			
Warranty	Lifetime Limited on Base Frame Three (3) Year Limited on Seat Frame One (1) Year Limited on Electronics One (1) Year Limited on Actuators Five (5) Year Limited Warranty on GB™ Motors Eighteen (18) Months Limited Warranty on Gear Motor			

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X
*Not available with Contoura Seat Cushion							
Back Heights	25" or 29.5"						

Battery Boxes	
*Standard Battery Tray	
TDX5	Group 24
TDX4	Group 24 22 NF (optional)
TDX3	22 NF Group 24 (Optional)
*Sealed Gel Cell Batteries Recommended	

TDX for Formula Tilt Only Features and Specifications
<ul style="list-style-type: none"> • Weight capacity of 300 lb. or 400 lb. (with heavy-duty motors) • Tilt range from 0° to 55° or 5° to 60° • Conventional or Contoura style seating • Adjustable seat frame with Conventional style • Seat-to-floor height of 18" at 0° tilt

Batteries
22NF batteries with terminal configuration (positive on the left and negative on the right) MUST be used. 22NF batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection.
GP24 batteries with terminal configuration (negative on the left and positive on the right) MUST be used. GP24 batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection. See Owner's Manual, part number 1114809. These recommendations MUST be followed otherwise injury and damage may occur.

Extended Warranty Information:
 Electrical Components: **2 Years**
 Motors (2 and 4-pole): **2 Years**
 GB Motors: **5 Years**
 * Warranty valid to original purchaser who purchases this product within Canada
 * For more product warranty information, please refer to the product owner's manual

*** UPDATED - 06/18/2007 - New Warranty Box**

Invacare Canada reserves the right to change prices and specifications without notice.
 HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4G 4Z4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)
 Regional Centre - 16769 Boul. Hymus, Kirkland, P.Q. H9H 3L4 Phones: (514) 630-6080 (800) 363-0137(P.Q.) Fax: (514) 630-6120
 Invacare®, and Yes, you can., are trademarks of Invacare Corporation.
 WWW.INVACARE.CA