



NAME \_\_\_\_\_ PHONE NO \_\_\_\_\_

ACCOUNT# \_\_\_\_\_ P.O. \_\_\_\_\_ DATE \_\_\_\_\_

ADDRESS \_\_\_\_\_

**FRAME**

**CHAIR TYPE/POWER SEATING BASE**

- TDX5-PS TDX5 for Formula w/MK6i..... \$8,974
- 1  TDX3-PS TDX3 for Formula w/MK6i...(includes Std MPJM6)..... \$8,164

**TRANSPORT TIE DOWN**

- TRBKTS Wheelchair Transport Brackets..... \$325

**SYSTEM TYPE**

- 2  TRE-TDX Formula Tilt/Recline for TDX..... \$9,810
- 3  E2-TDX Elevate to Assist for Transfers..... \$2,716

**USER WEIGHT LIMITS**

**TDX3, TDX5**

- U300 User Weight <300lbs..... No Charge

**TRUETRACK PACKAGE FOR TDX3**

- 4  TRUTK True Track Upgrade..... \$1,465

**WIDTH AND DEPTH**

**CONVENTIONAL SEAT WIDTH RANGE**

- 1620W 16" to 20" Width (Adjustable)..... No Charge
- 2024W 20" to 24" Width (Adjustable)..... \$429

**5 CONVENTIONAL & CONTOURA**

**SEAT WIDTH SETTING**

- W16 16" Seat Width..... No Charge
- W17 17" Seat Width (Conventional Only)..... No Charge
- W18 18" Seat Width..... No Charge
- W19 19" Seat Width (Conventional Only)..... No Charge
- W20 20" Seat Width..... No Charge
- W21 21" Seat Width (Conventional Only)..... No Charge
- W22 22" Seat Width..... No Charge
- W23 23" Seat Width (Conventional Only)..... No Charge
- W24 24" Seat Width..... No Charge

**ADJUSTABLE SEAT DEPTH RANGE**

**CONVENTIONAL**

- 1619D 16" to 19" Deep (Adjustable)..... No Charge
- 1922D 19" to 22" Deep (Adjustable)..... \$355

**6 CONTOURA**

- 1618D 16" to 18" Deep (Adjustable)..... No Charge
- 1921D 19" to 21" Deep (Adjustable)..... No Charge
- 22D 22" Deep..... No Charge

**\*Transport Tie Down**

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

\*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

**FOOTNOTES**

- 1 Package includes 4-pole motors and 6.5 MPH.
- 2 It is strongly recommended to upgrade to TrueTrack for improved battery charge life.
- 3 Must choose the E2-TDX option with the TE-TDX. If you want the Formula Tilt only (T-TDX), see form number 05-406C - Formula Tilt Only pricelist.
- 4 Package includes the TrueTrack driving technology, which enables driver to track straighter regardless of the terrain (7.5 MPH). Weight limit is 300 lbs. Must have upgraded MK6TT controller. When selecting this option, heavy-duty 4-pole motors picked in the user weight limit category will be upgraded to Gearless Brushless GB™ motors. TDX5 comes standard with TrueTrack Package.
- 5 Contoura only available in even widths - non adjustable in width. Please see chart on page 6.
- 6 Selected seat depth on Contoura assumes 1" of back cushion compression. For example: When choosing an 18" depth, the measurement from the front of the back cushion to the front of the seat equals 17". The Contoura back cushion is 2" thick.

**6 CONVENTIONAL & CONTOURA**

**SEAT DEPTH SETTING**

- 7  D16 16" Seat Depth..... No Charge
- 7  D17 17" Seat Depth..... No Charge
- D18 18" Seat Depth..... No Charge
- D19 19" Seat Depth..... No Charge
- D20 20" Seat Depth..... No Charge
- 8  D21 21" Seat Depth..... No Charge
- 8  D22 22" Seat Depth..... No Charge

**SEAT OPTIONS**

- CNB Conventional Back..... Standard
- CTS Contoura Seating..... \$553
- OMCU Omit Cushion..... Less Net (\$25)

**ARMS**

**9 ARM TYPE**

- RA19FR Reclining Full Length Adjustable Height Arm Right..... No Charge
- RA29FR Reclining Desk Length Adjustable Height Arm Right..... No Charge
- RA19FL Reclining Full Length Adjustable Height Arm Left..... No Charge
- RA29FL Reclining Desk Length Adjustable Height Arm Left..... No Charge

**BACKS**

**BACK HEIGHT**

**CONVENTIONAL**

- TSBH20 Back Height 20"..... No Charge
- TSBH21 Back Height 21"..... No Charge
- TSBH22 Back Height 22"..... No Charge
- TSBH23 Back Height 23"..... No Charge
- TSBH24 Back Height 24"..... No Charge
- TSBH25 Back Height 25"..... \$311
- TSBH26 Back Height 26"..... \$311

**CONTOURA**

- TSBH25 Back Height 25"..... No Charge
- TSBH295 Back Height 29.5"..... No Charge

**HEADREST TYPE (OPTIONAL)**

- TAR127700 Curved Headrest ..... Standard
- TAR128000 Two-Step Headrest ..... No Charge
- T0022400 Large Headrest ..... \$97
- T0022600 Small Headrest ..... \$97

**FOR ALL YOUR STEALTH NEEDS PLEASE REFER TO THE FOLLOWING  
STEALTH PRICE LIST FORM 05-329C**

**10 SEAT TILT**

- 0SA 0° (Tilt Range 0° - 55°)..... Standard
- 5SA 5° (Tilt Range 5° - 60°)..... No Charge

**HEADREST TYPE (OPTIONAL)**

- TAR127700 Curved Headrest ..... Standard
- TAR128000 Two-Step Headrest ..... No Charge
- T0022400 Large Headrest ..... \$97
- T0022600 Small Headrest ..... \$97

**11 OTTO BOCK ARM PADS (PAIR)**

- OB20400 One Piece Channel Armpads ..... \$165
- OB20500 Large Forearm Pads..... \$156
- OB20600 Medium Forearm Pads..... \$156
- OB20700 Small Forearm Pads..... \$156

**FOOTNOTES**

- 7 Not available with Contoura seat cushion on 22" wide.
- 8 Not available with Contoura seat cushion on 16" or 18" wide.
- 9 All arms are flip-up and detachable.
- 10 Indicates static tilt of seat. Impacts seat-to-floor height, every 3° raises the seat-to-floor height.
- 11 Not compatible with quad link retractable joystick mount.

**11 OTTO BOCK HAND PADS**

- OB23000 Large Flat Hand Pads - Pair..... \$128
- OB23300 Medium Flat Hand Pads - Pair..... \$128

**LATERALS (PAIR)**

**12 CONVENTIONAL BACK**

- TMSLM Swingaway Laterals Medium ..... \$315
- TMSLL Swingaway Laterals Large ..... \$315

**13 CONTOURA SEATING**

- CSLM Swingaway Contoura Laterals Medium ..... \$315
- CSLL Swingaway Contoura Laterals Large ..... \$315

**WHEELS AND COLOURS**

**WHEEL OPTION**

**TDX3 & TDX5**

- 1430-3 14" X 3" Pneumatic Tires..... **Standard**
- 1431-3 14" X 3" Tire w/Foam Filled Inserts..... \$184

**FRAME FINISH**

**CONTEMPORARY**

- 24P Wet Black..... No Charge
- 60P Silver Metallic..... No Charge
- 71P Silver Vein..... No Charge
- 105P Midnight Blue..... No Charge
- 115P Black Prism..... No Charge

**TOO COOL**

- 30P Sunny Yellow..... No Charge
- 61P Electric Red..... No Charge
- 119P Bubblegum Pink..... No Charge
- 120P Grape Madness..... No Charge

**WAY COOL**

- 62P Electric Blue..... No Charge
- 104P Electric Teal..... No Charge
- 121P Lolly Pop Blue..... No Charge
- 122P Cosmic Blue..... No Charge
- 125P Emerald Green..... No Charge
- 127P Grasshopper Green..... No Charge

**MK6i SELECTION**

**MK6i CONTROL**

**TDX3**

- MK690 90AMP Controller ..... Standard
- MK6TT True Track Controller..... No Charge

**TDX5**

- MK6TT True Track Controller..... No Charge

**STANDARD DRIVER CONTROL**

- MPJM6 MPJ+ Multiple Joystick Drive Joystick..... **Standard**

**ALTERNATIVE PROPORTIONAL CONTROL OPTIONS**

- 14  PSFM6 Personal Joystick, Inductive on Front..... No Charge
- 14  PSRM6 Personal Joystick, Inductive on Rear..... No Charge
- 16  PBOD Push Button On/Off w/PB Drive Select..... \$84
- 15,16  PBSS Push Button On/Off w/Speed Select..... \$84
- 16  OR On/Off Right..... No Charge
- 16  OL On/Off Left..... No Charge
- 17  ASLRDYM6 MK6 Display Only (No Driver Control)..... \$226

**FOR ALL YOUR ASL DRIVE SYSTEM NEEDS PLEASE REFER TO THE FOLLOWING ASL PRICE LIST FORM 05-220C**

- 1500M6 Proportional RIM Head Control..... \$2,263
- SMHD Small Headrest..... \$225
- LGHD Large Headrest..... \$225
- 1558MM6 Compact Joystick w/ Display..... \$1,024

**FOOTNOTES**

- 12 For use with channel canes only.
- 13 For use with Contoura Seating only.
- 14 TOD-SS: Toggle On/Off-Drive Select, Speed Potentiometer Standard.
- 15 Drive Select through mode switch, or mode switch and Joystick Commands.
- 16 PSRM6 Only.
- 17 Display only, no driver control. Allows choice from Price Sheet.

**ALTERNATIVE DIGITAL CONTROL SYSTEMS**

18	<input type="checkbox"/>	SNPM6	Sip-N-Puff.....	\$1,993
18	<input type="checkbox"/>	5018M6	Wafer Board w/SNP/DI.....	\$2,375
18	<input type="checkbox"/>	5020M6	Mini Tash Joystick & Digital Interface Box.....	\$2,506

**SIP N PUFF KIT**

	<input type="checkbox"/>	PKG32666	Therapin Sip-N-Puff Breath Tube Kit.....	\$411
--	--------------------------	----------	--	-------

**ADDITIONAL CONTROL CHOICES**

	<input type="checkbox"/>	1558M6	Compact Joystick - Less Std Mount.....	\$875
19	<input type="checkbox"/>	5018	Wafer Board.....	\$382
19	<input type="checkbox"/>	5020	Mini Joystick w/Push.....	\$513
	<input type="checkbox"/>	PACM6	Proportional Attendant Control.....	\$965
19	<input type="checkbox"/>	1552M	Attendant Control.....	\$529

**JOYSTICK MOUNTING HARDWARE**

20	<input type="checkbox"/>	QLAM6	Quad Link Retract Joystick Mount MPJ+.....	\$334
	<input type="checkbox"/>	LEFT	Left Hand Mounted Joystick.....	No Charge
	<input type="checkbox"/>	RIMHW	Rim Hardware For 1500.....	\$229
	<input type="checkbox"/>	PKG32669	Complete Bib Assembly (Tash/1558).....	\$379
	<input type="checkbox"/>	RCM	Rear Cane Mount for Compact Joystick.....	\$149

<input type="checkbox"/>	CGAT400-427I	Gatlin Midline Mount Horizontal - Invacare MK6i for MPJM6	\$914
<input type="checkbox"/>	CGAT400-426D	Gatlin Midline Mount Horizontal - Invacare for SPJPLUS	\$914
<input type="checkbox"/>	CGAT400V-427I	Gatlin Midline Mount Vertical (chin mount) - Invacare MMK6i for MPJM6	\$914
<input type="checkbox"/>	CGAT400V-426D	Gatlin Midline Mount Vertical (chin mount) - Invacare for SPJPLUS	\$914
<input type="checkbox"/>	CGAT400-480R	Gatlin Midline Mount Horizontal - Invacare 1812	\$914
<input type="checkbox"/>	CGAT400V-480R	Gatlin Midline Mount Vertical (chin mount) - Invacare 1812	\$914

**ELECTRONIC ACCESSORIES**

	<input type="checkbox"/>	1813M6	MK6i Programmer w/Pro Memory Card.....	\$575	<b>* UPDATED - 06/18/2007- Increase in Price</b>
	<input type="checkbox"/>	1813SD	Pro Memory Card.....	\$105	
	<input type="checkbox"/>	1560	"T" Handle Flexible Joystick Extension.....	\$86	
	<input type="checkbox"/>	1561	Straight Handle Flexible Joystick Extension.....	\$64	
	<input type="checkbox"/>	1826	Chin Cup.....	\$36	
	<input type="checkbox"/>	A24VPS	24V MK6i Auxiliary Power Source.....	\$235	
	<input type="checkbox"/>	AUX12M6	Auxiliary Module For 1 & 2.....	\$659	
	<input type="checkbox"/>	AUX34M6	Auxiliary Module For 3 & 4.....	\$659	
21	<input type="checkbox"/>	MPY	Mono Port "Y" Cable.....	\$20	

**POWER SEATING ACCESSORIES**

22	<input type="checkbox"/>	4WSB	4 Way Switch Box for Power Legs.....	No Charge
23	<input type="checkbox"/>	S4WSB	Multiple Actuator Interface Box.....	\$1,973
	<input type="checkbox"/>	FWT	4 Way Toggle.....	No Charge
	<input type="checkbox"/>	QPB	4 Push Buttons.....	No Charge

**POSITIONING STRAPS**

**SEAT POSITIONING STRAPS**

	<input type="checkbox"/>	1311BK	Airline Buckle Seat Positioning Strap.....	\$79
	<input type="checkbox"/>	1515	Auto Style Buckle Seat Positioning Strap.....	<b>Standard</b>

**CHEST POSITIONING STRAPS**

	<input type="checkbox"/>	7311BK	Airline Buckle Chest Positioning Strap.....	\$79
	<input type="checkbox"/>	9515	Auto Style Chest Positioning Strap.....	\$79
	<input type="checkbox"/>	7321BK	Hook & Loop Chest Pos. Strap.....	\$61

**BATTERY**

**24 BATTERY TRAY TYPE**

	<input type="checkbox"/>	22TRY	22NF Style Tray.....	<b>Standard</b>
--	--------------------------	-------	----------------------	-----------------

**25 BATTERIES**

	<input type="checkbox"/>	22NFBATTERY	22 NF Battery (On Chair).....	\$534
--	--------------------------	-------------	-------------------------------	-------

**FOOTNOTES**

- 18 Must also use either a Multiple Drive Joystick or MK6 Display.
- 19 Requires SMPM6 Sip-N-Puff interface Box.
- 20 QLA includes ARM250 standard.
- 21 For use with MPJ+ or MK6i. Display to provide additional switch site.
- 22 Required with all multiple actuator systems. Provides a D-9 port for the 4-Way Toggle or 4 push buttons.
- 23 Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control.
- 24 Due to the battery harness & tray, only batteries with terminal configuration & dimension per page 6 will work with the Invacare TDX chairs.
- 25 Sealed Gel Batteries highly recommended for use in Invacare power wheelchairs.

**ACCESSORIES**

**BACK ACCESSORIES**

- L8601 Omit Back Upholstery.....Less Net (\$20)
- RBP Reinforced Back Upholstery..... \$84

**MISCELLANEOUS ACCESSORIES**

- 1303 Wheel Locks..... \$81
- OMCU Omit Cushion.....Less Net (\$25)

**FRONT RIGGING OPTIONS**

**26 FRONT RIGGING**

**27 CENTER MOUNT**

- CTMT-F Centre Mount Adjustable Knee Angle 12 1/4"-15"..... \$262
- CTMT-FJ Junior Centre Mount Adj. Knee Angle 9 1/4"-12"..... \$262
- CTMT-FL Center Mount Rigging Extension 16"-18 3/4"..... \$416

**29 CENTRE PIVOT STYLE**

- 70NHD 70° Swingaway Footrests (Steel)..... \$225
- 70STAPER 70° Tapered Footrests (Aluminum)..... \$262
- AT5544 Manual Elevating Legrests (Aluminum)..... \$529

**29 PIN STYLE-HEMI (PAIR)**

- PHW93 70° Swingaway Footplates (Steel)..... \$163
- PHAL4A Manual Articulating/Elevating Legrests (Steel)..... \$366

**POWER ELEVATING LEGRESTS - ARTICULATING**

- 30 PCTMT Power Elev. Center Mount..... \$2,898
- 30 PCTMT-L Power Elev. Center Mount-Long..... \$2,898
- 28,31 ELRPW Power Elev. Legrests w/Articulation 20" Ext..... \$2,856
- 28,32 ELRPW22 Power Elev. Legrests w/Articulation 22" Ext..... \$2,856

**FOOTNOTES**

- 26 All riggings include a standard footplate with the exception of the 70NHD.
- 27 Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring.
- 28 Standard with impact guards.
- 29 Provides up to 1 1/2" of depth adjustment in 1/2" increments.
- 30 Provides 6.5" of articulation and 90° to 25° of elevation. Knee-to-heel is 10"-13" or 14"-17".
- 31 Independent powered elevating legrests include elevating and articulating legrests with angle adjustable footplates. **Extension range 13"-20"**. Sold as pair.
- 32 Independent powered elevating legrests include elevating and articulating legrests with angle adjustable footplates. **Extension range 15" - 22"**. Sold as pair.
- 33 One piece footboard for Center Mount (CTMT-F/CTMT-FL and PCTMT/CTMT-L) Riggings.

**FOOTPLATE OPTIONS (PAIR)**

- 1651 Flip-Up Composite Footplates..... **No Charge**
- 1350 Extra Large Aluminum Footplates..... \$99
- AT5543 Adjustable Angle Flip Up Footplates..... \$122
- 33 FBCTMT Footboard for Centre Mount..... \$262

**FOOTREST/LEGREST ACCESSORIES**

- 1337 Calf Strap..... \$59
- 1600BK Heel Loops w/Ankle Straps..... \$79
- LPT2 Longer Pivot Slide Tube For 1350..... \$123
- IMPACTG Impact Guards..... **No Charge**
- ALPT Longer Pivot And Slide Tube..... \$123
- ALPT4 Longer Pivot And Slide Tube 4"..... \$123

Battery Boxes	
*Standard Battery Tray	
TDX5	Group 24
TDX4	Group 24 22 NF (optional)
TDX3	22 NF Group 24 (Optional)
*Sealed Gel Cell Batteries Recommended	

**Batteries  
(not included)**

22NF batteries with terminal configuration (positive on the left and negative on the right) MUST be used. 22NF batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection.

Contoura Seat Sizes w/Elevate Combinations							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
*Not available with Contoura Seat Cushion.							
Back Heights	25" or 29.5"						

Conventional Seat Sizes	
Seat Widths	16" - 24"
Seat Depths	16" - 22"
Back Heights	20" - 26"

**FOR YOUR SEATING NEEDS  
PLEASE REFER TO THE FOLLOWING  
REHAB SEAT AND BACK CUSHIONS  
SEE PRICE LIST FORM # 05-338C.**

TDX for Formula Tilt/Recline/Elevate Specifications			
	TDX5	TDX4	TDX3
<b>Maximum Speed MPH</b>			
TrueTrack	7.5 (std)	7.5	7.5
4-pole	N/A	6.5	6.0
<b>Seat-To-Floor at 0°</b>			
Low	18"	18"	18"
High	19.25"	19.25"	19.25"
<b>at 5°</b>			
Low	19.75"	19.75"	19.75"
High	21"	21"	21"
<b>Weight Limits</b>			
<300 lb	N/A	X	X
<300 lb, True Track			
<b>Warranty</b>	Lifetime Limited on Base Frame Three (3) Year Limited on Seat Frame One (1) Year Limited on Electronics One (1) Year Limited on Actuators Five (5) Year Limited Warranty on GB™ Motors Eighteen (18) Months Limited Warranty on Gear Motor		

Formula Tilt/Recline/Elevate Features and Specifications	
•	Standard with 300 lb. weight capacity.
•	Tilt range from 0° to 55° or 5° to 60°.
•	Recline range of 170° Elevate range of 7".
•	Mechanical shear reduction of 3".
•	Conventional or Contoura style seating.
•	Adjustable seat frame with Conventional style.
•	Seat-to-floor height of 18" at 0° tilt.
•	Speed Inhibit feature of 20% of maximum speed while elevated.

**Extended Warranty Information:**

Electrical Components: **2 Years**  
Motors (2 and 4-pole): **2 Years**  
GB Motors: **5 Years**

\* Warranty valid to original purchaser who purchases this product within Canada  
\* For more product warranty information, please refer to the product owner's manual

**\* UPDATED - 06/18/2007 - New Warranty Box**

Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4G 4Z4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)  
Regional Centre - 16769 Boul. Hymus, Kirkland, P.Q. H9H 3L4 Phones: (514) 630-6080 (800) 363-0137(P.Q.) Fax: (514) 630-6120

Invacare®, and Yes, you can., are trademarks of Invacare Corporation.  
WWW.INVACARE.CA

NOTE: All specifications and dimensions are approximate.  
For the most current pricing information visit [www.invacare.ca](http://www.invacare.ca)