



NAME _____ PHONE NO _____

ACCOUNT# _____ P.O. _____ DATE _____

ADDRESS _____

FRAME

1 CHAIR TYPE/POWER SEATING BASE

- TDX5-PS TDX5 for Formula w/MK6i..... \$8,974
- TDX4-PS TDX4 for Formula w/MK6i..... \$7,791
- TDX3-PS TDX3 for Formula w/MK6i...(includes Std MPJM6)..... \$8,164

TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets..... \$325

SYSTEM TYPE

- TR-TDX Formula Tilt/Recline for TDX..... \$9,810

USER WEIGHT LIMITS

2 TDX5

- U400T5 True Track & User Weight <400lbs. No Charge

3 TDX4

- U300T4 User Weight <300lbs. No Charge
- U400T4 User Weight <400lbs. No Charge

3 TDX3

- U300 User Weight <300lbs. No Charge

2 TRUETRACK PACKAGE FOR TDX3 & TDX4

- TRUTK TrueTrack upgrade..... \$1,465

WIDTH AND DEPTH

CONVENTIONAL SEAT WIDTH RANGE

- 1620W 16" to 20" Width (Adjustable)..... No Charge
- 2024W 20" to 24" Width (Adjustable)..... \$429

6 CONVENTIONAL & CONTOURA

SEAT WIDTH SETTING

- W16 16" Seat Width..... No Charge
- W17 17" Seat Width (Conventional Only)..... No Charge
- W18 18" Seat Width..... No Charge
- W19 19" Seat Width (Conventional Only)..... No Charge
- W20 20" Seat Width..... No Charge
- W21 21" Seat Width (Conventional Only)..... No Charge
- W22 22" Seat Width..... No Charge
- W23 23" Seat Width (Conventional Only)..... No Charge
- W24 24" Seat Width..... No Charge

ADJUSTABLE SEAT DEPTH RANGE

CONVENTIONAL

- 1619D 16" to 19" Deep (Adjustable)..... No Charge
- 1922D 19" to 22" Deep (Adjustable)..... \$355

CONTOURA

- 1618D 16" to 18" Deep (Adjustable)..... No Charge
- 1921D 19" to 21" Deep (Adjustable)..... No Charge
- 22D 22" Deep..... No Charge

*Transport Tie Down
<p>TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.</p> <p>As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.</p> <p>*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.</p>

FOOTNOTES

- 1 With the vent tray option the maximum weight limit is less 50 lbs.
- 2 Package includes GB™ Motors and TrueTrack driving technology, which enables driver to track straighter regardless of the terrain. 7.5 MPH top speed. Must select MK6TT controller. TDX 5 comes standard with TRUTK package.
- 3 It is recommended to upgrade to Group 24 batteries and/or TrueTrack for improved battery charge life.
- 4 Package includes 4 pole motors and 6.5 MPH top speed unless TRUTK upgrade is selected.
- 5 Package includes heavy-duty 4 pole motors and 5 MPH top speed unless TRUTK upgrade is selected.
- 6 Contoura is only available in even widths - non adjustable in width. Please see chart on page 5.
- 7 Must pick heavy duty motor package for this option. U400T5, U400T4, HD400T4, U400T3, UTT4.

CONVENTIONAL & CONTOURA

SEAT DEPTH SETTING

8	<input type="checkbox"/> D16	16" Seat Depth.....	No Charge
8	<input type="checkbox"/> D17	17" Seat Depth.....	No Charge
	<input type="checkbox"/> D18	18" Seat Depth.....	No Charge
	<input type="checkbox"/> D19	19" Seat Depth.....	No Charge
	<input type="checkbox"/> D20	20" Seat Depth.....	No Charge
9	<input type="checkbox"/> D21	21" Seat Depth.....	No Charge
9	<input type="checkbox"/> D22	22" Seat Depth.....	No Charge

SEAT OPTIONS

	<input type="checkbox"/> CNB	Conventional Back.....	Standard
10	<input type="checkbox"/> CTS	Contoura Seating.....	\$553

ARMS

11 ARM TYPE

	<input type="checkbox"/> RA19FR	Reclining Full Length Adjustable Height Arm Right.....	No Charge
	<input type="checkbox"/> RA29FR	Reclining Desk Length Adjustable Height Arm Right.....	No Charge
	<input type="checkbox"/> RA19FL	Reclining Full Length Adjustable Height Arm Left.....	No Charge
	<input type="checkbox"/> RA29FL	Reclining Desk Length Adjustable Height Arm Left.....	No Charge

BACKS

BACK HEIGHT

CONVENTIONAL

	<input type="checkbox"/> TSBH20	Back Height 20".....	No Charge
	<input type="checkbox"/> TSBH21	Back Height 21".....	No Charge
	<input type="checkbox"/> TSBH22	Back Height 22".....	No Charge
	<input type="checkbox"/> TSBH23	Back Height 23".....	No Charge
	<input type="checkbox"/> TSBH24	Back Height 24".....	No Charge
	<input type="checkbox"/> TSBH25	Back Height 25".....	\$311
	<input type="checkbox"/> TSBH26	Back Height 26".....	\$311

CONTOURA

	<input type="checkbox"/> TSBH25	Back Height 25".....	No Charge
	<input type="checkbox"/> TSBH295	Back Height 29.5".....	No Charge

HEADREST TYPE (OPTIONAL)

	<input type="checkbox"/> TAR127700	Curved Headrest	Standard
	<input type="checkbox"/> TAR128000	Two-Step Headrest	No Charge
	<input type="checkbox"/> T0022400	Large Headrest	\$97
	<input type="checkbox"/> T0022600	Small Headrest	\$97

FOR ALL YOUR STEALTH NEEDS PLEASE REFER TO THE FOLLOWING STEALTH PRICE LIST FORM 05-329C

SEAT TILT

	<input type="checkbox"/> 0SA	0° (Tilt Range 0° - 55°).....	Standard
12	<input type="checkbox"/> 5SA	5° (Tilt Range 5° - 60°).....	No Charge

13 OTTO BOCK ARM PADS (PAIR)

	<input type="checkbox"/> OB20400	One Piece Channel Armpads	\$165
	<input type="checkbox"/> OB20500	Large Forearm Pads.....	\$156
	<input type="checkbox"/> OB20600	Medium Forearm Pads.....	\$156
	<input type="checkbox"/> OB20700	Small Forearm Pads.....	\$156

OTTO BOCK HAND PADS

	<input type="checkbox"/> OB23000	Large Flat Hand Pads - Pair.....	\$128
	<input type="checkbox"/> OB23300	Medium Flat Hand Pads - Pair.....	\$128

LATERALS (PAIR)

14 CONVENTIONAL BACK

	<input type="checkbox"/> TMSLM	Swingaway Laterals Medium	\$315
	<input type="checkbox"/> TMSLL	Swingaway Laterals Large	\$315

15 CONTOURA SEATING

	<input type="checkbox"/> CSLM	Swingaway Contoura Laterals Medium	\$315
	<input type="checkbox"/> CSLL	Swingaway Contoura Laterals Large	\$315

FOOTNOTES

- 8 Not available with Contoura seat cushion on 22" or 24" wide.
- 9 Not available with Contoura seat cushion on 16" or 18" wide.
- 10 Selected seat depth on Contoura assumes 1" of back cushion compression. For example: When choosing an 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick.
- 11 All arms are flip-up and detachable.
- 12 Indicates static tilt of seat. Impacts seat-to-floor height, every 5° raises the seat-to-floor height.
- 13 Not compatible with quad link retractable joystick mount.
- 14 For use with channel canes only.
- 15 For use with Contoura Seating only.

WHEELS AND COLOURS

SEAT TO FLOOR HEIGHTS

- LOW Seat-to-Floor Low..... No Charge
- TALL Seat-to-Floor Tall..... \$252

WHEEL OPTION

16 TDX3

- 1430-3T 14" X 3" Pneumatic Tires..... No Charge
- 1431-3T 14" X 3" Tire w/Foam Filled Inserts..... \$184

17 TDX4, TDX5

- 1430-3 14" X 3" Pneumatic Tires..... **Standard**
- 1431-3 14" X 3" Tire w/Foam Filled Inserts..... \$184

FRAME FINISH

CONTEMPORARY

- 24P Wet Black..... No Charge
- 60P Silver Metallic..... No Charge
- 71P Silver Vein..... No Charge
- 105P Midnight Blue..... No Charge
- 115P Black Prism..... No Charge

TOO COOL

- 30P Sunny Yellow..... No Charge
- 61P Electric Red..... No Charge
- 119P Bubblegum Pink..... No Charge
- 120P Grape Madness..... No Charge

WAY COOL

- 62P Electric Blue..... No Charge
- 104P Electric Teal..... No Charge
- 121P Lolly Pop Blue..... No Charge
- 122P Cosmic Blue..... No Charge
- 125P Emerald Green..... No Charge
- 127P Grasshopper Green..... No Charge

MK6i SELECTION

MK6i CONTROL

TDX3, TDX4

- MK690ACC Controller w/Actuator Control Standard
- MK6TT True Track Controller..... No Charge

TDX5

- MK6TT True Track Controller..... No Charge

STANDARD DRIVER CONTROL

- MPJM6 MPJ+ Multiple Joystick Drive Joystick..... **Standard**

ALTERNATIVE PROPORTIONAL CONTROL OPTIONS

- 18 PSFM6 Personal Joystick, Inductive on Front..... No Charge
- 18 PSRM6 Personal Joystick, Inductive on Rear..... No Charge
- 20 PBOD Push Button On/Off w/PB Drive Select..... \$84
- 19,20 PBSS Push Button On/Off w/Speed Select..... \$84
- 20 OR On/Off Right..... No Charge
- 20 OL On/Off Left..... No Charge
- 21 ASLRDYM6 MK6 Display Only (No Driver Control)..... \$226

FOR ALL YOUR ASL DRIVE SYSTEM NEEDS PLEASE REFER TO THE FOLLOWING ASL PRICE LIST FORM 05-220C

ALTERNATIVE PROPORTIONAL CONTROL OPTIONS CONTINUED

- 1500M6 Proportional RIM Head Control..... \$2,263
- SMHD Small Headrest..... \$225
- LGHD Large Headrest..... \$225
- 1558MM6 Compact Joystick w/ Display..... \$1,024

FOOTNOTES

- 16 The 12" pneumatic drive wheel is standard with the HD300T3 package. The 14" drive wheels are available as an option.
- 17 The 14" pneumatic drive wheel is standard with TDX3 motor packages; HD400T3, U400T3 and all TDX4/TDX5 motor options.
- 18 TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.
- 19 Drive Select through mode switch, or mode switch and Joystick Commands.
- 20 PSM6 Only.
- 21 Display Only, No Driver Control. Allows choice from order form.
- 22 Must also use either a Multiple Drive Joystick or MK6 Display.

ALTERNATIVE DIGITAL CONTROL SYSTEMS

- 22 SNPM6 Sip-N-Puff..... \$1,993
- 22 5018M6 Wafer Board w/SNP/DI..... \$2,375
- 22 5020M6 Mini Tash Joystick & Digital Interface Box..... \$2,506

SIP N PUFF KIT

23	<input type="checkbox"/> PKG32666	Therapin Sip-N-Puff Breath Tube Kit.....	\$411
----	-----------------------------------	--	-------

ADDITIONAL CONTROL CHOICES

	<input type="checkbox"/> 1558M6	Compact Joystick - Less Std Mount.....	\$875
23	<input type="checkbox"/> 5018	Wafer Board.....	\$382
23	<input type="checkbox"/> 5020	Mini Joystick w/Push.....	\$513
	<input type="checkbox"/> PACM6	Proportional Attendant Control.....	\$965
23	<input type="checkbox"/> 1552M	Attendant Control.....	\$529

JOYSTICK MOUNTING HARDWARE

24	<input type="checkbox"/> QLAM6	Quad Link Retract Joystick Mount MPJ+.....	\$334
	<input type="checkbox"/> LEFT	Left Hand Mounted Joystick.....	No Charge
	<input type="checkbox"/> RIMHW	Rim Hardware For 1500.....	\$229
	<input type="checkbox"/> PKG32669	Complete Bib Assembly (Tash/1558).....	\$379
	<input type="checkbox"/> RCM	Rear Cane Mount for Compact Joystick.....	\$149

<input type="checkbox"/> CGAT400-427I	Gatlin Midline Mount Horizontal - Invacare MK6i for MPJM6	\$914
<input type="checkbox"/> CGAT400-426D	Gatlin Midline Mount Horizontal - Invacare for SPJPLUS	\$914
<input type="checkbox"/> CGAT400V-427I	Gatlin Midline Mount Vertical (chin mount) - Invacare MMK6i for MPJM6	\$914
<input type="checkbox"/> CGAT400V-426D	Gatlin Midline Mount Vertical (chin mount) - Invacare for SPJPLUS	\$914
<input type="checkbox"/> CGAT400-480R	Gatlin Midline Mount Horizontal - Invacare 1812	\$914
<input type="checkbox"/> CGAT400V-480R	Gatlin Midline Mount Vertical (chin mount) - Invacare 1812	\$914

ELECTRONIC ACCESSORIES

	<input type="checkbox"/> 1813M6	MK6i Programmer w/Pro Memory Card.....	\$575
	<input type="checkbox"/> 1813SD	Pro Memory Card.....	\$105
	<input type="checkbox"/> 1560	"T" Handle Flexible Joystick Extension.....	\$86
	<input type="checkbox"/> 1561	Straight Handle Flexible Joystick Extension.....	\$64
	<input type="checkbox"/> 1826	Chin Cup.....	\$36
	<input type="checkbox"/> A24VPS	24V MK6i Auxiliary Power Source.....	\$235
	<input type="checkbox"/> AUX12M6	Auxiliary Module For 1 & 2.....	\$659
	<input type="checkbox"/> AUX34M6	Auxiliary Module For 3 & 4.....	\$659
25	<input type="checkbox"/> MPY	Mono Port "Y" Cable.....	\$20

*** UPDATED - 06/18/2007- Increase in Price**

POWER SEATING ACCESSORIES

26	<input type="checkbox"/> 4WSB	4 Way Switch Box for Power Legs.....	No Charge
27	<input type="checkbox"/> S4WSB	Multiple Actuator Interface Box.....	\$1,973
	<input type="checkbox"/> FWT	4 Way Toggle.....	No Charge
	<input type="checkbox"/> QPB	4 Push Buttons.....	No Charge

POSITIONING STRAPS

SEAT POSITIONING STRAPS

	<input type="checkbox"/> 1311BK	Airline Buckle Seat Positioning Strap.....	\$79
	<input type="checkbox"/> 1515	Auto Style Buckle Seat Positioning Strap.....	Standard

CHEST POSITIONING STRAPS

	<input type="checkbox"/> 7311BK	Airline Buckle Chest Positioning Strap.....	\$79
	<input type="checkbox"/> 9515	Auto Style Chest Positioning Strap.....	\$79
	<input type="checkbox"/> 7321BK	Hook & Loop Chest Pos. Strap.....	\$61

FOOTNOTES

- 23 Requires SNPM6 Sip-N-Puff Interface Box.
- 24 QLA includes ARM250.
- 25 For use with MPJ+ or MK6i Display to provide additional switch site.
- 26 Required with all multiple actuator systems. Provides a D-9 port for the 4-way toggle or 4 push buttons.
- 27 Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control.
- 28 Due to battery harness & tray, only batteries with terminal configuration & dimension per page 4 will work with the Invacare TDX chairs.
- 29 It is recommended to upgrade to Group 24 batteries and/or TrueTrack for improved battery charge life.
- 30 Standard with TDX3 and optional with TDX4.
- 31 Standard with TDX5 and TDX4, optional with TDX3.
- 32 Not available with vent tray.
- 33 Sealed Gel Batteries highly recommended for use in Invacare power wheelchairs.
- 34 With vent use it is strongly recommended to upgrade to TrueTrack for improved battery charge life.
- 35 Not available with TRBKTS.

BATTERY

28,29 BATTERY TRAY TYPE

30	<input type="checkbox"/> 22TRY	22NF Style Tray.....	Standard
31,32	<input type="checkbox"/> 24TRY	24 Group Style Tray.....	No Charge

29,33 BATTERIES

	<input type="checkbox"/> 24BATTERY	24 Gel Battery (On Chair).....	\$654
32	<input type="checkbox"/> 22NFBATTERY	22 NF Battery (On Chair).....	\$534

34,35 VENT TRAYS

	<input type="checkbox"/> VENT	Articulating Vent Tray.....	\$1,141
--	-------------------------------	-----------------------------	---------

INSTALLATION OF VENT BATTERIES

	<input type="checkbox"/> VENT22NF	Vent Back Up Battery - 22NF (On Chair).....	\$245
--	-----------------------------------	---	-------

ACCESSORIES

BACK ACCESSORIES

- L8601 Omit Back Upholstery..... Less Net (\$20)
- RBP Reinforced Back Upholstery..... \$84

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks..... \$81
- OMCU Omit Cushion.....Less Net (\$25)

FRONT RIGGING OPTIONS

36 FRONT RIGGING

37 CENTER MOUNT

- CTMT-F Centre Mount Adjustable Knee Angle 12 1/4"-15"..... \$262
- CTMT-FJ Junior Centre Mount Adj. Knee Angle 9 1/4"-12"..... \$262
- CTMT-FL Center Mount Rigging Extension 16"-18 3/4"..... \$416

CENTRE PIVOT STYLE

- 70NHD 70° Swingaway Footrests (Steel)..... \$225
- 38** 70STAPER 70° Tapered Footrests (Aluminum)..... \$262
- AT5544 Manual Elevating Legrests (Aluminum)..... \$529

PIN STYLE-HEMI (PAIR)

- PHW93 70° Swingaway Footplates (Steel)..... \$163
- PHAL4A Manual Articulating/Elevating Legrests (Steel)..... \$366

POWER ELEVATING LEGRESTS - ARTICULATING

- 39** PCTMT Power Elev. Center Mount..... \$2,898
- 39** PCTMT-L Power Elev. Center Mount-Long..... \$2,898
- 40** ELRPW Power Elev. Legrests w/Articulation 20" Ext..... \$2,856
- 40** ELRPW22 Power Elev. Legrests w/Articulation 22" Ext..... \$2,856

FOOTPLATE OPTIONS (PAIR)

- 1651 Composite Flip-Up Footplates..... No Charge
- 1350 Extra Large Aluminum Footplates..... \$99
- AT5543 Adjustable Angle Flip Up Footplates..... \$122

FOOTNOTES

- 36 All riggings include a standard footplate with the exception of the 70NHD.
- 37 Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring.
- 38 Standard with impact guards.
- 39 Provides 6.5" of articulation, 90° - 25° of elevation, 10"-13" or 14"-17" extension range.
- 40 Independent powered elevating legrests include elevating and articulating legrests with angle adjustable footplates. Extension range 13"-20" or 15"-22". Sold as pair.

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap..... \$59
- 1600BK Heel Loops w/Ankle Straps..... \$79
- LPT2 Longer Pivot Slide Tube For 1350..... \$123
- IMPACTG Impact Guards..... No Charge
- ALPT Longer Pivot And Slide Tube..... \$123
- ALPT4 Longer Pivot And Slide Tube 4"..... \$123

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING**

**REHAB SEAT AND BACK CUSHIONS
SEE PRICE LIST FORM # 05-338C.**

Batteries

22NF batteries with terminal configuration (positive on the left and negative on the right) MUST be used. 22NF batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection.

GP24 batteries with terminal configuration (negative on the left and positive on the right) MUST be used. GP24 batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection. See Owner's Manual, part number 1114809. These recommendations MUST be followed otherwise injury and damage may occur.

TDX for Formula Tilt/Recline Specifications

	TDX5	TDX4	TDX3
Maximum Speed MPH			
TrueTrack	7.5 (std)	7.5	7.5
4-pole	N/A	6.5	6.0
4-pole HD	N/A	5	5
Seat-To-Floor at 0°			
Low	18"	18"	18"
High	19.25"	19.25"	19.25"
at 5°			
Low	19.75"	19.75"	19.75"
High	21"	21"	21"
Weight Limits			
<300 lb. (250 lb. w/vent)	N/A	X	X
<400 lb. (350 lb. w/vent)	N/A	X	X
<400 lb. True Track (350 lb. w/vent)	X	X	X
Warranty	Lifetime Limited on Base Frame Three (3) Year Limited on Seat Frame One (1) Year Limited on Electronics One (1) Year Limited on Actuators Five (5) Year Limited Warranty on GB™ Motors Eighteen (18) Months Limited Warranty on Gear Motor		

Battery Boxes	
	*Standard Battery Tray
TDX5	Group 24
TDX4	Group 24 22 NF (optional)
TDX3	22 NF Group 24 (Optional)
*Sealed Gel Cell Batteries Recommended	

Contoura Seat Sizes

	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X

*Not available with Contoura Seat Cushion

Back Heights | 25" or 29.5"

Conventional Seat Sizes

Seat Widths	12" - 24"
Seat Depths	12" - 22"
Back Heights	16" - 26" (16"-20" only available with 12"-15" seat widths)

Formula Tilt/Recline Features and Specifications

- Weight capacity of 300 lb. or 400 lb. (with heavy-duty motors).
- Tilt range from 0° to 55° or 5° to 60°.
- Recline range of 170°.
- Mechanical shear reduction of 3".
- Conventional or Contoura style seating.
- Adjustable seat frame with Conventional style.
- Seat-to-floor height 18" at 0° tilt.

Extended Warranty Information:

Electrical Components: **2 Years**

Motors (2 and 4-pole): **2 Years**

GB Motors: **5 Years**

* Warranty valid to original purchaser who purchases this product within Canada

* For more product warranty information, please refer to the product owner's manual

*** UPDATED - 06/18/2007 - New Warranty Box**

Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4G 4Z4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)
Regional Centre - 16769 Boul. Hymus, Kirkland, P.Q. H9H 3L4 Phones: (514) 630-6080 (800) 363-0137(P.Q.) Fax: (514) 630-6120

Invacare®, and Yes, you can., are trademarks of Invacare Corporation.

WWW.INVACARE.CA

NOTE: All specifications and dimensions are approximate.
For the most current pricing information visit www.invacare.ca