



NAME _____ PHONE NO _____

ACCOUNT # _____ P.O. _____ DATE _____

ADDRESS _____

1 SYSTEM TYPE

- TILTSYS Tilt Only System..... \$6,547
- RECLSYS Recline Only System..... \$9,660
- TLRLSYS Tilt/Recline System..... \$11,541
- PTOKIT Power Tilt Only..... \$2,995

EXISTING CHAIR TYPE (PTO ONLY)

- STMRWD Existing Storm® RWD ASBA Seat..... No Charge
- PRTR2M4 Existing Pronto™ R2 w/ASBA Seat w/MKIV™..... No Charge
- PRTR2RII Existing Pronto R2 w/ASBA Seat w/MKRII..... \$406

REAR/FRONT WHEEL DRIVE

- RWD Rear Wheel Drive..... No Charge
- FWD Front Wheel Drive..... No Charge

USER WEIGHT LIMITS

- STDWT Under 250 Lbs (185 lbs. w/Vent)..... **Standard**
- 2 HDTB Heavy Duty Tarsys Ready Base..... \$512

VENT TRAYS

- 3 VENT Vent Tray..... \$1,141

5 BACK STYLE

- 4,6 CURABS Curved ABS Back..... **Standard**
- 4,7 PROFILE Profile Back..... \$257
- 8 ROUND Round 7/8" Back Canes Only No Charge
- 9 CURVABS Curved ABS Back- P.T.O \$806

SEAT WIDTH

- TD14 Seat Width 14" PTO Only..... No Charge
- TR15 Seat Width 15" PTO Only..... No Charge
- TR16 Seat Width 16"..... No Charge
- TR17 Seat Width 17"..... No Charge
- TR18 Seat Width 18"..... No Charge
- TR19 Seat Width 19"..... No Charge
- TR20 Seat Width 20"..... No Charge
- TD21 Seat Width 21"..... \$323
- TR22 Seat Width 22"..... \$323

SEAT DEPTH

- TAR16 Seat Depth 16"..... No Charge
- TAR17 Seat Depth 17"..... No Charge
- TAR18 Seat Depth 18"..... No Charge
- TAR19 Seat Depth 19"..... No Charge
- TAR20 Seat Depth 20"..... No Charge
- TAR21 Seat Depth 21"..... \$323
- TAR22 Seat Depth 22"..... \$323

POWER TILT ONLY (PTO) SEAT DEPTH

- PSU17 Powered Seat Depth 17" or Less..... No Charge
- PSU19 Powered Seat Depth 18" - 19"..... No Charge
- PSU22 Powered Seat Depth 20" - 22"..... \$247

SEAT TILT

- 0SA 0 Degree..... No Charge
- 5SA 5 Degree..... No Charge
- 10 10SA 10 Degree..... No Charge
- 10 15SA 15 Degree..... No Charge

FOOTNOTES

- 1 These choices include your powered seating system only. For retrofit onto Storm Arrow, Torque SP and RX RWD bases
- 2 Must have Heavy Duty Package for Storm Arrow or Torque SP. (Not with PTO)
- 3 Available only with long frame
- 4 When ordering a curved or profile back, user seat depth is reduced by 1".
- 5 See page 5 for back height and back type ordering information.
- 6 CURABS back style available in back heights 20"-25" on Tilt and 20"-28" on Recline and Tilt/Recline.
- 7 PROFILE back style available in heights 23" and 24" on Tilt and 25" and 26" on Recline and Tilt/Recline.
- 8 ROUND back style available in back heights of 20"-28" For Tilt systems only.
- 9 CURVABS back style available in back heights of 15"-24" on Power Tilt Only.
- 10 For Recline Systems only

5 BACK HEIGHT

<input type="checkbox"/>	TSBH15	Back Height 15".....	No Charge
<input type="checkbox"/>	TSBH16	Back Height 16".....	No Charge
<input type="checkbox"/>	TSBH17	Back Height 17".....	No Charge
<input type="checkbox"/>	TSBH18	Back Height 18".....	No Charge
<input type="checkbox"/>	TSBH19	Back Height 19".....	No Charge
<input type="checkbox"/>	TSBH20	Back Height 20".....	No Charge
<input type="checkbox"/>	TSBH21	Back Height 21".....	No Charge
<input type="checkbox"/>	TSBH22	Back Height 22".....	No Charge
<input type="checkbox"/>	TSBH23	Back Height 23".....	No Charge
<input type="checkbox"/>	TSBH24	Back Height 24".....	No Charge
<input type="checkbox"/>	TSBH25	Back Height 25".....	\$311
<input type="checkbox"/>	TSBH26	Back Height 26".....	\$311
<input type="checkbox"/>	TSBH27	Back Height 27".....	\$311
<input type="checkbox"/>	TSBH28	Back Height 28".....	\$311

TILT/RECLINER HEADREST

<input type="checkbox"/>	TAR127700	Curved Headrest.....	Standard
<input type="checkbox"/>	TAR128000	Two-Step Headrest.....	No Charge
<input type="checkbox"/>	T0022400	Large Headrest.....	\$97
<input type="checkbox"/>	T0022600	Small Headrest.....	\$97

ARM TYPE

11 CONVENTIONAL ARM STYLES

<input type="checkbox"/>	19FL	Full Length Adj. Hgt. Conv. Flip Up Left.....	\$174
<input type="checkbox"/>	19FR	Full Length Adj. Hgt. Conv. Flip-Up Right.....	\$174
<input type="checkbox"/>	29FL	Desk Length Adj. Hgt. Conv. Flip-Up Left.....	\$174
<input type="checkbox"/>	29FR	Desk Length Adj. Hgt. Conv. Flip-Up Right....	\$174

12,13 RECLINING ARM STYLES

<input type="checkbox"/>	RA19FL	Reclining Full Length Arms Left.....	No Charge
<input type="checkbox"/>	RA19FR	Reclining Full Length Arms Right.....	No Charge
<input type="checkbox"/>	RA29FL	Reclining Desk Length Arms Left.....	No Charge
<input type="checkbox"/>	RA29FR	Reclining Desk Length Arms Right.....	No Charge

POWER TILT ONLY (PTO) ARM STYLES

<input type="checkbox"/>	EAFULLT	Existing Arm - Full Length Left.....	No Charge
<input type="checkbox"/>	EAFULRT	Existing Arm - Full Length Right.....	No Charge
<input type="checkbox"/>	EADSKLT	Existing Arm - Desk Length Left.....	No Charge
<input type="checkbox"/>	EADSKRT	Existing Arm - Desk Length Right.....	No Charge

14 LATERALS

<input type="checkbox"/>	TMSLM	Multibar Swingaway Laterals Med. Pr.....	\$315
<input type="checkbox"/>	TMSLL	Multibar Swingaway Laterals Lg. Pr.....	\$315
<input type="checkbox"/>	TCSLM	Clamp on Swingaway Laterals Medium Pr....	\$315
<input type="checkbox"/>	TCSLL	Clamp on Swingaway Laterals Large Pr.....	\$315

OTTO BOCK ARM PADS

15	<input type="checkbox"/>	OB20400	One Piece Channel Armpads Pr.....	\$165
15	<input type="checkbox"/>	OB20500	Large Forearm Pads Pr.....	\$156
	<input type="checkbox"/>	OB20600	Medium Forearm Pads Pr.....	\$156
	<input type="checkbox"/>	OB20700	Small Forearm Pads Pr.....	\$156

OTTO BOCK HAND PADS

<input type="checkbox"/>	OB23000	Large Flat Hand Pads Pr.....	\$128
<input type="checkbox"/>	OB23300	Medium Flat Hand Pads Pr.....	\$128
<input type="checkbox"/>	OB16480	Cone Pads-Left.....	\$141
<input type="checkbox"/>	OB16490	Cone Pad-Right.....	\$141
<input type="checkbox"/>	OB21500	Palm Extensor-Left.....	\$152
<input type="checkbox"/>	OB21600	Palm Extensor-Right.....	\$152

OTTO BOCK MISCELLANEOUS ACCESSORIES

<input type="checkbox"/>	OB21700	Hand/Forearm Pad Strap Kit Pr.....	\$65
<input type="checkbox"/>	OB21300	Elevating Swivel Unit-Left.....	\$161
<input type="checkbox"/>	OB21400	Elevating Swivel Unit-Right.....	\$161
<input type="checkbox"/>	OB20200	Trough Calf Pad Pr.....	\$190
<input type="checkbox"/>	OB20300	Calf Pad Adapter Plate Pr.....	\$50

FOOTNOTES

- 11 For Tilt systems only.
- 12 Arms recline with the system.
- 13 For Recline and Tilt/Recline systems only.
- 14 Not available with Tilt or Round Canes
- 15 Must order full length arms.

Invacare® TARSYS 2G Power Tilt/Recline "Off Chair" Systems Only Accessories & Options

HANGER BRACKET TYPE

- 17 ASB Center Pivot Style Bracket..... No Charge
- 18 PSB Pin Style Bracket..... No Charge
- 19 GHB Genius Hanger Bracket..... No Charge
- NHB No Hanger Bracket..... No Charge

FRONT RIGGINGS

CENTRE PIVOT STYLE

- 20 60HD 60 ° Swingaway Footrests..... \$225
- 20 SWING914 60 ° Swingaway Footrests MFX..... \$351
- 20 70NHD 70 ° Swingaway Footrests..... \$225
- 20 90SWING 90 ° Swingaway Footrests..... \$231
- 20,21 70STAPER 70 ° Tapered Footrests..... \$262
- 20 70MFX 70 ° Swingaway Footrests MFX..... \$351
- AT5544 Elevating Legrests..... \$529

PIN STYLE

- PW93 Swingaway Footrests 13 1/2"-19"..... \$163
- 22 P93M Minimum Swingaway Footrests 6 1/4" - 10"..... \$214
- 20,22 PW93E Footrests 9 3/4" - 15 1/4"..... \$214
- 22 P904A Swingaway Elevating Legrests 16" - 19 3/4"..... \$333
- PW93ST Footrests 14 3/4" - 19 3/4"..... \$333
- PAL4A Smartleg Art. Legrests 17 1/4" - 21 1/2"..... \$366

ELEVATING LEGRESTS

- 13 GL18S Genius Legrests 18" Ext. Std. Footplates..... \$575
- GL18L Genius Legrests 18" Ext. Lg. Footplates..... \$615
- GL20S Genius Legrests 20" Ext. Std. Footplates..... \$575
- GL20L Genius Legrests 20" Ext. Lg. Footplates..... \$615
- AT5544P Elevating Legrests \$529

POWER ELEVATING LEGRESTS

- 24 ELRPW Power ELR with Articulation..... \$2,856

POWER LEGS CONTROL MODULE

- 25 TRCM Tilt Recline Control Module..... \$1,973
- 26 TRCM-R Tilt Recline Control Module - Recline..... \$1,973

27 FRAME FINISH

CONTEMPORARY

- 24P Wet Black..... No Charge
- 60P Silver Metallic..... No Charge
- 71P Silver Vein..... No Charge
- 105P Midnight Blue..... No Charge
- 115P Black Prism..... No Charge

TOO HOT

- 30P Sunny Yellow..... No Charge
- 61P Electric Red..... No Charge
- 117P Tangerine Red..... No Charge
- 119P Bubblegum Pink..... No Charge
- 120P Grape Madness..... No Charge

WAY COOL

- 62P Electric Blue..... No Charge
- 104P Electric Teal..... No Charge
- 121P Lolly Pop Blue..... No Charge
- 122P Cosmic Blue..... No Charge
- 127P Grasshopper Green..... No Charge

28 FOOTPLATE OPTIONS

- 29 1651 Flip-Up Composite Footplates..... No Charge
- 30 1350 Extra Large Aluminum Footplates..... \$99
- 1450 Aluminum Footplates-Deluxe..... \$96
- 31,32 AT5543 Adjustable Angle Flip-Up Footplates..... \$122
- 32 AT5043 Flip-Up Footplates..... No Charge
- 33 S1263 One Piece Wooden Footboard..... \$131

FOOTNOTES

- 16 Pricing reflects standard sizes. Seating Discounts apply. Check the Invacare Seating Price List for complete listing.
- 17 For Centre Pivot style front rigging only.
- 18 For Pin & Saddle front rigging only.
- 19 For Mechanical elevating front rigging only. Not available with Tilt only systems.
- 20 All footrests are standard with heel loops.
- 21 70STAPER not available with 9" casters and flip-up composite footplates on 19" or 20" widths, Impact guards are standard.
- 22 Only available on 16" or wider.
- 23 Mechanically Elevating Legrests elevate in conjunction with recline cycle. Available in Black only.
- 24 Independent powered elevating legrests includes elevating and Articulating legrests with angle adjustable footplates.: Extension range is 13"-20", Centre Pivot Hanger Bracket
- 25 Tilt/Recline and Recline only come standard with TRCM, which provides for operation through a 4-way toggle or driver control. TRCM upgrade needed with power leg pick on Tilt/Recline systems.
- 26 Needed if combined with TRCM elevating legs, Power upgrade for recline
- 27 Indicates colour of front rigging only. All elevating legrests are wet black only.
- 28 For FWD bases.
- 29 Standard on 60HD, SWING914, 70NHD, 70TAPER, and AT5544.
- 30 Not available w/70STAPER, 60MFX, and 70MFX.
- 31 Standard on 70MFX.
- 32 Only available on 90SWING and 70MFX.
- 33 Not available on 90SWING and 70MFX.

Back Height Specifications Chart

	20"	21"	22"	23"	24"	25"	26"	27"	28"
Back Heights									
Tilt Only Systems									
-Back Cane Height	20"	20"	20"	20"	20"	20"	N/A	N/A	N/A
-ABS Back Pan/Cushion Height	16"	16"	18"	19"	20"	21"	22"	N/A	N/A
-Profile Back	N/A	N/A	N/A	16"	18"	N/A	N/A	N/A	N/A
Recline Only and Tilt/Recline Systems									
-Back Cane Height	19"	19"	19 3/4"	20 3/4"	21 3/4"	22 3/4"	23 3/4"	24 3/4"	25 3/4"
-ABS Back Pan/Cushion Height	11"	12"	13"	14"	15"	16"	17"	18"	19"
-Profile Back	N/A	N/A	N/A	N/A	N/A	16"	18"	N/A	N/A

NOTE: All specifications and dimensions are approximate. For the most current pricing information visit www.invacare.ca.

Invacare® TARSYS 2G Power Tilt/Recline "Off Chair" Systems Only Accessories & Options

FOOTREST/LEGREST ACCESSORIES

34	<input type="checkbox"/>	1820	Non-Removable Device.....	\$36
34	<input type="checkbox"/>	1337	Calf Strap.....	\$59
34	<input type="checkbox"/>	1600BK	Heel Loops With Ankle Straps.....	\$79
34	<input type="checkbox"/>	1345	Toe Loops.....	\$61
34	<input type="checkbox"/>	1348	Legrest Panel.....	\$86
34	<input type="checkbox"/>	9342	"H" Strap.....	\$87
35	<input type="checkbox"/>	LPT2	Longer Pivot Slide Tube For 1350.....	\$122
35	<input type="checkbox"/>	IMPACTG	Impact Guards.....	No Charge
34	<input type="checkbox"/>	FPPP	Footplate Parallel With Floor.....	\$199
35	<input type="checkbox"/>	ALPT	Longer Pivot And Slide Tube.....	\$123
35	<input type="checkbox"/>	ALPT4	Longer Pivot And Slide Tube 4".....	\$123
34	<input type="checkbox"/>	1346	Heel Loops.....	\$59

FOOTNOTES

- 34 Only available on Pin & Saddle front riggings.
- 35 Only available on Centre Pivot style front riggings.
- 36 Not available with Recline or Tilt/Recline systems.
- 37 Informational only; to send correct wiring harness with the Seat.

ELECTRONIC ACCESSORIES

36	<input type="checkbox"/>	TAC	Two Actuator Controller.....	\$994
	<input type="checkbox"/>	TAR014800	Interface Cable Invacare Rim.....	\$246
	<input type="checkbox"/>	TAR015400	Interface Cable Peachtree.....	\$188
	<input type="checkbox"/>	SFTSTD	Single Function Toggle Std Mount.....	No Charge
	<input type="checkbox"/>	SFTQL	Single Function Toggle Quadlink.....	No Charge

SEAT POSITIONING STRAPS

<input type="checkbox"/>	1311BK	Airline Buckle Seat Positioning Strap.....	\$79
<input type="checkbox"/>	1515	Auto Style Buckle Seat Positioning Strap.....	Standard

CHEST POSITIONING STRAPS

<input type="checkbox"/>	7311BK	Airline Buckle Chest Positioning Strap.....	\$79
<input type="checkbox"/>	9515	Auto Style Chest Positioning Strap.....	\$79
<input type="checkbox"/>	7321BK	Hook & Loop Chest Positioning Strap.....	\$61

37 BATTERY BOXES

<input type="checkbox"/>	NN372	Group 24 Battery Boxes.....	No Charge
<input type="checkbox"/>	NN373	22NF Battery Boxes.....	No Charge

** Dimensions calculated on Invacare Storm Series® Bases				
Dimensions	FWD	RWD		
Seat Width*	16" - 22"	16" - 22"		
Seat Depth	16" - 22"	16" - 22"		
Overall Width	25.25"	24.75"		
Overall Length (from headtubes to active anti-tipper)		Short	Long	
With Vent Tray	37"	38.5"	41.5"	
Without Vent Tray	34"	31.5"	34.5"	
Seat-to-Floor Height (based on 16" depth)				
At 0°	19.5"	18.5"		
At 5°	20.25"	19.5"		
Weight Capacity - User				
With Vent Tray	185 lb. / 285 lb. with Heavy-Duty pkg.	185 lb. / 285 lb. with Heavy-Duty pkg.		
Without Vent Tray	250 lb. / 350 lb. with Heavy-Duty pkg.	250 lb. / 350 lb. with Heavy-Duty pkg.		
Product Weight				
Without Base	105 lb.	105 lb.		
With Base	300 lb.	260 lb.		
Electronics	Standard	Optional		
2GT™	Single Toggle	TAC		
2GR™	TAC	TRCM		
2GTR™	TRCM	n/a		
		Actual Back Material Height		
	Tilt Only	Profile	Recline/VSR	
Back Height	Curved ABS†	Profile	Curved ABS‡	Profile
20"	16.25"	n/a	11.25"	n/a
21"	17.25"	n/a	12.25"	n/a
22"	18.25"	n/a	13.25"	n/a
23"	19.25"	16"	14.25"	n/a
24"	20.25"	18"	15.25"	n/a
25"	21.25"	n/a	16.25"	16"
26"	22.25"	n/a	17.25"	18"

* 23" and 24" available through special order.
† On the 2GT™, there is a 3.75" gap from the seat pan to the bottom of the back. This allows for the cushion and cannot be reduced.
‡ On the 2GTR™, there is an 8.75" gap from the seat pan to the bottom of the back cushion. This allows for 5" of VSR movement plus 3.75" to the back-cane pivot point. This can be reduced by ordering a special longer back but will limit the VSR movement.

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEAT AND BACK CUSHIONS
SEE PRICE LIST FORM # 05-338C.**

Extended Warranty Information:

Electrical Components: **2 Years**
Motors (2 and 4-pole): **2 Years**
GB Motors: **5 Years**

* Warranty valid to original purchaser who purchases this product within Canada

* For more product warranty information, please refer to the product owner's manual

*** UPDATED - 06/18/2007 - New Warranty Box**

TARSYS TECH TIPS

Tilt Only Systems

Standard Electronics:

Single function Toggle-direct operation of tilt

Through Driver Controls:

Option 1

Upgrade to Two Actuator Control (TAC)

Provides seamless integration and Many critical programming parameters including) drive lockout

Option 2

Upgrade to TRCM

Required if using powered elevating legs

Recline Only and Tilt/Recline Systems

Option 1

Standard Electronics:

Two Actuator Control(TRCM) (TRCM-R)

Can be Operated through 4 Way Toggle or Driver Control w/o additional pieces.

(MKIVA electronics needed)

Option 2

TRCM-R Needed if combined with powered elevating leg rests.

With a recline only system.

Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4G 4Z4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)

Regional Centre - 16769 Boul. Hymus, Kirkland, P.Q. H9H 3L4 Phones: (514) 630-6080 (800) 363-0137(P.Q.) Fax: (514) 630-6120

Invacare®, and Yes, you can., are trademarks of Invacare Corporation.

WWW.INVACARE.CA