

JAZZY® 1103 ULTRA

Jazzy® 1103 Ultra shown in Candy Apple Red with gray contoured medium-back seat, foot platform and VSI controller.




The Jazzy® 1103 Ultra brings sporty style and high-performance features, like a standard power elevating seat and Active-Trac® Suspension, together into one perfectly balanced power chair. Yet, even fully loaded, the Jazzy 1103 Ultra remains highly compact and maneuverable.

SPECIFICATIONS

Weight capacity:	300 lbs.
Maximum speed:	Up to 4 mph
Ground clearance:	2.5"
Turning radius:	18.5"
Overall length:	35.5"
Overall width:	23"
Wheels:	Drive: 10" pneumatic Casters: 8" solid Anti-tips: 5" solid
Suspension:	Active-Trac®
Drive train:	(2) motor, mid-wheel drive
Braking system:	Intelligent braking (electronic, regenerative disc brakes)
Std. Electronics:	50 amp, 6-key VSI
Specialty controls:	Yes
Battery charger:	On-board, 4 amp
Per-charge range:	Up to 25 miles
Battery requirements:	(2) 12 volt deep cycle
Battery size:	U-1, 32AH
Transport weight:	
Batteries:	24.5 lbs. each
Base:	97 lbs.
Standard seat:	33 lbs. (medium-back)
Seat dimensions:	20" x 20" (max. contoured seating)
Warranty:	Lifetime limited warranty on frame 2-year limited warranty on electronics 18-months limited warranty on drive motors

Colors:



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

FEATURES

- 300 lbs. weight capacity
- Standard power elevating seat provides greater access to objects out of reach
- Powerful two-motor, patented mid-wheel drive design offers outstanding maneuverability
- Low center of gravity design delivers a very smooth and stable ride
- 10" knobby drive wheels and optimum weight distribution provide superior traction
- Speed range adjustable up to 4 mph
- Standard lap belt
- On-board battery charger

OPTIONS

- Contoured high-back seat
- Synergy™ Seating System
- Solid seat pan
- Solid tires
- TRU-Balance™ Power Positioning System*
- Cup holder
- Heavy-duty dust cover
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Rear basket
- Angle-adjustable foot plates
- Saddlebag
- Off-board charger

* Requires off-board charger

Pride
Mobility Products Corp