

Jazzy 1113[®]

and 1113 ATS

The Jazzy[®] 1113 offers you everything you need—all of the “big chair” options and accessories, even full rehab capabilities—on the most compact frame available. You’ll take on tight spaces and narrow hallways without sacrificing a strong outdoor performance never before seen in a chair of this size.

For sharp ramp transitions as well as a smooth and stable ride even over exceptionally rough and uneven terrain, the Jazzy 1113 ATS is equipped with Active-Trac[®] Suspension. As an added plus, both the Jazzy 1113 and 1113 ATS fully disassemble for transport.




Jazzy 1113
shown in Candy Apple Red
with gray contoured high-back
seat, foot platform and PG VSI
controller.



Weight capacity:	300 lbs.*	
Maximum speed:	Up to 4 mph	
Ground clearance:	1113: 2"	1113 ATS: 2.5"
Turning radius:	1113: 16.5"	1113 ATS: 18"
Overall length:	1113: 32"	1113 ATS: 36.5"
Overall width:	23"	
Wheels:	Drive: 10" pneumatic	
	Casters: 8" solid	
	Anti-tips: 6" solid	
Suspension:	1113: Limited	
	1113 ATS: Active-Trac®	
Drive train:	Two-motor, mid-wheel drive	
Braking system:	Intelligent braking (electronic, regenerative disc brakes)	
Std. Electronics:	50 amp, PG VSI controller	
Specialty controls:	Yes	
Battery charger:	On-board, 4 amp (standard)	
Per-charge range:	Up to 25 miles	
Battery requirements:	(2) 12 volt deep cycle	
Battery size:	32AH, U-1	
Transport weight		
Battery weight:	24.5 lbs. each	
Base weight:	1113: 79.5 lbs. 1113 ATS: 90.5 lbs.	
Std. seat weight:	33 lbs. (medium-back)	
Seat dimensions:	20" x 20" (max. contoured seating)	
Warranty:	Lifetime limited warranty on frame 2-year limited warranty on electronics 18-months limited warranty on drive motors	

Colors:



FEATURES

- Powerful two-motor, mid-wheel drive design offers outstanding maneuverability
- Low center of gravity design delivers a very smooth and stable ride
- 10" knobby drive wheels and optimum weight distribution provide superior traction
- Speed range adjustable up to 4 mph
- On-board battery charger
- Fully programmable controller
- Angle-adjustable footplate
- 6" resilient front anti-tips

OPTIONS

- TRU-Balance™ Power Tilt/ Super Low Tilt
- TRU-Balance™ Power Positioning Systems
- Contoured high-back seat
- Synergy™ seating system
- Contoured medium/high-back with solid seat pan
- Lap belt
- Cup holder
- Heavy-duty dust cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Solid tires
- Adjustable-height armrests (optional on standard seats)
- Removable height-adjustable arms (optional on Synergy Seat)
- Swing-away joystick mount
- Safety flag

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

* Maximum weight capacity is only 200 lbs. when used with the TRU-Balance Power Positioning Systems

