

# JAZZY<sup>®</sup> 1121



**JAZZY 1121** shown in Viper Blue with optional contoured high-back seat, PG VSI controller and foot platform

The Jazzy<sup>®</sup> 1121 with enhanced Active-Trac<sup>®</sup> Suspension (ATS) makes easy work of rough and uneven outdoor terrain. Plus, an articulating rear caster beam and patented resilient, adjustable front anti-tips further serve to add to the 1121's comfort and stability. Even tight-space indoor maneuverability is made easy with the 1121's compact overall size and mid-wheel drive design. The Jazzy 1121 with ATS is the perfect power chair for your diverse and demanding lifestyle.



800-800-8586 (US) • 888-570-1113 (Canada) • [pridemobility.com](http://pridemobility.com)

<b>Weight capacity:</b>	300 lbs.
<b>Maximum speed:</b>	Up to 5 mph
<b>Ground clearance:</b>	3.5"
<b>Turning radius:</b>	20.5"
<b>Overall length:</b>	38.5"
<b>Overall width:</b>	23.8"
<b>Wheels:</b>	Drive: 14" pneumatic
	Casters: 8" solid
	Anti-tips: 6" solid
<b>Suspension:</b>	Active-Trac®
<b>Drive train:</b>	Two-motor, mid-wheel drive
<b>Braking system:</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Std. Electronics:</b>	50 amp, PG VSI controller
<b>Specialty controls:</b>	Yes
<b>Battery charger:</b>	On-board, 5 amp (standard)
<b>Per-charge range:</b>	Up to 25 miles
<b>Battery requirements:</b>	(2) 12 volt, deep cycle
<b>Battery size:</b>	NF-22
<b>Transport weight</b>	
<b>Battery weight:</b>	38 lbs. each
<b>Base weight:</b>	135.5 lbs.
<b>Std. seat weight:</b>	33 lbs. (medium-back)
<b>Seat dimensions:</b>	20" x 22" (max. contoured seating)
<b>Warranty:</b>	Lifetime limited warranty on frame 2-year limited warranty on electronics 18-months limited warranty on drive motors

**Colors:**



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

# FEATURES

- Active-Trac® Suspension for enhanced performance on rough and uneven terrain, creating a smooth and stable ride
- Rear articulating caster beam for great performance and durability
- One-piece, 2" diameter tubular frame for maximum strength and enhanced performance
- Easy front access to batteries
- Easy accessible side charger port
- PG VSI integral controller

# OPTIONS

- Contoured high-back
- Solid seat pan
- Synergy™ Seat
- Synergy™ Manual Recline
- Synergy™ Manual Tilt
- TRU-Balance™ Power Positioning System
- Solid tires
- Lap belt
- Cup holder
- Heavy-duty dust cover
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Power elevating seat
- Rear basket
- Angle-adjustable foot plates
- 70° and 80° heavy-duty drop-in footrests

**HD PACKAGE:**

- 400 lbs. weight capacity
- Solid tires (required)
- 70 amp, PG controller
- Max. speed up to 4 mph
- 24.2" turning radius
- Base weight of 144.5 lbs.
- 41" length