

JAZZY[®] 1122



JAZZY 1122 shown in Candy Apple Red with deluxe contoured high-back seat, PG Remote Plus Controller and foot platform



At first glance, the Jazzy[®] 1122 catches your eye with its sleek styling and sporty edge. Yet, the 1122 is so much more. The 1122 features standard Active-Trac[®] and rear caster suspension and a high-speed/high-torque motor package that allows an optimal speed of 6 mph and a weight capacity of 300 lbs. In addition to these impressive performance-enhancing features, the 1122 also offers easy battery access through the front of the chair without the removal of the shroud or seating system and a height and depth-adjustable solid foot platform. The Jazzy 1122 is truly the most advanced power chair available today!



800-800-8586 (US) • 888-570-1113 (Canada) • pridemobility.com

SPECIFICATIONS

Weight capacity:	300 lbs.
Maximum speed:	Up to 6 mph
Ground clearance:	3.25"
Turning radius:	21"
Overall length:	40.5"
Overall width:	25.4"
Wheels:	Drive: 14" pneumatic Casters: 8" solid Anti-tips: 6" solid
Suspension:	Active-Trac® and rear caster
Drive train:	Two-motor, mid-wheel drive; high-speed motor package (std.)
Braking system:	Intelligent braking (electronic, regenerative disc brakes)
Std. Electronics:	70 amp PG Pilot controller
Specialty controls:	Yes
Battery charger:	On-board, 5 amp (standard)
Per-charge range:	Up to 25 miles
Battery requirements:	(2) 12 volt, deep cycle
Battery size:	70AH, Group 24
Transport weight	
Battery weight:	53.5 lbs. each
Base weight:	129 lbs.
Std. seat weight:	40 lbs. (high-back)
Seat dimensions:	22" x 22" (max. contoured seating)
Warranty:	Lifetime limited warranty on frame 2-year limited warranty on electronics 18-months limited warranty on drive motors

Colors:



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

FEATURES

- High-speed motor package—top speed of 6 mph
- Easy front-access batteries
- Powerful, patented, two-motor, mid-wheel drive design provides superb power and acceleration
- Active-Trac® and rear caster suspension for outstanding outdoor performance
- Adjustable front anti-tips to meet your preferences and terrain challenges
- On-board battery charger
- Battery-conserving sleep mode available with Remote Plus controller

OPTIONS

- Synergy™ Seat
- Synergy™ Manual Recline
- Synergy™ Manual Tilt
- TRU-Balance™ Super Low¹
- TRU-Balance™ Power Tilt
- Solid seat pan
- Power elevating seat
- Solid tires
- Lap belt
- Cup holder
- Heavy-duty dust cover
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Rear basket
- Specialty controls
- Angle-adjustable footplates
- 70° and 80° heavy-duty drop-in footrests
- Saddlebag
- Off-board charger

¹ Requires off-board charger