

# Jazzy®

## 1133

The sporty Jazzy® 1133 delivers high-performance, hot style, rehab capability, and patented features, like Active-Trac® Suspension, in a full-size, yet easily transportable package that conveniently disassembles into five manageable base pieces. All these features combine with true rigid-frame stability to bring innovation to the world of full-sized transportable power chairs. Simply put, the Jazzy 1133 is the perfect full-size, take-along power chair choice for people on the go.



### DISASSEMBLY

*Quickly and easily disassembles*



### VERSATILITY

*Rehab capabilities*



### PERFORMANCE

*Active-Trac® Suspension for excellent outdoor performance*



### STABILITY


*Full-sized rigid frame stability*

**Pride**<sup>®</sup>  
Mobility Products Corp.

800-800-8586 (US)  
888-570-1113 (Canada)  
pridemobility.com

<b>Weight capacity:</b>	300 lbs.
<b>Maximum speed:</b>	Up to 5 mph
<b>Ground clearance:</b>	3.5"
<b>Turning radius:</b>	21.5"
<b>Overall length:</b>	41"
<b>Overall width:</b>	25.5"
<b>Wheels:</b>	<b>Drive:</b> 14" pneumatic (standard)
	<b>Casters:</b> 8" solid
	<b>Anti-tips:</b> 6" solid
<b>Suspension:</b>	Active-Trac®
<b>Drive train:</b>	Two-motor, mid-wheel drive
<b>Braking system:</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Std. Electronics:</b>	50 amp, PG VSI controller
<b>Specialty controls:</b>	N/A
<b>Battery charger:</b>	On-board, 4 amp (standard)
<b>Per-charge range:</b>	Up to 25 miles
<b>Battery requirements:</b>	(2) 12 volt deep cycle
<b>Battery size:</b>	NF-22
<b>Transport weight</b>	
<b>Battery weight:</b>	38 lbs. each
<b>Base weight:</b>	95 lbs.
<b>Std. seat weight:</b>	37 lbs. (medium-back)
<b>Seat dimensions:</b>	20" x 22" (max. contoured seating)
<b>Warranty:</b>	Lifetime limited warranty on frame 2-year limited warranty on electronics 18-months limited warranty on drive motors

**Colors:**



Easy disassembly for quick and convenient transport



# FEATURES

- Easily disassembles for transport
- Powerful two-motor, mid-wheel drive design offers outstanding maneuverability
- Active-Trac® Suspension for outstanding performance
- Adjustable front anti-tips for personal preferences and terrain challenges
- Speed range adjustable up to 5 mph
- On-board battery charger
- Fully programmable controller
- Easy rear battery access

# OPTIONS

- Synergy® seating
- Contoured high-back seat
- Lap belt
- Cup holder
- Dust cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Flat-free inserts
- Adjustable-height armrests
- Flat-free tires
- Saddlebag

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.